

TimeTagger

2.9.0.0

Generated by Doxygen 1.8.11



# Contents

<b>1</b>	<b>TimeTagger</b>	<b>1</b>
<b>2</b>	<b>Deprecated List</b>	<b>3</b>
<b>3</b>	<b>Module Index</b>	<b>5</b>
3.1	Modules . . . . .	5
<b>4</b>	<b>Hierarchical Index</b>	<b>7</b>
4.1	Class Hierarchy . . . . .	7
<b>5</b>	<b>Class Index</b>	<b>9</b>
5.1	Class List . . . . .	9
<b>6</b>	<b>File Index</b>	<b>11</b>
6.1	File List . . . . .	11
<b>7</b>	<b>Module Documentation</b>	<b>13</b>
7.1	base iterators . . . . .	13
7.1.1	Detailed Description . . . . .	14

<b>8</b>	<b>Class Documentation</b>	<b>15</b>
8.1	Coincidence Class Reference	15
8.1.1	Detailed Description	16
8.1.2	Constructor & Destructor Documentation	16
8.1.2.1	Coincidence(TimeTaggerBase *tagger, std::vector< channel_t > channels, timestamp_t coincidenceWindow=1000, CoincidenceTimestamp timestamp=CoincidenceTimestamp::Last)	16
8.1.3	Member Function Documentation	16
8.1.3.1	getChannel()	16
8.2	Coincidences Class Reference	16
8.2.1	Detailed Description	17
8.2.2	Constructor & Destructor Documentation	18
8.2.2.1	Coincidences(TimeTaggerBase *tagger, std::vector< std::vector< channel_t >> coincidenceGroups, timestamp_t coincidenceWindow, CoincidenceTimestamp timestamp=CoincidenceTimestamp::Last)	18
8.2.2.2	~Coincidences()	18
8.2.3	Member Function Documentation	18
8.2.3.1	getChannels()	18
8.2.3.2	next_impl(std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override	18
8.2.3.3	setCoincidenceWindow(timestamp_t coincidenceWindow)	18
8.2.4	Friends And Related Function Documentation	19
8.2.4.1	CoincidencesImpl	19
8.3	Combiner Class Reference	19
8.3.1	Detailed Description	20
8.3.2	Constructor & Destructor Documentation	20
8.3.2.1	Combiner(TimeTaggerBase *tagger, std::vector< channel_t > channels)	20
8.3.2.2	~Combiner()	20
8.3.3	Member Function Documentation	20
8.3.3.1	clear_impl() override	20
8.3.3.2	getChannel()	20
8.3.3.3	getData(std::function< int64_t *(size_t)> array_out)	21

8.3.3.4	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time) override</code>	21
8.3.4	Friends And Related Function Documentation	21
8.3.4.1	<code>CombinerImpl</code>	21
8.4	<code>ConstantFractionDiscriminator</code> Class Reference	21
8.4.1	Detailed Description	22
8.4.2	Constructor & Destructor Documentation	22
8.4.2.1	<code>ConstantFractionDiscriminator(TimeTaggerBase *tagger, std::vector&lt; channel_t &gt; channels, timestamp_t search_window)</code>	22
8.4.2.2	<code>~ConstantFractionDiscriminator()</code>	23
8.4.3	Member Function Documentation	23
8.4.3.1	<code>getChannels()</code>	23
8.4.3.2	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time) override</code>	23
8.4.3.3	<code>on_start() override</code>	23
8.4.4	Friends And Related Function Documentation	24
8.4.4.1	<code>ConstantFractionDiscriminatorImpl</code>	24
8.5	<code>Correlation</code> Class Reference	24
8.5.1	Detailed Description	25
8.5.2	Constructor & Destructor Documentation	25
8.5.2.1	<code>Correlation(TimeTaggerBase *tagger, channel_t channel_1, channel_t channel_2=CHANNEL_UNUSED, timestamp_t binwidth=1000, int n_bins=1000)</code>	25
8.5.2.2	<code>~Correlation()</code>	25
8.5.3	Member Function Documentation	25
8.5.3.1	<code>clear_impl() override</code>	25
8.5.3.2	<code>getData(std::function&lt; int32_t *(size_t)&gt; array_out)</code>	26
8.5.3.3	<code>getDataNormalized(std::function&lt; double *(size_t)&gt; array_out)</code>	26
8.5.3.4	<code>getIndex(std::function&lt; long long *(size_t)&gt; array_out)</code>	26
8.5.3.5	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time) override</code>	26
8.5.4	Friends And Related Function Documentation	26
8.5.4.1	<code>CorrelationImpl</code>	26
8.6	<code>CountBetweenMarkers</code> Class Reference	27

8.6.1	Detailed Description	28
8.6.2	Constructor & Destructor Documentation	28
8.6.2.1	CountBetweenMarkers(TimeTaggerBase *tagger, channel_t click_channel, channel_t begin_channel, channel_t end_channel=CHANNEL_UNUSED, int32_t n_values=1000)	28
8.6.2.2	~CountBetweenMarkers()	28
8.6.3	Member Function Documentation	28
8.6.3.1	clear_impl() override	28
8.6.3.2	getBinWidths(std::function< long long *(size_t)> array_out)	28
8.6.3.3	getData(std::function< int32_t *(size_t)> array_out)	28
8.6.3.4	getIndex(std::function< long long *(size_t)> array_out)	29
8.6.3.5	next_impl(std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override	29
8.6.3.6	ready()	29
8.6.4	Friends And Related Function Documentation	29
8.6.4.1	CountBetweenMarkersImpl	29
8.7	Counter Class Reference	29
8.7.1	Detailed Description	30
8.7.2	Constructor & Destructor Documentation	31
8.7.2.1	Counter(TimeTaggerBase *tagger, std::vector< channel_t > channels, timestamp_t binwidth=1000000000, int32_t n_values=1)	31
8.7.2.2	~Counter()	31
8.7.3	Member Function Documentation	31
8.7.3.1	clear_impl() override	31
8.7.3.2	getData(std::function< int32_t *(size_t, size_t)> array_out)	31
8.7.3.3	getIndex(std::function< long long *(size_t)> array_out)	31
8.7.3.4	next_impl(std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override	31
8.7.3.5	on_start() override	32
8.7.4	Friends And Related Function Documentation	32
8.7.4.1	CounterImpl	32
8.8	Countrate Class Reference	32
8.8.1	Detailed Description	33

8.8.2	Constructor & Destructor Documentation . . . . .	33
8.8.2.1	Countrate(TimeTaggerBase *tagger, std::vector< channel_t > channels) . . . . .	33
8.8.2.2	~Countrate() . . . . .	34
8.8.3	Member Function Documentation . . . . .	34
8.8.3.1	clear_impl() override . . . . .	34
8.8.3.2	getCountsTotal(std::function< int64_t *(size_t)> array_out) . . . . .	34
8.8.3.3	getData(std::function< double *(size_t)> array_out) . . . . .	34
8.8.3.4	next_impl(std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override . . . . .	34
8.8.3.5	on_start() override . . . . .	35
8.8.4	Friends And Related Function Documentation . . . . .	35
8.8.4.1	CountrateImpl . . . . .	35
8.9	CustomLogger Class Reference . . . . .	35
8.9.1	Constructor & Destructor Documentation . . . . .	35
8.9.1.1	CustomLogger() . . . . .	35
8.9.1.2	~CustomLogger() . . . . .	35
8.9.2	Member Function Documentation . . . . .	35
8.9.2.1	disable() . . . . .	35
8.9.2.2	enable() . . . . .	35
8.9.2.3	Log(int level, const std::string &msg)=0 . . . . .	35
8.10	CustomMeasurementBase Class Reference . . . . .	36
8.10.1	Constructor & Destructor Documentation . . . . .	37
8.10.1.1	CustomMeasurementBase(TimeTaggerBase *tagger) . . . . .	37
8.10.1.2	~CustomMeasurementBase() override . . . . .	37
8.10.2	Member Function Documentation . . . . .	37
8.10.2.1	_lock() . . . . .	37
8.10.2.2	_unlock() . . . . .	37
8.10.2.3	clear_impl() override . . . . .	37
8.10.2.4	finalize_init() . . . . .	37
8.10.2.5	is_running() const . . . . .	37

8.10.2.6	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time) override</code>	37
8.10.2.7	<code>next_impl_cs(void *tags_ptr, uint64_t num_tags, timestamp_t begin_time, timestamp_t end_time)</code>	38
8.10.2.8	<code>on_start() override</code>	38
8.10.2.9	<code>on_stop() override</code>	38
8.10.2.10	<code>register_channel(channel_t channel)</code>	38
8.10.2.11	<code>stop_all_custom_measurements()</code>	38
8.10.2.12	<code>unregister_channel(channel_t channel)</code>	38
8.11	DelayedChannel Class Reference	38
8.11.1	Detailed Description	39
8.11.2	Constructor & Destructor Documentation	39
8.11.2.1	<code>DelayedChannel(TimeTaggerBase *tagger, channel_t input_channel, timestamp_t delay)</code>	39
8.11.2.2	<code>DelayedChannel(TimeTaggerBase *tagger, std::vector&lt; channel_t &gt; input_channels, timestamp_t delay)</code>	40
8.11.2.3	<code>~DelayedChannel()</code>	40
8.11.3	Member Function Documentation	40
8.11.3.1	<code>getChannel()</code>	40
8.11.3.2	<code>getChannels()</code>	40
8.11.3.3	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time) override</code>	40
8.11.3.4	<code>on_start() override</code>	41
8.11.3.5	<code>setDelay(timestamp_t delay)</code>	41
8.11.4	Friends And Related Function Documentation	41
8.11.4.1	<code>DelayedChannelImpl</code>	41
8.12	Dump Class Reference	41
8.12.1	Detailed Description	42
8.12.2	Constructor & Destructor Documentation	42
8.12.2.1	<code>Dump(TimeTaggerBase *tagger, std::string filename, int64_t max_tags, std::vector&lt; channel_t &gt; channels=std::vector&lt; channel_t &gt;())</code>	42
8.12.2.2	<code>~Dump()</code>	43
8.12.3	Member Function Documentation	43



8.12.3.1	<code>clear_impl()</code> override	43
8.12.3.2	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time)</code> override	43
8.12.3.3	<code>on_start()</code> override	43
8.12.3.4	<code>on_stop()</code> override	43
8.12.4	Friends And Related Function Documentation	44
8.12.4.1	<code>DumpImpl</code>	44
8.13	Event Struct Reference	44
8.13.1	Member Data Documentation	44
8.13.1.1	<code>state</code>	44
8.13.1.2	<code>time</code>	44
8.14	EventGenerator Class Reference	44
8.14.1	Detailed Description	45
8.14.2	Constructor & Destructor Documentation	45
8.14.2.1	<code>EventGenerator(TimeTaggerBase *tagger, channel_t trigger_channel, std::vector&lt; timestamp_t &gt; pattern, uint64_t trigger_divider=1, uint64_t divider_offset=0, channel_t stop_channel=CHANNEL_UNUSED)</code>	45
8.14.2.2	<code>~EventGenerator()</code>	46
8.14.3	Member Function Documentation	46
8.14.3.1	<code>clear_impl()</code> override	46
8.14.3.2	<code>getChannel()</code>	46
8.14.3.3	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time)</code> override	46
8.14.3.4	<code>on_start()</code> override	46
8.14.4	Friends And Related Function Documentation	47
8.14.4.1	<code>EventGeneratorImpl</code>	47
8.15	FastBinning Class Reference	47
8.15.1	Detailed Description	47
8.15.2	Member Enumeration Documentation	47
8.15.2.1	<code>Mode</code>	47
8.15.3	Constructor & Destructor Documentation	48
8.15.3.1	<code>FastBinning()</code>	48

8.15.3.2	FastBinning(uint64_t divisor, uint64_t max_duration_)	48
8.15.4	Member Function Documentation	48
8.15.4.1	divide(uint64_t duration) const	48
8.15.4.2	getMode() const	48
8.16	FileReader Class Reference	48
8.16.1	Detailed Description	48
8.16.2	Constructor & Destructor Documentation	48
8.16.2.1	FileReader(std::vector< std::string > filenames)	48
8.16.2.2	FileReader(const std::string &filename)	49
8.16.2.3	~FileReader()	49
8.16.3	Member Function Documentation	49
8.16.3.1	getConfiguration()	49
8.16.3.2	getData(uint64_t n_events)	49
8.16.3.3	getDataRaw(std::vector< Tag > &tag_buffer)	49
8.16.3.4	getLastMarker()	50
8.16.3.5	hasData()	50
8.16.4	Friends And Related Function Documentation	50
8.16.4.1	FileReaderImpl	50
8.17	FileWriter Class Reference	50
8.17.1	Detailed Description	51
8.17.2	Constructor & Destructor Documentation	51
8.17.2.1	FileWriter(TimeTaggerBase *tagger, const std::string &filename, std::vector< channel_t > channels)	51
8.17.2.2	~FileWriter()	52
8.17.3	Member Function Documentation	52
8.17.3.1	clear_impl() override	52
8.17.3.2	getMaxFileSize()	52
8.17.3.3	getTotalEvents()	52
8.17.3.4	getTotalSize()	52
8.17.3.5	next_impl(std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override	52

8.17.3.6	<code>on_start()</code> override	53
8.17.3.7	<code>on_stop()</code> override	53
8.17.3.8	<code>setMarker(const std::string &amp;marker)</code>	53
8.17.3.9	<code>setMaxFileSize(uint64_t max_file_size)</code>	53
8.17.3.10	<code>split(const std::string &amp;new_filename="")</code>	54
8.17.4	Friends And Related Function Documentation	54
8.17.4.1	<code>FileWriterImpl</code>	54
8.18	Flim Class Reference	54
8.18.1	Detailed Description	56
8.18.2	Constructor & Destructor Documentation	56
8.18.2.1	<code>Flim(TimeTaggerBase *tagger, channel_t start_channel, channel_t click_← channel, channel_t pixel_begin_channel, uint32_t n_pixels, uint32_t n_bins, timestamp_t binwidth, channel_t pixel_end_channel=CHANNEL_UNUSED, channel_t frame_begin_channel=CHANNEL_UNUSED, uint32_t finish_after_← _outputframe=0, uint32_t n_frame_average=1, bool pre_initialize=true)</code>	56
8.18.2.2	<code>~Flim()</code>	57
8.18.3	Member Function Documentation	57
8.18.3.1	<code>clear_impl()</code> override	57
8.18.3.2	<code>frameReady(uint32_t frame_number, std::vector&lt; uint32_t &gt; &amp;data, std::vector&lt; timestamp_t &gt; &amp;pixel_begin_times, std::vector&lt; timestamp_t &gt; &amp;pixel_end_← times, timestamp_t frame_begin_time, timestamp_t frame_end_time)</code>	57
8.18.3.3	<code>get_ready_index(int32_t index)</code>	57
8.18.3.4	<code>getCurrentFrame(std::function&lt; uint32_t *(size_t, size_t)&gt; array_out)</code>	57
8.18.3.5	<code>getCurrentFrameEx()</code>	57
8.18.3.6	<code>getCurrentFrameIntensity(std::function&lt; float *(size_t)&gt; array_out)</code>	57
8.18.3.7	<code>getFramesAcquired() const</code>	57
8.18.3.8	<code>getIndex(std::function&lt; long long *(size_t)&gt; array_out)</code>	57
8.18.3.9	<code>getReadyFrame(std::function&lt; uint32_t *(size_t, size_t)&gt; array_out, int32_← t index=-1)</code>	57
8.18.3.10	<code>getReadyFrameEx(int32_t index=-1)</code>	58
8.18.3.11	<code>getReadyFrameIntensity(std::function&lt; float *(size_t)&gt; array_out, int32_← t index=-1)</code>	58
8.18.3.12	<code>getSummedFrames(std::function&lt; uint32_t *(size_t, size_t)&gt; array_out, bool only_ready_frames=true, bool clear_summed=false)</code>	58

8.18.3.13	getSummedFramesEx(bool only_ready_frames=true, bool clear_summed=false)	59
8.18.3.14	getSummedFramesIntensity(std::function< float *(size_t)> array_out, bool only_ready_frames=true, bool clear_summed=false)	59
8.18.3.15	initialize()	59
8.18.3.16	on_frame_end() override	59
8.18.4	Member Data Documentation	59
8.18.4.1	accum_diffs	59
8.18.4.2	back_frames	59
8.18.4.3	captured_frames	59
8.18.4.4	frame_begins	60
8.18.4.5	frame_ends	60
8.18.4.6	last_frame	60
8.18.4.7	pixels_completed	60
8.18.4.8	summed_frames	60
8.18.4.9	swap_chain_lock	60
8.18.4.10	total_frames	60
8.19	FlimAbstract Class Reference	60
8.19.1	Constructor & Destructor Documentation	61
8.19.1.1	FlimAbstract(TimeTaggerBase *tagger, channel_t start_channel, channel_t click_channel, channel_t pixel_begin_channel, uint32_t n_pixels, uint32_t n_bins, timestamp_t binwidth, channel_t pixel_end_channel=CHANNEL_UNUSED, channel_t frame_begin_channel=CHANNEL_UNUSED, uint32_t finish_after_outputframe=0, uint32_t n_frame_average=1, bool pre_initialize=true)	61
8.19.1.2	~FlimAbstract()	62
8.19.2	Member Function Documentation	62
8.19.2.1	clear_impl() override	62
8.19.2.2	isAcquiring() const	62
8.19.2.3	next_impl(std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override	62
8.19.2.4	on_frame_end()=0	63
8.19.2.5	on_start() override	63
8.19.2.6	process_tags(const std::vector< Tag > &incoming_tags)	63
8.19.3	Member Data Documentation	63

8.19.3.1	acquiring	63
8.19.3.2	acquisition_lock	63
8.19.3.3	binner	63
8.19.3.4	binwidth	63
8.19.3.5	click_channel	63
8.19.3.6	current_frame_begin	63
8.19.3.7	current_frame_end	63
8.19.3.8	data_base	63
8.19.3.9	finish_after_outputframe	64
8.19.3.10	frame	64
8.19.3.11	frame_acquisition	64
8.19.3.12	frame_begin_channel	64
8.19.3.13	frames_completed	64
8.19.3.14	initialized	64
8.19.3.15	n_bins	64
8.19.3.16	n_frame_average	64
8.19.3.17	n_pixels	64
8.19.3.18	pixel_acquisition	64
8.19.3.19	pixel_begin_channel	64
8.19.3.20	pixel_begins	64
8.19.3.21	pixel_end_channel	64
8.19.3.22	pixel_ends	64
8.19.3.23	pixels_processed	64
8.19.3.24	previous_starts	64
8.19.3.25	start_channel	64
8.19.3.26	ticks	64
8.19.3.27	time_window	64
8.20	FlimBase Class Reference	65
8.20.1	Constructor & Destructor Documentation	65

8.20.1.1	FlimBase(TimeTaggerBase *tagger, channel_t start_channel, channel_t click↔ _channel, channel_t pixel_begin_channel, uint32_t n_pixels, uint32_t n_bins, timestamp_t binwidth, channel_t pixel_end_channel=CHANNEL_UNUSED, channel_t frame_begin_channel=CHANNEL_UNUSED, uint32_t finish_after↔ _outputframe=0, uint32_t n_frame_average=1, bool pre_initialize=true)	65
8.20.1.2	~FlimBase()	66
8.20.2	Member Function Documentation	66
8.20.2.1	frameReady(uint32_t frame_number, std::vector< uint32_t > &data, std::vector< timestamp_t > &pixel_begin_times, std::vector< timestamp_t > &pixel_end_↔ times, timestamp_t frame_begin_time, timestamp_t frame_end_time)	66
8.20.2.2	initialize()	66
8.20.2.3	on_frame_end() override	66
8.20.3	Member Data Documentation	66
8.20.3.1	total_frames	66
8.21	FlimFrameInfo Class Reference	67
8.21.1	Member Function Documentation	67
8.21.1.1	getFrameNumber() const	67
8.21.1.2	getHistograms(std::function< uint32_t *(size_t, size_t)> array_out)	67
8.21.1.3	getIntensities(std::function< float *(size_t)> array_out)	67
8.21.1.4	getPixelBegins(std::function< long long *(size_t)> array_out)	67
8.21.1.5	getPixelEnds(std::function< long long *(size_t)> array_out)	67
8.21.1.6	getPixelPosition() const	67
8.21.1.7	getSummedCounts(std::function< uint64_t *(size_t)> array_out)	68
8.21.1.8	isValid() const	68
8.21.2	Friends And Related Function Documentation	68
8.21.2.1	Flim	68
8.21.3	Member Data Documentation	68
8.21.3.1	bins	68
8.21.3.2	frame_number	68
8.21.3.3	pixel_position	68
8.21.3.4	pixels	68
8.22	FrequencyMultiplier Class Reference	68
8.22.1	Detailed Description	69

8.22.2	Constructor & Destructor Documentation	69
8.22.2.1	FrequencyMultiplier(TimeTaggerBase *tagger, channel_t input_channel, int32_t multiplier)	69
8.22.2.2	~FrequencyMultiplier()	70
8.22.3	Member Function Documentation	70
8.22.3.1	getChannel()	70
8.22.3.2	getMultiplier()	70
8.22.3.3	next_impl(std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override	70
8.22.4	Friends And Related Function Documentation	70
8.22.4.1	FrequencyMultiplierImpl	70
8.23	GatedChannel Class Reference	70
8.23.1	Detailed Description	71
8.23.2	Constructor & Destructor Documentation	71
8.23.2.1	GatedChannel(TimeTaggerBase *tagger, channel_t input_channel, channel_t gate_start_channel, channel_t gate_stop_channel)	71
8.23.2.2	~GatedChannel()	72
8.23.3	Member Function Documentation	72
8.23.3.1	getChannel()	72
8.23.3.2	next_impl(std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override	72
8.23.4	Friends And Related Function Documentation	72
8.23.4.1	GatedChannelImpl	72
8.24	Histogram Class Reference	73
8.24.1	Detailed Description	74
8.24.2	Constructor & Destructor Documentation	74
8.24.2.1	Histogram(TimeTaggerBase *tagger, channel_t click_channel, channel_t start_channel=CHANNEL_UNUSED, timestamp_t binwidth=1000, int32_t n_bins=1000)	74
8.24.2.2	~Histogram()	74
8.24.3	Member Function Documentation	74
8.24.3.1	clear_impl() override	74
8.24.3.2	getData(std::function< int32_t *(size_t)> array_out)	75
8.24.3.3	getIndex(std::function< long long *(size_t)> array_out)	75

8.24.3.4	<a href="#">next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time) override</a>	75
8.24.3.5	<a href="#">on_start() override</a>	76
8.24.4	<a href="#">Friends And Related Function Documentation</a>	76
8.24.4.1	<a href="#">TimeDifferencesImpl&lt; Histogram &gt;</a>	76
8.25	<a href="#">Histogram2D Class Reference</a>	76
8.25.1	<a href="#">Detailed Description</a>	77
8.25.2	<a href="#">Constructor &amp; Destructor Documentation</a>	77
8.25.2.1	<a href="#">Histogram2D(TimeTaggerBase *tagger, channel_t start_channel, channel_t stop_channel_1, channel_t stop_channel_2, timestamp_t binwidth_1, timestamp_t binwidth_2, int32_t n_bins_1, int32_t n_bins_2)</a>	77
8.25.2.2	<a href="#">~Histogram2D()</a>	78
8.25.3	<a href="#">Member Function Documentation</a>	78
8.25.3.1	<a href="#">clear_impl() override</a>	78
8.25.3.2	<a href="#">getData(std::function&lt; int32_t *(size_t, size_t)&gt; array_out)</a>	78
8.25.3.3	<a href="#">getIndex(std::function&lt; long long *(size_t, size_t, size_t)&gt; array_out)</a>	78
8.25.3.4	<a href="#">getIndex_1(std::function&lt; long long *(size_t)&gt; array_out)</a>	78
8.25.3.5	<a href="#">getIndex_2(std::function&lt; long long *(size_t)&gt; array_out)</a>	78
8.25.3.6	<a href="#">next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time) override</a>	78
8.25.4	<a href="#">Friends And Related Function Documentation</a>	79
8.25.4.1	<a href="#">Histogram2DImpl</a>	79
8.26	<a href="#">HistogramLogBins Class Reference</a>	79
8.26.1	<a href="#">Detailed Description</a>	80
8.26.2	<a href="#">Constructor &amp; Destructor Documentation</a>	80
8.26.2.1	<a href="#">HistogramLogBins(TimeTaggerBase *tagger, channel_t click_channel, channel_t start_channel, double exp_start, double exp_stop, int32_t n_bins)</a>	80
8.26.2.2	<a href="#">~HistogramLogBins()</a>	80
8.26.3	<a href="#">Member Function Documentation</a>	80
8.26.3.1	<a href="#">clear_impl() override</a>	80
8.26.3.2	<a href="#">getBinEdges(std::function&lt; long long *(size_t)&gt; array_out)</a>	81
8.26.3.3	<a href="#">getData(std::function&lt; uint64_t *(size_t)&gt; array_out)</a>	81
8.26.3.4	<a href="#">getDataNormalizedCountsPerPs(std::function&lt; double *(size_t)&gt; array_out)</a>	81



8.26.3.5	<code>getDataNormalizedG2(std::function&lt; double *(size_t)&gt; array_out)</code>	81
8.26.3.6	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time) override</code>	81
8.26.4	Friends And Related Function Documentation	81
8.26.4.1	<code>HistogramLogBinsImpl</code>	82
8.27	Iterator Class Reference	82
8.27.1	Detailed Description	83
8.27.2	Constructor & Destructor Documentation	83
8.27.2.1	<code>Iterator(TimeTaggerBase *tagger, channel_t channel)</code>	83
8.27.2.2	<code>~Iterator()</code>	83
8.27.3	Member Function Documentation	83
8.27.3.1	<code>clear_impl()</code> override	83
8.27.3.2	<code>next()</code>	83
8.27.3.3	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time) override</code>	83
8.27.3.4	<code>size()</code>	84
8.27.4	Friends And Related Function Documentation	84
8.27.4.1	<code>IteratorImpl</code>	84
8.28	IteratorBase Class Reference	84
8.28.1	Detailed Description	87
8.28.2	Constructor & Destructor Documentation	87
8.28.2.1	<code>IteratorBase(TimeTaggerBase *tagger, std::string base_type_=""IteratorBase"", std::string extra_info_="")</code>	87
8.28.2.2	<code>~IteratorBase()</code>	87
8.28.3	Member Function Documentation	87
8.28.3.1	<code>clear()</code>	87
8.28.3.2	<code>clear_impl()</code>	87
8.28.3.3	<code>finish_running()</code>	87
8.28.3.4	<code>finishInitialization()</code>	87
8.28.3.5	<code>getCaptureDuration()</code>	88
8.28.3.6	<code>getLock()</code>	88
8.28.3.7	<code>getNewVirtualChannel()</code>	88

8.28.3.8	<code>isRunning()</code> . . . . .	88
8.28.3.9	<code>lock()</code> . . . . .	88
8.28.3.10	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time)=0</code> . . . . .	88
8.28.3.11	<code>on_start()</code> . . . . .	89
8.28.3.12	<code>on_stop()</code> . . . . .	89
8.28.3.13	<code>parallelize(OrderedPipeline &amp;pipeline)</code> . . . . .	89
8.28.3.14	<code>registerChannel(channel_t channel)</code> . . . . .	89
8.28.3.15	<code>start()</code> . . . . .	90
8.28.3.16	<code>startFor(timestamp_t capture_duration, bool clear=true)</code> . . . . .	90
8.28.3.17	<code>stop()</code> . . . . .	90
8.28.3.18	<code>unlock()</code> . . . . .	90
8.28.3.19	<code>unregisterChannel(channel_t channel)</code> . . . . .	90
8.28.3.20	<code>waitUntilFinished(int64_t timeout=-1)</code> . . . . .	91
8.28.4	Friends And Related Function Documentation . . . . .	91
8.28.4.1	<code>SynchronizedMeasurements</code> . . . . .	91
8.28.4.2	<code>TimeTaggerProxy</code> . . . . .	91
8.28.4.3	<code>TimeTaggerRunner</code> . . . . .	91
8.28.5	Member Data Documentation . . . . .	91
8.28.5.1	<code>autostart</code> . . . . .	91
8.28.5.2	<code>capture_duration</code> . . . . .	91
8.28.5.3	<code>channels_registered</code> . . . . .	91
8.28.5.4	<code>running</code> . . . . .	91
8.28.5.5	<code>tagger</code> . . . . .	91
8.29	OrderedBarrier Class Reference . . . . .	91
8.29.1	Constructor & Destructor Documentation . . . . .	92
8.29.1.1	<code>OrderedBarrier()</code> . . . . .	92
8.29.1.2	<code>~OrderedBarrier()</code> . . . . .	92
8.29.2	Member Function Documentation . . . . .	92
8.29.2.1	<code>queue()</code> . . . . .	92
8.29.2.2	<code>waitUntilFinished()</code> . . . . .	92

8.29.3 Friends And Related Function Documentation . . . . .	92
8.29.3.1 OrderInstance . . . . .	92
8.30 OrderedPipeline Class Reference . . . . .	92
8.30.1 Constructor & Destructor Documentation . . . . .	93
8.30.1.1 OrderedPipeline() . . . . .	93
8.30.1.2 ~OrderedPipeline() . . . . .	93
8.30.2 Friends And Related Function Documentation . . . . .	93
8.30.2.1 IteratorBase . . . . .	93
8.31 OrderedBarrier::OrderInstance Class Reference . . . . .	93
8.31.1 Constructor & Destructor Documentation . . . . .	93
8.31.1.1 OrderInstance() . . . . .	93
8.31.1.2 OrderInstance(OrderedBarrier *parent, uint64_t instance_id) . . . . .	93
8.31.1.3 ~OrderInstance() . . . . .	93
8.31.2 Member Function Documentation . . . . .	93
8.31.2.1 release() . . . . .	93
8.31.2.2 sync() . . . . .	93
8.31.3 Friends And Related Function Documentation . . . . .	93
8.31.3.1 OrderedBarrier . . . . .	93
8.32 Scope Class Reference . . . . .	94
8.32.1 Constructor & Destructor Documentation . . . . .	94
8.32.1.1 Scope(TimeTaggerBase *tagger, std::vector< channel_t > event_channels, channel_t trigger_channel, timestamp_t window_size=1000000000, int32_t n_↵ traces=1, int32_t n_max_events=1000) . . . . .	94
8.32.1.2 ~Scope() . . . . .	95
8.32.2 Member Function Documentation . . . . .	95
8.32.2.1 clear_impl() override . . . . .	95
8.32.2.2 getData() . . . . .	95
8.32.2.3 getWindowSize() . . . . .	95
8.32.2.4 next_impl(std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override . . . . .	95
8.32.2.5 ready() . . . . .	96
8.32.2.6 triggered() . . . . .	96

8.32.3	Friends And Related Function Documentation	96
8.32.3.1	ScopImpl	96
8.33	StartStop Class Reference	96
8.33.1	Detailed Description	97
8.33.2	Constructor & Destructor Documentation	97
8.33.2.1	StartStop(TimeTaggerBase *tagger, channel_t click_channel, channel_t start_↵ channel=CHANNEL_UNUSED, timestamp_t binwidth=1000)	97
8.33.2.2	~StartStop()	97
8.33.3	Member Function Documentation	97
8.33.3.1	clear_impl() override	97
8.33.3.2	getData(std::function< long long *(size_t, size_t)> array_out)	98
8.33.3.3	next_impl(std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time) override	98
8.33.3.4	on_start() override	98
8.33.4	Friends And Related Function Documentation	98
8.33.4.1	StartStopImpl	98
8.34	SynchronizedMeasurements Class Reference	98
8.34.1	Detailed Description	99
8.34.2	Constructor & Destructor Documentation	99
8.34.2.1	SynchronizedMeasurements(TimeTaggerBase *tagger)	99
8.34.2.2	~SynchronizedMeasurements()	100
8.34.3	Member Function Documentation	100
8.34.3.1	clear()	100
8.34.3.2	getTagger()	100
8.34.3.3	isRunning()	100
8.34.3.4	registerMeasurement(IteratorBase *measurement)	100
8.34.3.5	runCallback(TimeTaggerBase::IteratorCallback callback, std::unique_lock< std↵ ::mutex > &lk, bool block=true)	100
8.34.3.6	start()	100
8.34.3.7	startFor(timestamp_t capture_duration, bool clear=true)	100
8.34.3.8	stop()	101
8.34.3.9	unregisterMeasurement(IteratorBase *measurement)	101

8.34.3.10	<code>waitUntilFinished(int64_t timeout=-1)</code>	101
8.34.4	Friends And Related Function Documentation	101
8.34.4.1	<code>TimeTaggerProxy</code>	101
8.35	<code>SyntheticSingleTag</code> Class Reference	101
8.35.1	Detailed Description	102
8.35.2	Constructor & Destructor Documentation	102
8.35.2.1	<code>SyntheticSingleTag(TimeTaggerBase *tagger, channel_t base_channel=CHANNEL_UNUSED)</code>	102
8.35.2.2	<code>~SyntheticSingleTag()</code>	102
8.35.3	Member Function Documentation	102
8.35.3.1	<code>getChannel() const</code>	103
8.35.3.2	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time) override</code>	103
8.35.3.3	<code>trigger()</code>	103
8.35.4	Friends And Related Function Documentation	103
8.35.4.1	<code>SyntheticSingleTagImpl</code>	103
8.36	<code>Tag</code> Struct Reference	103
8.36.1	Detailed Description	104
8.36.2	Member Enumeration Documentation	104
8.36.2.1	Type	104
8.36.3	Member Data Documentation	104
8.36.3.1	<code>channel</code>	104
8.36.3.2	<code>missed_events</code>	104
8.36.3.3	<code>reserved</code>	104
8.36.3.4	<code>time</code>	104
8.36.3.5	<code>type</code>	104
8.37	<code>TimeDifferences</code> Class Reference	105
8.37.1	Detailed Description	106
8.37.2	Constructor & Destructor Documentation	106
8.37.2.1	<code>TimeDifferences(TimeTaggerBase *tagger, channel_t click_channel, channel_t start_channel=CHANNEL_UNUSED, channel_t next_channel=CHANNEL_UNUSED, channel_t sync_channel=CHANNEL_UNUSED, timestamp_t binwidth=1000, int32_t n_bins=1000, int32_t n_histograms=1)</code>	106

8.37.2.2	<code>~TimeDifferences()</code>	107
8.37.3	Member Function Documentation	107
8.37.3.1	<code>clear_impl()</code> override	107
8.37.3.2	<code>getCounts()</code>	107
8.37.3.3	<code>getData(std::function&lt; int32_t *(size_t, size_t)&gt; array_out)</code>	107
8.37.3.4	<code>getIndex(std::function&lt; long long *(size_t)&gt; array_out)</code>	107
8.37.3.5	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time)</code> override	107
8.37.3.6	<code>on_start()</code> override	108
8.37.3.7	<code>ready()</code>	108
8.37.3.8	<code>setMaxCounts(uint64_t max_counts)</code>	108
8.37.4	Friends And Related Function Documentation	108
8.37.4.1	<code>TimeDifferencesImpl&lt; TimeDifferences &gt;</code>	108
8.38	<code>TimeDifferencesImpl&lt; T &gt;</code> Class Template Reference	108
8.39	<code>TimeDifferencesND</code> Class Reference	109
8.39.1	Detailed Description	110
8.39.2	Constructor & Destructor Documentation	110
8.39.2.1	<code>TimeDifferencesND(TimeTaggerBase *tagger, channel_t click_channel, channel_t start_channel, std::vector&lt; channel_t &gt; next_channels, std::vector&lt; channel_t &gt; sync_channels, std::vector&lt; int32_t &gt; n_histograms, timestamp_t binwidth, int32_t n_bins)</code>	110
8.39.2.2	<code>~TimeDifferencesND()</code>	110
8.39.3	Member Function Documentation	110
8.39.3.1	<code>clear_impl()</code> override	111
8.39.3.2	<code>getData(std::function&lt; int32_t *(size_t, size_t)&gt; array_out)</code>	111
8.39.3.3	<code>getIndex(std::function&lt; long long *(size_t)&gt; array_out)</code>	111
8.39.3.4	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time)</code> override	111
8.39.3.5	<code>on_start()</code> override	111
8.39.4	Friends And Related Function Documentation	112
8.39.4.1	<code>TimeDifferencesNDImpl</code>	112
8.40	<code>TimeTagger</code> Class Reference	112
8.40.1	Detailed Description	114

8.40.2 Member Function Documentation . . . . .	114
8.40.2.1 autoCalibration(std::function< double *(size_t)> array_out)=0 . . . . .	114
8.40.2.2 clearConditionalFilter()=0 . . . . .	114
8.40.2.3 factoryAccess(uint32_t pw, uint32_t addr, uint32_t data, uint32_t mask)=0 . . . . .	114
8.40.2.4 getChannelList(ChannelEdge type=ChannelEdge::All)=0 . . . . .	114
8.40.2.5 getChannelNumberScheme()=0 . . . . .	115
8.40.2.6 getConditionalFilterFiltered()=0 . . . . .	115
8.40.2.7 getConditionalFilterTrigger()=0 . . . . .	115
8.40.2.8 getDACRange()=0 . . . . .	115
8.40.2.9 getDistributionCount(std::function< uint64_t *(size_t, size_t)> array_out)=0 . . . . .	115
8.40.2.10 getDistributionPSecs(std::function< long long *(size_t, size_t)> array_out)=0 . . . . .	115
8.40.2.11 getEventDivider(channel_t channel)=0 . . . . .	115
8.40.2.12 getFirmwareVersion()=0 . . . . .	116
8.40.2.13 getHardwareBufferSize()=0 . . . . .	116
8.40.2.14 getHardwareDelayCompensation(channel_t channel)=0 . . . . .	116
8.40.2.15 getInputMux(channel_t channel)=0 . . . . .	116
8.40.2.16 getLicenseInfo()=0 . . . . .	117
8.40.2.17 getModel()=0 . . . . .	117
8.40.2.18 getNormalization(channel_t channel)=0 . . . . .	117
8.40.2.19 getPcbVersion()=0 . . . . .	117
8.40.2.20 getPsPerClock()=0 . . . . .	117
8.40.2.21 getSensorData()=0 . . . . .	118
8.40.2.22 getSerial()=0 . . . . .	118
8.40.2.23 getStreamBlockSizeEvents()=0 . . . . .	118
8.40.2.24 getStreamBlockSizeLatency()=0 . . . . .	118
8.40.2.25 getTestSignalDivider()=0 . . . . .	118
8.40.2.26 getTriggerLevel(channel_t channel)=0 . . . . .	118
8.40.2.27 reset()=0 . . . . .	118
8.40.2.28 setConditionalFilter(std::vector< channel_t > trigger, std::vector< channel_t > filtered, bool hardwareDelayCompensation=true)=0 . . . . .	118
8.40.2.29 setEventDivider(channel_t channel, unsigned int divider)=0 . . . . .	119

8.40.2.30	setHardwareBufferSize(int size)=0 . . . . .	119
8.40.2.31	setInputMux(channel_t channel, int mux_mode)=0 . . . . .	119
8.40.2.32	setLED(uint32_t bitmask)=0 . . . . .	119
8.40.2.33	setNormalization(std::vector< channel_t > channel, bool state)=0 . . . . .	120
8.40.2.34	setSoundFrequency(uint32_t freq_hz)=0 . . . . .	120
8.40.2.35	setStreamBlockSize(int max_events, int max_latency)=0 . . . . .	120
8.40.2.36	setTestSignalDivider(int divider)=0 . . . . .	120
8.40.2.37	setTriggerLevel(channel_t channel, double voltage)=0 . . . . .	121
8.41	TimeTaggerBase Class Reference . . . . .	121
8.41.1	Member Typedef Documentation . . . . .	123
8.41.1.1	IteratorCallback . . . . .	123
8.41.1.2	IteratorCallbackMap . . . . .	123
8.41.2	Constructor & Destructor Documentation . . . . .	123
8.41.2.1	TimeTaggerBase() . . . . .	123
8.41.2.2	~TimeTaggerBase() . . . . .	123
8.41.2.3	TimeTaggerBase(const TimeTaggerBase &)=delete . . . . .	123
8.41.3	Member Function Documentation . . . . .	123
8.41.3.1	addChild(TimeTaggerBase *child)=0 . . . . .	123
8.41.3.2	addIterator(IteratorBase *it)=0 . . . . .	123
8.41.3.3	clearOverflows()=0 . . . . .	123
8.41.3.4	freeIterator(IteratorBase *it)=0 . . . . .	123
8.41.3.5	freeVirtualChannel(channel_t channel)=0 . . . . .	123
8.41.3.6	getConfiguration()=0 . . . . .	123
8.41.3.7	getDeadtime(channel_t channel)=0 . . . . .	123
8.41.3.8	getDelayHardware(channel_t channel)=0 . . . . .	124
8.41.3.9	getDelaySoftware(channel_t channel)=0 . . . . .	124
8.41.3.10	getFence(bool alloc_fence=true)=0 . . . . .	124
8.41.3.11	getInputDelay(channel_t channel)=0 . . . . .	125
8.41.3.12	getInvertedChannel(channel_t channel)=0 . . . . .	125
8.41.3.13	getNewVirtualChannel()=0 . . . . .	125



8.41.3.14	<code>getOverflows()</code>	125
8.41.3.15	<code>getOverflowsAndClear()</code>	125
8.41.3.16	<code>getTestSignal(channel_t channel)</code>	125
8.41.3.17	<code>isUnusedChannel(channel_t channel)</code>	125
8.41.3.18	<code>operator=(const TimeTaggerBase &amp;)=delete</code>	126
8.41.3.19	<code>registerChannel(channel_t channel)</code>	126
8.41.3.20	<code>release()</code>	126
8.41.3.21	<code>removeChild(TimeTaggerBase *child)</code>	126
8.41.3.22	<code>runSynchronized(const IteratorCallbackMap &amp;callbacks, bool block=true)</code>	126
8.41.3.23	<code>setDeadtime(channel_t channel, timestamp_t deadtime)</code>	126
8.41.3.24	<code>setDelayHardware(channel_t channel, timestamp_t delay)</code>	127
8.41.3.25	<code>setDelaySoftware(channel_t channel, timestamp_t delay)</code>	127
8.41.3.26	<code>setInputDelay(channel_t channel, timestamp_t delay)</code>	127
8.41.3.27	<code>setTestSignal(channel_t channel, bool enabled)</code>	128
8.41.3.28	<code>setTestSignal(std::vector&lt; channel_t &gt; channel, bool enabled)</code>	128
8.41.3.29	<code>sync(int64_t timeout=-1)</code>	128
8.41.3.30	<code>unregisterChannel(channel_t channel)</code>	128
8.41.3.31	<code>waitForFence(unsigned int fence, int64_t timeout=-1)</code>	128
8.41.4	Friends And Related Function Documentation	129
8.41.4.1	IteratorBase	129
8.41.4.2	TimeTaggerProxy	129
8.41.4.3	TimeTaggerRunner	129
8.42	TimeTaggerVirtual Class Reference	129
8.42.1	Detailed Description	130
8.42.2	Member Function Documentation	130
8.42.2.1	<code>clearConditionalFilter()</code>	130
8.42.2.2	<code>getConditionalFilterFiltered()</code>	130
8.42.2.3	<code>getConditionalFilterTrigger()</code>	131
8.42.2.4	<code>getReplaySpeed()</code>	131
8.42.2.5	<code>replay(const std::string &amp;file, timestamp_t begin=0, timestamp_t duration=-1, bool queue=true)</code>	131

8.42.2.6	<code>reset()</code> =0	131
8.42.2.7	<code>setConditionalFilter(std::vector&lt; channel_t &gt; trigger, std::vector&lt; channel_t &gt; filtered)</code> =0	131
8.42.2.8	<code>setReplaySpeed(double speed)</code> =0	132
8.42.2.9	<code>stop()</code> =0	132
8.42.2.10	<code>waitForCompletion(uint64_t ID=0, int64_t timeout=-1)</code> =0	132
8.43	TimeTagStream Class Reference	133
8.43.1	Detailed Description	133
8.43.2	Constructor & Destructor Documentation	134
8.43.2.1	<code>TimeTagStream(TimeTaggerBase *tagger, uint64_t n_max_events, std::vector&lt; channel_t &gt; channels=std::vector&lt; channel_t &gt;())</code>	134
8.43.2.2	<code>~TimeTagStream()</code>	134
8.43.3	Member Function Documentation	134
8.43.3.1	<code>clear_impl()</code> override	134
8.43.3.2	<code>getCounts()</code>	134
8.43.3.3	<code>getData()</code>	134
8.43.3.4	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time)</code> override	134
8.43.4	Friends And Related Function Documentation	135
8.43.4.1	<code>TimeTagStreamImpl</code>	135
8.44	TimeTagStreamBuffer Class Reference	135
8.44.1	Member Function Documentation	136
8.44.1.1	<code>getChannels(std::function&lt; int *(size_t)&gt; array_out)</code>	136
8.44.1.2	<code>getEventTypes(std::function&lt; unsigned char *(size_t)&gt; array_out)</code>	136
8.44.1.3	<code>getMissedEvents(std::function&lt; unsigned short *(size_t)&gt; array_out)</code>	136
8.44.1.4	<code>getOverflows(std::function&lt; unsigned char *(size_t)&gt; array_out)</code>	136
8.44.1.5	<code>getTimestamps(std::function&lt; long long *(size_t)&gt; array_out)</code>	136
8.44.2	Friends And Related Function Documentation	136
8.44.2.1	<code>FileReaderImpl</code>	136
8.44.2.2	<code>TimeTagStreamImpl</code>	136
8.44.3	Member Data Documentation	136
8.44.3.1	<code>hasOverflows</code>	136

8.44.3.2	<code>size</code>	136
8.44.3.3	<code>tGetData</code>	136
8.44.3.4	<code>tStart</code>	136
8.45	TriggerOnCountrate Class Reference	136
8.45.1	Detailed Description	137
8.45.2	Constructor & Destructor Documentation	138
8.45.2.1	<code>TriggerOnCountrate(TimeTaggerBase *tagger, channel_t input_channel, double reference_countrate, double hysteresis, timestamp_t time_window)</code>	138
8.45.2.2	<code>~TriggerOnCountrate()</code>	138
8.45.3	Member Function Documentation	138
8.45.3.1	<code>getChannelAbove()</code>	138
8.45.3.2	<code>getChannelBelow()</code>	138
8.45.3.3	<code>getChannels()</code>	138
8.45.3.4	<code>getCurrentCountrate()</code>	138
8.45.3.5	<code>injectCurrentState()</code>	138
8.45.3.6	<code>isAbove()</code>	139
8.45.3.7	<code>isBelow()</code>	139
8.45.3.8	<code>next_impl(std::vector&lt; Tag &gt; &amp;incoming_tags, timestamp_t begin_time, timestamp_t end_time) override</code>	139
8.45.3.9	<code>on_start() override</code>	139
8.45.4	Friends And Related Function Documentation	139
8.45.4.1	<code>TriggerOnCountrateImpl</code>	139

<b>9 File Documentation</b>	<b>141</b>
9.1 Iterators.h File Reference	141
9.1.1 Macro Definition Documentation	143
9.1.1.1 BINNING_TEMPLATE_HELPER	143
9.1.2 Enumeration Type Documentation	144
9.1.2.1 CoincidenceTimestamp	144
9.1.2.2 State	144
9.2 TimeTagger.h File Reference	144
9.2.1 Macro Definition Documentation	147
9.2.1.1 channel_t	147
9.2.1.2 ErrorLog	147
9.2.1.3 ErrorLogCensored	147
9.2.1.4 GET_DATA_1D	147
9.2.1.5 GET_DATA_1D_OP1	148
9.2.1.6 GET_DATA_1D_OP2	148
9.2.1.7 GET_DATA_2D	148
9.2.1.8 GET_DATA_2D_OP1	148
9.2.1.9 GET_DATA_2D_OP2	148
9.2.1.10 GET_DATA_3D	148
9.2.1.11 InfoLog	148
9.2.1.12 InfoLogCensored	148
9.2.1.13 LogMessage	148
9.2.1.14 LogMessageCensored	148
9.2.1.15 timestamp_t	148
9.2.1.16 TIMETAGGER_VERSION	149
9.2.1.17 TT_API	149
9.2.1.18 WarningLog	149
9.2.1.19 WarningLogCensored	149
9.2.2 Typedef Documentation	149
9.2.2.1 _Iterator	149

9.2.2.2	<code>logger_callback</code>	149
9.2.3	Enumeration Type Documentation	149
9.2.3.1	<code>ChannelEdge</code>	149
9.2.3.2	<code>FrontendType</code>	149
9.2.3.3	<code>LanguageUsed</code>	150
9.2.3.4	<code>LogLevel</code>	150
9.2.3.5	<code>Resolution</code>	150
9.2.3.6	<code>UsageStatisticsStatus</code>	150
9.2.4	Function Documentation	150
9.2.4.1	<code>createTimeTagger(std::string serial="", Resolution resolution=Resolution::↔ Standard)</code>	150
9.2.4.2	<code>createTimeTaggerVirtual()</code>	151
9.2.4.3	<code>extractLicenseInfo(const std::string &amp;license)</code>	151
9.2.4.4	<code>flashLicense(const std::string &amp;serial, const std::string &amp;license)</code>	151
9.2.4.5	<code>freeTimeTagger(TimeTaggerBase *tagger)</code>	151
9.2.4.6	<code>getTimeTaggerChannelNumberScheme()</code>	152
9.2.4.7	<code>getTimeTaggerModel(const std::string &amp;serial)</code>	152
9.2.4.8	<code>getUsageStatisticsReport()</code>	152
9.2.4.9	<code>getUsageStatisticsStatus()</code>	152
9.2.4.10	<code>getVersion()</code>	152
9.2.4.11	<code>hasTimeTaggerVirtualLicense()</code>	152
9.2.4.12	<code>LogBase(LogLevel level, const char *file, int line, bool censored, const char *fmt,...)</code>	152
9.2.4.13	<code>scanTimeTagger()</code>	152
9.2.4.14	<code>setCustomBitFileName(const std::string &amp;bitFileName)</code>	152
9.2.4.15	<code>setFrontend(FrontendType frontend)</code>	153
9.2.4.16	<code>setLanguageInfo(std::uint32_t pw, LanguageUsed language, std::string version)</code>	153
9.2.4.17	<code>setLogger(logger_callback callback)</code>	153
9.2.4.18	<code>setTimeTaggerChannelNumberScheme(int scheme)</code>	153
9.2.4.19	<code>setUsageStatisticsStatus(UsageStatisticsStatus new_status)</code>	154
9.2.5	Variable Documentation	154
9.2.5.1	<code>CHANNEL_UNUSED</code>	154
9.2.5.2	<code>CHANNEL_UNUSED_OLD</code>	154
9.2.5.3	<code>TT_CHANNEL_FALLING_EDGES</code>	154
9.2.5.4	<code>TT_CHANNEL_NUMBER_SCHEME_AUTO</code>	154
9.2.5.5	<code>TT_CHANNEL_NUMBER_SCHEME_ONE</code>	154
9.2.5.6	<code>TT_CHANNEL_NUMBER_SCHEME_ZERO</code>	154
9.2.5.7	<code>TT_CHANNEL_RISING_AND_FALLING_EDGES</code>	154
9.2.5.8	<code>TT_CHANNEL_RISING_EDGES</code>	154



# Chapter 1

## TimeTagger

backend for [TimeTagger](#), an OpalKelly based single photon counting library

### Author

Markus Wick [markus@swabianinstruments.com](mailto:markus@swabianinstruments.com)

Helmut Fedder [helmut@swabianinstruments.com](mailto:helmut@swabianinstruments.com)

Michael Schlagmüller [michael@swabianinstruments.com](mailto:michael@swabianinstruments.com)

[TimeTagger](#) provides an easy to use and cost effective hardware solution for time-resolved single photon counting applications.

This document describes the C++ native interface to the [TimeTagger](#) device.





## Chapter 2

# Deprecated List

### Class **Dump**

use [FileWriter](#)

### Class **Iterator**

use [TimeTagStream](#)

### Member **IteratorBase::lock** ()

use `getLock`

### Member **IteratorBase::unlock** ()

use `getLock`

### Member **TimeTagger::getDistributionPSecs** (std::function< long long \*(size\_t, size\_t)> array\_out)=0



## Chapter 3

# Module Index

### 3.1 Modules

Here is a list of all modules:

base iterators . . . . .	13
--------------------------	----



## Chapter 4

# Hierarchical Index

### 4.1 Class Hierarchy

This inheritance list is sorted roughly, but not completely, alphabetically:

CustomLogger . . . . .	35
Event . . . . .	44
FastBinning . . . . .	47
FileReader . . . . .	48
FlimFrameInfo . . . . .	67
IteratorBase . . . . .	84
Coincidences . . . . .	16
Coincidence . . . . .	15
Combiner . . . . .	19
ConstantFractionDiscriminator . . . . .	21
Correlation . . . . .	24
CountBetweenMarkers . . . . .	27
Counter . . . . .	29
Countrate . . . . .	32
CustomMeasurementBase . . . . .	36
DelayedChannel . . . . .	38
Dump . . . . .	41
EventGenerator . . . . .	44
FileWriter . . . . .	50
FlimAbstract . . . . .	60
Flim . . . . .	54
FlimBase . . . . .	65
FrequencyMultiplier . . . . .	68
GatedChannel . . . . .	70
Histogram . . . . .	73
Histogram2D . . . . .	76
HistogramLogBins . . . . .	79
Iterator . . . . .	82
Scope . . . . .	94
StartStop . . . . .	96
SyntheticSingleTag . . . . .	101
TimeDifferences . . . . .	105
TimeDifferencesND . . . . .	109
TimeTagStream . . . . .	133
TriggerOnCountrate . . . . .	136

OrderedBarrier . . . . .	91
OrderedPipeline . . . . .	92
OrderedBarrier::OrderInstance . . . . .	93
SynchronizedMeasurements . . . . .	98
Tag . . . . .	103
TimeDifferencesImpl< T > . . . . .	108
TimeTaggerBase . . . . .	121
TimeTagger . . . . .	112
TimeTaggerVirtual . . . . .	129
TimeTagStreamBuffer . . . . .	135

## Chapter 5

# Class Index

### 5.1 Class List

Here are the classes, structs, unions and interfaces with brief descriptions:

Coincidence	
Coincidence	monitor for one or more channel groups . . . . . 15
Coincidences	
Coincidence	monitor for one or more channel groups . . . . . 16
Combiner	
Combine	some channels in a virtual channel which has a tick for each tick in the input channels 19
ConstantFractionDiscriminator	
Virtual CFD	implementation which returns the mean time between a raising and a falling pair of edges . . . . . 21
Correlation	
Cross-correlation	between two channels . . . . . 24
CountBetweenMarkers	
Simple counter	where external marker signals determine the bins . . . . . 27
Counter	
Simple counter	on one or more channels . . . . . 29
Countrate	
Count rate	on one or more channels . . . . . 32
CustomLogger	. . . . . 35
CustomMeasurementBase	. . . . . 36
DelayedChannel	
Simple delayed queue	. . . . . 38
Dump	
Dump	all time tags to a file . . . . . 41
Event	. . . . . 44
EventGenerator	
Generate predefined events	in a virtual channel relative to a trigger event . . . . . 44
FastBinning	. . . . . 47
FileReader	. . . . . 48
FileWriter	
Compresses and stores	all time tags to a file . . . . . 50
Flim	
Fluorescence lifetime imaging	. . . . . 54
FlimAbstract	. . . . . 60
FlimBase	. . . . . 65
FlimFrameInfo	. . . . . 67

<a href="#">FrequencyMultiplier</a>	
The signal of an input channel is scaled up to a higher frequency according to the multiplier passed as a parameter . . . . .	68
<a href="#">GatedChannel</a>	
An input channel is gated by a gate channel . . . . .	70
<a href="#">Histogram</a>	
Accumulate time differences into a histogram . . . . .	73
<a href="#">Histogram2D</a>	
A 2-dimensional histogram of time differences. This can be used in measurements similar to 2D NRM spectroscopy . . . . .	76
<a href="#">HistogramLogBins</a>	
Accumulate time differences into a histogram with logarithmic increasing bin sizes . . . . .	79
<a href="#">Iterator</a>	
Simple event queue . . . . .	82
<a href="#">IteratorBase</a>	
Base class for all iterators . . . . .	84
<a href="#">OrderedBarrier</a>	91
<a href="#">OrderedPipeline</a>	92
<a href="#">OrderedBarrier::OrderInstance</a>	93
<a href="#">Scope</a>	94
<a href="#">StartStop</a>	
Simple start-stop measurement . . . . .	96
<a href="#">SynchronizedMeasurements</a>	
Start, stop and clear several measurements synchronized . . . . .	98
<a href="#">SyntheticSingleTag</a>	
Synthetic trigger timetag generator . . . . .	101
<a href="#">Tag</a>	
Single event on a channel . . . . .	103
<a href="#">TimeDifferences</a>	
Accumulates the time differences between clicks on two channels in one or more histograms .	105
<a href="#">TimeDifferencesImpl&lt; T &gt;</a>	108
<a href="#">TimeDifferencesND</a>	
Accumulates the time differences between clicks on two channels in a multi-dimensional histogram . . . . .	109
<a href="#">TimeTagger</a>	
Backend for the <a href="#">TimeTagger</a> . . . . .	112
<a href="#">TimeTaggerBase</a>	121
<a href="#">TimeTaggerVirtual</a>	
Virtual <a href="#">TimeTagger</a> based on dump files . . . . .	129
<a href="#">TimeTagStream</a>	
Access the time tag stream . . . . .	133
<a href="#">TimeTagStreamBuffer</a>	135
<a href="#">TriggerOnCountrate</a>	
Inject trigger events when exceeding or falling below a given count rate within a rolling time window . . . . .	136



## Chapter 6

# File Index

### 6.1 File List

Here is a list of all files with brief descriptions:

<a href="#">Iterators.h</a> . . . . .	141
<a href="#">TimeTagger.h</a> . . . . .	144



# Chapter 7

## Module Documentation

### 7.1 base iterators

base iterators for photon counting applications

#### Classes

- class [Combiner](#)  
*Combine some channels in a virtual channel which has a tick for each tick in the input channels.*
- class [CountBetweenMarkers](#)  
*a simple counter where external marker signals determine the bins*
- class [Counter](#)  
*a simple counter on one or more channels*
- class [Coincidences](#)  
*a coincidence monitor for one or more channel groups*
- class [Coincidence](#)  
*a coincidence monitor for one or more channel groups*
- class [Countrate](#)  
*count rate on one or more channels*
- class [DelayedChannel](#)  
*a simple delayed queue*
- class [TriggerOnCountrate](#)  
*Inject trigger events when exceeding or falling below a given count rate within a rolling time window.*
- class [GatedChannel](#)  
*An input channel is gated by a gate channel.*
- class [FrequencyMultiplier](#)  
*The signal of an input channel is scaled up to a higher frequency according to the multiplier passed as a parameter.*
- class [Iterator](#)  
*a simple event queue*
- class [TimeTagStream](#)  
*access the time tag stream*
- class [Dump](#)  
*dump all time tags to a file*
- class [StartStop](#)  
*simple start-stop measurement*

- class [TimeDifferences](#)  
*Accumulates the time differences between clicks on two channels in one or more histograms.*
- class [Histogram2D](#)  
*A 2-dimensional histogram of time differences. This can be used in measurements similar to 2D NRM spectroscopy.*
- class [TimeDifferencesND](#)  
*Accumulates the time differences between clicks on two channels in a multi-dimensional histogram.*
- class [Histogram](#)  
*Accumulate time differences into a histogram.*
- class [HistogramLogBins](#)  
*Accumulate time differences into a histogram with logarithmic increasing bin sizes.*
- class [Correlation](#)  
*cross-correlation between two channels*
- class [Scope](#)
- class [SynchronizedMeasurements](#)  
*start, stop and clear several measurements synchronized*
- class [ConstantFractionDiscriminator](#)  
*a virtual CFD implementation which returns the mean time between a raising and a falling pair of edges*
- class [FileWriter](#)  
*compresses and stores all time tags to a file*
- class [EventGenerator](#)  
*Generate predefined events in a virtual channel relative to a trigger event.*
- class [Flim](#)  
*Fluorescence lifetime imaging.*
- class [SyntheticSingleTag](#)  
*synthetic trigger timetag generator.*

### 7.1.1 Detailed Description

base iterators for photon counting applications

## Chapter 8


# Class Documentation

### 8.1 Coincidence Class Reference

a coincidence monitor for one or more channel groups

```
#include <Iterators.h>
```

Inheritance diagram for Coincidence:



#### Public Member Functions

- `Coincidence (TimeTaggerBase *tagger, std::vector< channel_t > channels, timestamp_t coincidence←Window=1000, CoincidenceTimestamp timestamp=CoincidenceTimestamp::Last)`  
*construct a coincidence*
- `channel_t getChannel ()`  
*virtual channel which contains the coincidences*

## Additional Inherited Members

### 8.1.1 Detailed Description

a coincidence monitor for one or more channel groups

Monitor coincidences for a given channel groups passed by the constructor. A coincidence is event is detected when all slected channels have a click within the given coincidenceWindow [ps] The coincidence will create a virtual events on a virtual channel with the channel number provided by [getChannel\(\)](#). For multiple coincidence channel combinations use the class [Coincidences](#) which outperforms multiple instances of Conincidence.

### 8.1.2 Constructor & Destructor Documentation

**8.1.2.1** `Coincidence::Coincidence ( TimeTaggerBase * tagger, std::vector< channel_t > channels, timestamp_t coincidenceWindow = 1000, CoincidenceTimestamp timestamp = CoincidenceTimestamp::Last )`  
`[inline]`

construct a coincidence

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>channels</i>	vector of channels to match
<i>coincidenceWindow</i>	max distance between all clicks for a coincidence [ps]
<i>timestamp</i>	type of timestamp for virtual channel (Last, Average, First, ListedFirst)

### 8.1.3 Member Function Documentation

**8.1.3.1** `channel_t Coincidence::getChannel ( )` `[inline]`

virtual channel which contains the coincidences

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.2 Coincidences Class Reference

a coincidence monitor for one or more channel groups

```
#include <Iterators.h>
```

Inheritance diagram for Coincidences:



### Public Member Functions

- `Coincidences` (`TimeTaggerBase` \*tagger, `std::vector`< `std::vector`< `channel_t` >> coincidenceGroups, `timestamp_t` coincidenceWindow, `CoincidenceTimestamp` timestamp=`CoincidenceTimestamp::Last`)  
construct a `Coincidences`
- `~Coincidences` ()
- `std::vector`< `channel_t` > `getChannels` ()  
fetches the block of virtual channels for those coincidence groups
- void `setCoincidenceWindow` (`timestamp_t` coincidenceWindow)

### Protected Member Functions

- bool `next_impl` (`std::vector`< `Tag` > &incoming\_tags, `timestamp_t` begin\_time, `timestamp_t` end\_time) over-ride  
update iterator state

### Friends

- class `CoincidencesImpl`

### Additional Inherited Members

#### 8.2.1 Detailed Description

a coincidence monitor for one or more channel groups

Monitor coincidences for given coincidence groups passed by the constructor. A coincidence is hereby defined as for a given coincidence group a) the incoming is part of this group b) at least tag arrived within the coincidence↔ Window [ps] for all other channels of this coincidence group Each coincidence will create a virtual event. The block of event IDs for those coincidence group can be fetched.

## 8.2.2 Constructor & Destructor Documentation

**8.2.2.1** `Coincidences::Coincidences ( TimeTaggerBase * tagger, std::vector< std::vector< channel_t >> coincidenceGroups, timestamp_t coincidenceWindow, CoincidenceTimestamp timestamp = CoincidenceTimestamp::Last )`

construct a [Coincidences](#)

### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>coincidenceGroups</i>	a vector of channels defining the coincidences
<i>coincidenceWindow</i>	the size of the coincidence window in picoseconds
<i>timestamp</i>	type of timestamp for virtual channel (Last, Average, First, ListedFirst)

**8.2.2.2** `Coincidences::~~Coincidences ( )`

## 8.2.3 Member Function Documentation

**8.2.3.1** `std::vector<channel_t> Coincidences::getChannels ( )`

fetches the block of virtual channels for those coincidence groups

**8.2.3.2** `bool Coincidences::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time )` `[override]`, `[protected]`, `[virtual]`

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

### Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

### Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

**8.2.3.3** `void Coincidences::setCoincidenceWindow ( timestamp_t coincidenceWindow )`



## 8.2.4 Friends And Related Function Documentation

### 8.2.4.1 friend class CoincidencesImpl [friend]

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.3 Combiner Class Reference

Combine some channels in a virtual channel which has a tick for each tick in the input channels.

```
#include <Iterators.h>
```

Inheritance diagram for Combiner:



### Public Member Functions

- `Combiner` (`TimeTaggerBase` \*`tagger`, `std::vector`< `channel_t` > `channels`)  
*construct a combiner*
- `~Combiner` ()
- `void` `getData` (`std::function`< `int64_t` \*(`size_t`)> `array_out`)  
*get sum of counts*
- `channel_t` `getChannel` ()  
*the new virtual channel*

### Protected Member Functions

- `bool` `next_impl` (`std::vector`< `Tag` > &`incoming_tags`, `timestamp_t` `begin_time`, `timestamp_t` `end_time`) override  
*update iterator state*
- `void` `clear_impl` () override  
*clear `Iterator` state.*

## Friends

- class [CombinerImpl](#)

## Additional Inherited Members

### 8.3.1 Detailed Description

Combine some channels in a virtual channel which has a tick for each tick in the input channels.

This iterator can be used to get aggregation channels, eg if you want to monitor the countrate of the sum of two channels.

### 8.3.2 Constructor & Destructor Documentation

#### 8.3.2.1 `Combiner::Combiner ( TimeTaggerBase * tagger, std::vector< channel_t > channels )`

construct a combiner

##### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>channels</i>	vector of channels to combine

#### 8.3.2.2 `Combiner::~~Combiner ( )`

### 8.3.3 Member Function Documentation

#### 8.3.3.1 `void Combiner::clear_impl ( ) [override], [protected], [virtual]`

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

#### 8.3.3.2 `channel_t Combiner::getChannel ( )`

the new virtual channel

This function returns the new allocated virtual channel. It can be used now in any new iterator.

8.3.3.3 `void Combiner::getData ( std::function< int64_t *(size_t)> array_out )`

get sum of counts

For reference, this iterators sums up how much ticks are generated because of which input channel. So this functions returns an array with one value per input channel.

8.3.3.4 `bool Combiner::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time )` `[override]`, `[protected]`, `[virtual]`

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

## 8.3.4 Friends And Related Function Documentation

8.3.4.1 `friend class CombinerImpl` `[friend]`

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.4 ConstantFractionDiscriminator Class Reference

a virtual CFD implementation which returns the mean time between a raising and a falling pair of edges

```
#include <Iterators.h>
```

Inheritance diagram for ConstantFractionDiscriminator:



## Public Member Functions

- `ConstantFractionDiscriminator` (`TimeTaggerBase *tagger`, `std::vector< channel_t > channels`, `timestamp_t search_window`)  
*constructor of a `ConstantFractionDiscriminator`*
- `~ConstantFractionDiscriminator` ()
- `std::vector< channel_t > getChannels` ()  
*the list of new virtual channels*

## Protected Member Functions

- `bool next_impl` (`std::vector< Tag > &incoming_tags`, `timestamp_t begin_time`, `timestamp_t end_time`) override  
*update iterator state*
- `void on_start` () override  
*callback when the measurement class is started*

## Friends

- class `ConstantFractionDiscriminatorImpl`

## Additional Inherited Members

### 8.4.1 Detailed Description

a virtual CFD implementation which returns the mean time between a raising and a falling pair of edges

### 8.4.2 Constructor & Destructor Documentation

#### 8.4.2.1 `ConstantFractionDiscriminator::ConstantFractionDiscriminator ( TimeTaggerBase * tagger, std::vector< channel_t > channels, timestamp_t search_window )`

constructor of a `ConstantFractionDiscriminator`

## Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>channels</i>	list of channels for the CFD, the formers of the raising+falling pairs must be given
<i>search_window</i>	interval for the CFD window, must be positive

## 8.4.2.2 ConstantFractionDiscriminator::~~ConstantFractionDiscriminator ( )

## 8.4.3 Member Function Documentation

## 8.4.3.1 std::vector&lt;channel\_t&gt; ConstantFractionDiscriminator::getChannels ( )

the list of new virtual channels

This function returns the list of new allocated virtual channels. It can be used now in any new measurement class.

## 8.4.3.2 bool ConstantFractionDiscriminator::next\_impl ( std::vector&lt; Tag &gt; &amp; incoming\_tags, timestamp\_t begin\_time, timestamp\_t end\_time ) [override], [protected], [virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

## Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

## Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

## 8.4.3.3 void ConstantFractionDiscriminator::on\_start ( ) [override], [protected], [virtual]

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

## 8.4.4 Friends And Related Function Documentation

### 8.4.4.1 friend class ConstantFractionDiscriminatorImpl [friend]

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.5 Correlation Class Reference

cross-correlation between two channels

```
#include <Iterators.h>
```

Inheritance diagram for Correlation:



### Public Member Functions

- [Correlation](#) ([TimeTaggerBase](#) \*tagger, [channel\\_t](#) channel\_1, [channel\\_t](#) channel\_2=CHANNEL\_UNUSED, [timestamp\\_t](#) binwidth=1000, int n\_bins=1000)  
*constructor of a correlation measurement*
- [~Correlation](#) ()
- void [getData](#) (std::function< int32\_t \*(size\_t)> array\_out)  
*returns a one-dimensional array of size n\_bins containing the histogram*
- void [getDataNormalized](#) (std::function< double \*(size\_t)> array\_out)  
*get the histogram - normalized such that a perfectly uncorrelated signals would be flat at a height of one*
- void [getIndex](#) (std::function< long long \*(size\_t)> array\_out)  
*returns a vector of size n\_bins containing the time bins in ps*

### Protected Member Functions

- bool [next\\_impl](#) (std::vector< [Tag](#) > &incoming\_tags, [timestamp\\_t](#) begin\_time, [timestamp\\_t](#) end\_time) override  
*update iterator state*
- void [clear\\_impl](#) () override  
*clear [Iterator](#) state.*

## Friends

- class [CorrelationImpl](#)

## Additional Inherited Members

### 8.5.1 Detailed Description

cross-correlation between two channels

Accumulates time differences between clicks on two channels into a histogram, where all ticks are considered both as start and stop clicks and both positive and negative time differences are considered. The histogram is determined by the number of total bins and the binwidth.

### 8.5.2 Constructor & Destructor Documentation

**8.5.2.1** `Correlation::Correlation ( TimeTaggerBase * tagger, channel_t channel_1, channel_t channel_2 = CHANNEL_UNUSED, timestamp_t binwidth = 1000, int n_bins = 1000 )`

constructor of a correlation measurement

If *channel\_2* is left empty or set to CHANNEL\_UNUSED, an auto-correlation measurement is performed. This is the same as setting *channel\_2* = *channel\_1*.

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>channel_1</i>	first click channel
<i>channel_2</i>	second click channel
<i>binwidth</i>	width of one histogram bin in ps
<i>n_bins</i>	the number of bins in the resulting histogram

**8.5.2.2** `Correlation::~~Correlation ( )`

### 8.5.3 Member Function Documentation

**8.5.3.1** `void Correlation::clear_impl ( )` `[override]`, `[protected]`, `[virtual]`

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

**8.5.3.2** void Correlation::getData ( std::function< int32\_t \*(size\_t)> array\_out )

returns a one-dimensional array of size n\_bins containing the histogram

**8.5.3.3** void Correlation::getDataNormalized ( std::function< double \*(size\_t)> array\_out )

get the histogram - normalized such that a perfectly uncorrelated signals would be flat at a height of one

**8.5.3.4** void Correlation::getIndex ( std::function< long long \*(size\_t)> array\_out )

returns a vector of size n\_bins containing the time bins in ps

**8.5.3.5** bool Correlation::next\_impl ( std::vector< Tag > & incoming\_tags, timestamp\_t begin\_time, timestamp\_t end\_time ) [override], [protected], [virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

#### Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

#### Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

## 8.5.4 Friends And Related Function Documentation

**8.5.4.1** friend class CorrelationImpl [friend]

The documentation for this class was generated from the following file:

- [Iterators.h](#)




## 8.6 CountBetweenMarkers Class Reference

a simple counter where external marker signals determine the bins

```
#include <Iterators.h>
```

Inheritance diagram for CountBetweenMarkers:



### Public Member Functions

- `CountBetweenMarkers` (`TimeTaggerBase` \*tagger, `channel_t` click\_channel, `channel_t` begin\_channel, `channel_t` end\_channel=`CHANNEL_UNUSED`, `int32_t` n\_values=1000)  
*constructor of CountBetweenMarkers*
- `~CountBetweenMarkers` ()
- `bool` ready ()  
*tbd*
- `void` getData (`std::function`< `int32_t` \*(`size_t`)> array\_out)  
*tbd*
- `void` getBinWidths (`std::function`< `long long` \*(`size_t`)> array\_out)  
*fetches the widths of each bins*
- `void` getIndex (`std::function`< `long long` \*(`size_t`)> array\_out)  
*fetches the starting time of each bin*

### Protected Member Functions

- `bool` next\_impl (`std::vector`< `Tag` > &incoming\_tags, `timestamp_t` begin\_time, `timestamp_t` end\_time) override  
*update iterator state*
- `void` clear\_impl () override  
*clear iterator state.*

### Friends

- class `CountBetweenMarkersImpl`

## Additional Inherited Members

### 8.6.1 Detailed Description

a simple counter where external marker signals determine the bins

[Counter](#) with external signals that trigger beginning and end of each counter accumulation. This can be used to implement counting triggered by a pixel clock and gated counting. The thread waits for the first time tag on the 'begin\_channel', then begins counting time tags on the 'click\_channel'. It ends counting when a tag on the 'end\_channel' is detected.

### 8.6.2 Constructor & Destructor Documentation

**8.6.2.1** `CountBetweenMarkers::CountBetweenMarkers ( TimeTaggerBase * tagger, channel_t click_channel, channel_t begin_channel, channel_t end_channel = CHANNEL_UNUSED, int32_t n_values = 1000 )`

constructor of [CountBetweenMarkers](#)

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>click_channel</i>	channel that increases the count
<i>begin_channel</i>	channel that triggers beginning of counting and stepping to the next value
<i>end_channel</i>	channel that triggers end of counting
<i>n_values</i>	the number of counter values to be stored

**8.6.2.2** `CountBetweenMarkers::~~CountBetweenMarkers ( )`

### 8.6.3 Member Function Documentation

**8.6.3.1** `void CountBetweenMarkers::clear_impl ( ) [override], [protected], [virtual]`

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

**8.6.3.2** `void CountBetweenMarkers::getBinWidths ( std::function< long long *(size_t)> array_out )`

fetches the widths of each bins

**8.6.3.3** `void CountBetweenMarkers::getData ( std::function< int32_t *(size_t)> array_out )`

tbd

8.6.3.4 `void CountBetweenMarkers::getIndex ( std::function< long long *(size_t)> array_out )`

fetches the starting time of each bin

8.6.3.5 `bool CountBetweenMarkers::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time )` `[override]`, `[protected]`, `[virtual]`

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

8.6.3.6 `bool CountBetweenMarkers::ready ( )`

tbd

## 8.6.4 Friends And Related Function Documentation

8.6.4.1 `friend class CountBetweenMarkersImpl` `[friend]`

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.7 Counter Class Reference

a simple counter on one or more channels

```
#include <Iterators.h>
```

Inheritance diagram for Counter:



## Public Member Functions

- `Counter` (`TimeTaggerBase *tagger`, `std::vector< channel_t > channels`, `timestamp_t binwidth=1000000000`, `int32_t n_values=1`)  
*construct a counter*
- `~Counter` ()
- `void getData` (`std::function< int32_t *(size_t, size_t)> array_out`)  
*get counts*
- `void getIndex` (`std::function< long long *(size_t)> array_out`)

## Protected Member Functions

- `bool next_impl` (`std::vector< Tag > &incoming_tags`, `timestamp_t begin_time`, `timestamp_t end_time`) override  
*update iterator state*
- `void clear_impl` () override  
*clear `Iterator` state.*
- `void on_start` () override  
*callback when the measurement class is started*

## Friends

- class `CounterImpl`

## Additional Inherited Members

### 8.7.1 Detailed Description

a simple counter on one or more channels

`Counter` with fixed binwidth and circular buffer output. This class is suitable to generate a time trace of the count rate on one or more channels. The thread repeatedly counts clicks on a single channel over a given time interval and stores the results in a two-dimensional array. The array is treated as a circular buffer. I.e., once the array is full, each new value shifts all previous values one element to the left.

## 8.7.2 Constructor & Destructor Documentation

**8.7.2.1** `Counter::Counter ( TimeTaggerBase * tagger, std::vector< channel_t > channels, timestamp_t binwidth = 1000000000, int32_t n_values = 1 )`

construct a counter

Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>channels</i>	channels to count on
<i>binwidth</i>	counts are accumulated for binwidth picoseconds
<i>n_values</i>	number of counter values stored (for each channel)

**8.7.2.2** `Counter::~Counter ( )`

## 8.7.3 Member Function Documentation

**8.7.3.1** `void Counter::clear_impl ( )` `[override]`, `[protected]`, `[virtual]`

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear\\_impl\(\)](#) method to reset its internal state. The [clear\\_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

**8.7.3.2** `void Counter::getData ( std::function< int32_t *(size_t, size_t)> array_out )`

get counts

the counts are copied to a newly allocated allocated memory, an the pointer to this location is returned.

**8.7.3.3** `void Counter::getIndex ( std::function< long long *(size_t)> array_out )`

**8.7.3.4** `bool Counter::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time )` `[override]`, `[protected]`, `[virtual]`

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

#### Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

**8.7.3.5** `void Counter::on_start( )` `[override]`, `[protected]`, `[virtual]`

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

### 8.7.4 Friends And Related Function Documentation

**8.7.4.1** `friend class CounterImpl` `[friend]`

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.8 Countrate Class Reference

count rate on one or more channels

```
#include <Iterators.h>
```

Inheritance diagram for Countrate:



## Public Member Functions

- [Countrate](#) ([TimeTaggerBase](#) \*[tagger](#), [std::vector](#)< [channel\\_t](#) > [channels](#))  
*constructor of [Countrate](#)*
- [~Countrate](#) ()
- void [getData](#) ([std::function](#)< double \*([size\\_t](#))> [array\\_out](#))  
*get the count rates*
- void [getCountsTotal](#) ([std::function](#)< [int64\\_t](#) \*([size\\_t](#))> [array\\_out](#))  
*get the total amount of events*

## Protected Member Functions

- bool [next\\_impl](#) ([std::vector](#)< [Tag](#) > &[incoming\\_tags](#), [timestamp\\_t](#) [begin\\_time](#), [timestamp\\_t](#) [end\\_time](#)) override  
*update iterator state*
- void [clear\\_impl](#) () override  
*clear [Iterator](#) state.*
- void [on\\_start](#) () override  
*callback when the measurement class is started*

## Friends

- class [CountrateImpl](#)

## Additional Inherited Members

### 8.8.1 Detailed Description

count rate on one or more channels

Measures the average count rate on one or more channels. Specifically, it counts incoming clicks and determines the time between the initial click and the latest click. The number of clicks divided by the time corresponds to the average countrate since the initial click.

### 8.8.2 Constructor & Destructor Documentation

#### 8.8.2.1 [Countrate::Countrate](#) ( [TimeTaggerBase](#) \* [tagger](#), [std::vector](#)< [channel\\_t](#) > [channels](#) )

constructor of [Countrate](#)

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>channels</i>	the channels to count on

### 8.8.2.2 Countrate::~~Countrate ( )

## 8.8.3 Member Function Documentation

### 8.8.3.1 void Countrate::clear\_impl ( ) [override],[protected],[virtual]

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear\\_impl\(\)](#) method to reset its internal state. The [clear\\_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

### 8.8.3.2 void Countrate::getCountsTotal ( std::function< int64\_t \*(size\_t)> array\_out )

get the total amount of events

Returns the total amount of events per channel as an array.

### 8.8.3.3 void Countrate::getData ( std::function< double \*(size\_t)> array\_out )

get the count rates

Returns the average rate of events per second per channel as an array.

### 8.8.3.4 bool Countrate::next\_impl ( std::vector< Tag > & incoming\_tags, timestamp\_t begin\_time, timestamp\_t end\_time ) [override],[protected],[virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

#### Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

#### Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).



8.8.3.5 `void Countrate::on_start ( ) [override],[protected],[virtual]`

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

## 8.8.4 Friends And Related Function Documentation

8.8.4.1 `friend class CountrateImpl [friend]`

The documentation for this class was generated from the following file:

- [Iterators.h](#)

## 8.9 CustomLogger Class Reference

```
#include <TimeTagger.h>
```

### Public Member Functions

- [CustomLogger](#) ()
- virtual [~CustomLogger](#) ()
- void [enable](#) ()
- void [disable](#) ()
- virtual void [Log](#) (int level, const std::string &msg)=0

### 8.9.1 Constructor & Destructor Documentation

8.9.1.1 `CustomLogger::CustomLogger ( )`

8.9.1.2 `virtual CustomLogger::~~CustomLogger ( ) [virtual]`

### 8.9.2 Member Function Documentation

8.9.2.1 `void CustomLogger::disable ( )`

8.9.2.2 `void CustomLogger::enable ( )`

8.9.2.3 `virtual void CustomLogger::Log ( int level, const std::string & msg ) [pure virtual]`


The documentation for this class was generated from the following file:

- [TimeTagger.h](#)

## 8.10 CustomMeasurementBase Class Reference

```
#include <Iterators.h>
```

Inheritance diagram for CustomMeasurementBase:



### Public Member Functions

- [~CustomMeasurementBase](#) () override
- void [register\\_channel](#) ([channel\\_t](#) channel)
- void [unregister\\_channel](#) ([channel\\_t](#) channel)
- void [finalize\\_init](#) ()
- bool [is\\_running](#) () const
- void [\\_lock](#) ()
- void [\\_unlock](#) ()

### Static Public Member Functions

- static void [stop\\_all\\_custom\\_measurements](#) ()

### Protected Member Functions

- [CustomMeasurementBase](#) ([TimeTaggerBase](#) \*tagger)
- virtual bool [next\\_impl](#) (std::vector< [Tag](#) > &incoming\_tags, [timestamp\\_t](#) begin\_time, [timestamp\\_t](#) end\_time) override  
*update iterator state*
- virtual void [next\\_impl\\_cs](#) (void \*tags\_ptr, uint64\_t num\_tags, [timestamp\\_t](#) begin\_time, [timestamp\\_t](#) end\_time)
- virtual void [clear\\_impl](#) () override  
*clear [Iterator](#) state.*
- virtual void [on\\_start](#) () override  
*callback when the measurement class is started*
- virtual void [on\\_stop](#) () override  
*callback when the measurement class is stopped*

## Additional Inherited Members

### 8.10.1 Constructor & Destructor Documentation

8.10.1.1 CustomMeasurementBase::CustomMeasurementBase ( TimeTaggerBase \* *tagger* ) [protected]

8.10.1.2 CustomMeasurementBase::~CustomMeasurementBase ( ) [override]

### 8.10.2 Member Function Documentation

8.10.2.1 void CustomMeasurementBase::\_lock ( )

8.10.2.2 void CustomMeasurementBase::\_unlock ( )

8.10.2.3 virtual void CustomMeasurementBase::clear\_impl ( ) [override],[protected],[virtual]

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear\\_impl\(\)](#) method to reset its internal state. The [clear\\_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

8.10.2.4 void CustomMeasurementBase::finalize\_init ( )

8.10.2.5 bool CustomMeasurementBase::is\_running ( ) const

8.10.2.6 virtual bool CustomMeasurementBase::next\_impl ( std::vector< Tag > & *incoming\_tags*, timestamp\_t *begin\_time*, timestamp\_t *end\_time* ) [override],[protected],[virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

#### Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

#### Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

8.10.2.7 `virtual void CustomMeasurementBase::next_impl_cs ( void * tags_ptr, uint64_t num_tags, timestamp_t begin_time, timestamp_t end_time )` [protected],[virtual]

8.10.2.8 `virtual void CustomMeasurementBase::on_start ( )` [override],[protected],[virtual]

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

8.10.2.9 `virtual void CustomMeasurementBase::on_stop ( )` [override],[protected],[virtual]

callback when the measurement class is stopped

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

8.10.2.10 `void CustomMeasurementBase::register_channel ( channel_t channel )`

8.10.2.11 `static void CustomMeasurementBase::stop_all_custom_measurements ( )` [static]

8.10.2.12 `void CustomMeasurementBase::unregister_channel ( channel_t channel )`

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.11 DelayedChannel Class Reference

a simple delayed queue

```
#include <Iterators.h>
```

Inheritance diagram for DelayedChannel:



## Public Member Functions

- [DelayedChannel](#) ([TimeTaggerBase](#) \*tagger, [channel\\_t](#) input\_channel, [timestamp\\_t](#) delay)  
*constructor of a [DelayedChannel](#)*
- [DelayedChannel](#) ([TimeTaggerBase](#) \*tagger, [std::vector](#)< [channel\\_t](#) > input\_channels, [timestamp\\_t](#) delay)  
*constructor of a [DelayedChannel](#) for delaying many channels at once*
- [~DelayedChannel](#) ()
- [channel\\_t](#) getChannel ()  
*the first new virtual channel*
- [std::vector](#)< [channel\\_t](#) > getChannels ()  
*the new virtual channels*
- void [setDelay](#) ([timestamp\\_t](#) delay)  
*set the delay time delay for the cloned tags in the virtual channels. A negative delay will delay all other events.*

## Protected Member Functions

- bool [next\\_impl](#) ([std::vector](#)< [Tag](#) > &incoming\_tags, [timestamp\\_t](#) begin\_time, [timestamp\\_t](#) end\_time) override  
*update iterator state*
- void [on\\_start](#) () override  
*callback when the measurement class is started*

## Friends

- class [DelayedChannelImpl](#)

## Additional Inherited Members

### 8.11.1 Detailed Description

a simple delayed queue

A simple first-in first-out queue of delayed event timestamps.

### 8.11.2 Constructor & Destructor Documentation

#### 8.11.2.1 DelayedChannel::DelayedChannel ( [TimeTaggerBase](#) \* tagger, [channel\\_t](#) input\_channel, [timestamp\\_t](#) delay )

constructor of a [DelayedChannel](#)

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>input_channel</i>	channel which is delayed
<i>delay</i>	amount of time to delay

### 8.11.2.2 DelayedChannel::DelayedChannel ( TimeTaggerBase \* *tagger*, std::vector< channel\_t > *input\_channels*, timestamp\_t *delay* )

constructor of a [DelayedChannel](#) for delaying many channels at once

This function is not exposed to Python/C#/Matlab/Labview

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>input_channels</i>	channels which will be delayed
<i>delay</i>	amount of time to delay

### 8.11.2.3 DelayedChannel::~~DelayedChannel ( )

## 8.11.3 Member Function Documentation

### 8.11.3.1 channel\_t DelayedChannel::getChannel ( )

the first new virtual channel

This function returns the first of the new allocated virtual channels. It can be used now in any new iterator.

### 8.11.3.2 std::vector<channel\_t> DelayedChannel::getChannels ( )

the new virtual channels

This function returns the new allocated virtual channels. It can be used now in any new iterator.

### 8.11.3.3 bool DelayedChannel::next\_impl ( std::vector< Tag > & *incoming\_tags*, timestamp\_t *begin\_time*, timestamp\_t *end\_time* ) [override],[protected],[virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

#### Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

**Returns**

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

**8.11.3.4** `void DelayedChannel::on_start( )` `[override]`, `[protected]`, `[virtual]`

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

**8.11.3.5** `void DelayedChannel::setDelay( timestamp_t delay )`

set the delay time delay for the cloned tags in the virtual channels. A negative delay will delay all other events.

Note: When the delay is the same or greater than the previous value all incoming tags will be visible at virtual channel. By applying a shorter delay time, the tags stored in the local buffer will be flushed and won't be visible in the virtual channel.

## 8.11.4 Friends And Related Function Documentation

**8.11.4.1** `friend class DelayedChannelImpl` `[friend]`

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.12 Dump Class Reference

dump all time tags to a file

```
#include <Iterators.h>
```

Inheritance diagram for Dump:



## Public Member Functions

- [Dump](#) ([TimeTaggerBase](#) \**tagger*, `std::string` *filename*, `int64_t` *max\_tags*, `std::vector`< [channel\\_t](#) > *channels*=`std::vector`< [channel\\_t](#) >())  
*constructor of a [Dump](#) thread*
- [~Dump](#) ()  
*tbd*

## Protected Member Functions

- `bool` [next\\_impl](#) (`std::vector`< [Tag](#) > &*incoming\_tags*, [timestamp\\_t](#) *begin\_time*, [timestamp\\_t](#) *end\_time*) override  
*update iterator state*
- `void` [clear\\_impl](#) () override  
*clear [Iterator](#) state.*
- `void` [on\\_start](#) () override  
*callback when the measurement class is started*
- `void` [on\\_stop](#) () override  
*callback when the measurement class is stopped*

## Friends

- class [DumpImpl](#)

## Additional Inherited Members

### 8.12.1 Detailed Description

dump all time tags to a file

**Deprecated** use [FileWriter](#)

### 8.12.2 Constructor & Destructor Documentation

8.12.2.1 `Dump::Dump ( TimeTaggerBase * tagger, std::string filename, int64_t max_tags, std::vector< channel\_t > channels = std::vector< channel\_t > () )`

constructor of a [Dump](#) thread

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>filename</i>	name of the file to dump to
<i>max_tags</i>	stop after this number of tags has been dumped. Negative values will dump forever
<i>channels</i>	channels which are dumped to the file (when empty or not passed all active channels are dumped)



## 8.12.2.2 Dump::~Dump ( )

tbd

## 8.12.3 Member Function Documentation

## 8.12.3.1 void Dump::clear\_impl ( ) [override],[protected],[virtual]

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear\\_impl\(\)](#) method to reset its internal state. The [clear\\_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

## 8.12.3.2 bool Dump::next\_impl ( std::vector&lt; Tag &gt; &amp; incoming\_tags, timestamp\_t begin\_time, timestamp\_t end\_time ) [override],[protected],[virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

## Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

## Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

## 8.12.3.3 void Dump::on\_start ( ) [override],[protected],[virtual]

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

## 8.12.3.4 void Dump::on\_stop ( ) [override],[protected],[virtual]

callback when the measurement class is stopped

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

## 8.12.4 Friends And Related Function Documentation

### 8.12.4.1 friend class DumpImpl [friend]

The documentation for this class was generated from the following file:

- [Iterators.h](#)

## 8.13 Event Struct Reference

```
#include <Iterators.h>
```

### Public Attributes

- [timestamp\\_t time](#)
- [State state](#)

### 8.13.1 Member Data Documentation

#### 8.13.1.1 State Event::state

#### 8.13.1.2 timestamp\_t Event::time

The documentation for this struct was generated from the following file:


- [Iterators.h](#)

## 8.14 EventGenerator Class Reference

Generate predefined events in a virtual channel relative to a trigger event.

```
#include <Iterators.h>
```

Inheritance diagram for EventGenerator:



## Public Member Functions

- [EventGenerator](#) ([TimeTaggerBase](#) \**tagger*, [channel\\_t](#) *trigger\_channel*, [std::vector](#)< [timestamp\\_t](#) > *pattern*, [uint64\\_t](#) *trigger\_divider*=1, [uint64\\_t](#) *divider\_offset*=0, [channel\\_t](#) *stop\_channel*=[CHANNEL\\_UNUSED](#))  
*construct a event generator*
- [~EventGenerator](#) ()
- [channel\\_t](#) *getChannel* ()  
*the new virtual channel*

## Protected Member Functions

- [bool](#) *next\_impl* ([std::vector](#)< [Tag](#) > &*incoming\_tags*, [timestamp\\_t](#) *begin\_time*, [timestamp\\_t](#) *end\_time*) *override*  
*update iterator state*
- [void](#) *clear\_impl* () *override*  
*clear [Iterator](#) state.*
- [void](#) *on\_start* () *override*  
*callback when the measurement class is started*

## Friends

- class [EventGeneratorImpl](#)

## Additional Inherited Members

### 8.14.1 Detailed Description

Generate predefined events in a virtual channel relative to a trigger event.

This iterator can be used to generate a predefined series of events, the pattern, relative to a trigger event on a defined channel. A *trigger\_divider* can be used to fire the pattern not on every, but on every n'th trigger received. The *trigger\_offset* can be used to select on which of the triggers the pattern will be generated when *trigger\_divider* is greater than 1. To abort the pattern being generated, a *stop\_channel* can be defined. In case it is the very same as the *trigger\_channel*, the subsequent generated patterns will not overlap.

### 8.14.2 Constructor & Destructor Documentation

- 8.14.2.1 [EventGenerator::EventGenerator](#) ( [TimeTaggerBase](#) \* *tagger*, [channel\\_t](#) *trigger\_channel*, [std::vector](#)< [timestamp\\_t](#) > *pattern*, [uint64\\_t](#) *trigger\_divider* = 1, [uint64\\_t](#) *divider\_offset* = 0, [channel\\_t](#) *stop\_channel* = [CHANNEL\\_UNUSED](#) )

construct a event generator

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>trigger_channel</i>	trigger for generating the pattern
<i>pattern</i>	vector of time stamp generated relativ to the trigger event
<i>trigger_divider</i>	establishes every how many trigger events a pattern is generated
<i>divider_offset</i>	the offset of the divided trigger when the pattern shall be emitted
<i>stop_channel</i>	channel on which a received event will stop all pending patterns from being generated

#### 8.14.2.2 EventGenerator::~~EventGenerator ( )

### 8.14.3 Member Function Documentation

#### 8.14.3.1 void EventGenerator::clear\_impl ( ) [override],[protected],[virtual]

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear\\_impl\(\)](#) method to reset its internal state. The [clear\\_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

#### 8.14.3.2 channel\_t EventGenerator::getChannel ( )

the new virtual channel

This function returns the new allocated virtual channel. It can be used now in any new iterator.

#### 8.14.3.3 bool EventGenerator::next\_impl ( std::vector< Tag > & incoming\_tags, timestamp\_t begin\_time, timestamp\_t end\_time ) [override],[protected],[virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

##### Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

##### Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

#### 8.14.3.4 void EventGenerator::on\_start ( ) [override],[protected],[virtual]

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

### 8.14.4 Friends And Related Function Documentation

#### 8.14.4.1 friend class EventGeneratorImpl [friend]

The documentation for this class was generated from the following file:

- [Iterators.h](#)

## 8.15 FastBinning Class Reference

```
#include <Iterators.h>
```

### Public Types

- enum [Mode](#) {  
    [Mode::ConstZero](#), [Mode::Dividend](#), [Mode::PowerOfTwo](#), [Mode::FixedPoint\\_32](#),  
    [Mode::FixedPoint\\_64](#), [Mode::Divide\\_32](#), [Mode::Divide\\_64](#) }

### Public Member Functions

- [FastBinning](#) ()
- [FastBinning](#) (uint64\_t divisor, uint64\_t max\_duration\_)
- template<Mode mode>  
    uint64\_t [divide](#) (uint64\_t duration) const
- [Mode](#) [getMode](#) () const

### 8.15.1 Detailed Description

Helper class for fast division with a constant divisor. It chooses the method on initialization time and precompile the evaluation functions for all methods.

### 8.15.2 Member Enumeration Documentation

#### 8.15.2.1 enum FastBinning::Mode [strong]

Enumerator

***ConstZero***

***Dividend***

***PowerOfTwo***

***FixedPoint\_32***

***FixedPoint\_64***

***Divide\_32***

***Divide\_64***

### 8.15.3 Constructor & Destructor Documentation

8.15.3.1 `FastBinning::FastBinning ( )` `[inline]`

8.15.3.2 `FastBinning::FastBinning ( uint64_t divisor, uint64_t max_duration_ )`

### 8.15.4 Member Function Documentation

8.15.4.1 `template<Mode mode> uint64_t FastBinning::divide ( uint64_t duration ) const` `[inline]`

8.15.4.2 `Mode FastBinning::getMode ( ) const` `[inline]`

The documentation for this class was generated from the following file:

- [Iterators.h](#)

## 8.16 FileReader Class Reference

```
#include <Iterators.h>
```

### Public Member Functions

- [FileReader](#) (std::vector< std::string > filenames)
- [FileReader](#) (const std::string &filename)
- [~FileReader](#) ()
- bool [hasData](#) ()
- [TimeTagStreamBuffer](#) [getData](#) (uint64\_t n\_events)
- bool [getDataRaw](#) (std::vector< [Tag](#) > &tag\_buffer)
- std::string [getConfiguration](#) ()
- std::string [getLastMarker](#) ()

### Friends

- class [FileReaderImpl](#)

### 8.16.1 Detailed Description

Reads tags from the disk files, which has been created by [FileWriter](#). Its usage is compatible with the [TimeTagStream](#).

### 8.16.2 Constructor & Destructor Documentation

8.16.2.1 `FileReader::FileReader ( std::vector< std::string > filenames )`

Creates a file reader with the given filename. The file reader automatically continues to read split [FileWriter](#) Streams. In case multiple filenames are given, the files will be read in successively.

## Parameters

<i>filenames</i>	list of files to read
------------------	-----------------------

8.16.2.2 `FileReader::FileReader ( const std::string & filename )`

Creates a file reader with the given filename. The file reader automatically continues to read split [FileWriter](#) Streams

## Parameters

<i>filename</i>	file to read
-----------------	--------------

8.16.2.3 `FileReader::~FileReader ( )`

## 8.16.3 Member Function Documentation

8.16.3.1 `std::string FileReader::getConfiguration ( )`

Fetches the overall configuration status of the Time Tagger object, which was serialized in the current file.

## Returns

a JSON serialized string with all configuration and status flags.

8.16.3.2 `TimeTagStreamBuffer FileReader::getData ( uint64_t n_events )`

Fetches and delete the next tags from the internal buffer. Every tag is returned exactly once. If less than `n_events` are returned, the reader is at the end-of-files.

## Parameters

<i>n_events</i>	maximum amount of elements to fetch
-----------------	-------------------------------------

## Returns

a [TimeTagStreamBuffer](#) with up to `n_events` events

8.16.3.3 `bool FileReader::getDataRaw ( std::vector< Tag > & tag_buffer )`

Low level file reading. This function will return the next non-empty buffer in a raw format.

## Parameters

<i>tag_buffer</i>	a buffer, which will be filled with the new events
-------------------	--

**Returns**

true if fetching the data was successfully

**8.16.3.4** `std::string FileReader::getLastMarker ( )`

return the last processed marker from the file.

**Returns**

the last marker from the file

**8.16.3.5** `bool FileReader::hasData ( )`

Checks if there are still events in the [FileReader](#)

**Returns**

false if no more events can be read from this [FileReader](#)

**8.16.4 Friends And Related Function Documentation****8.16.4.1** `friend class FileReaderImpl [friend]`

The documentation for this class was generated from the following file:


- [Iterators.h](#)

**8.17 FileWriter Class Reference**

compresses and stores all time tags to a file

```
#include <Iterators.h>
```

Inheritance diagram for FileWriter:





## Public Member Functions

- [FileWriter](#) ([TimeTaggerBase](#) \**tagger*, const std::string &filename, std::vector< [channel\\_t](#) > channels)  
*constructor of a [FileWriter](#)*
- [~FileWriter](#) ()
- void [split](#) (const std::string &new\_filename="")
- void [setMaxFileSize](#) (uint64\_t max\_file\_size)
- uint64\_t [getMaxFileSize](#) ()
- uint64\_t [getTotalEvents](#) ()
- uint64\_t [getTotalSize](#) ()
- void [setMarker](#) (const std::string &marker)

## Protected Member Functions

- bool [next\\_impl](#) (std::vector< [Tag](#) > &incoming\_tags, [timestamp\\_t](#) begin\_time, [timestamp\\_t](#) end\_time) override  
*update iterator state*
- void [clear\\_impl](#) () override  
*clear [Iterator](#) state.*
- void [on\\_start](#) () override  
*callback when the measurement class is started*
- void [on\\_stop](#) () override  
*callback when the measurement class is stopped*

## Friends

- class [FileWriterImpl](#)

## Additional Inherited Members

### 8.17.1 Detailed Description

compresses and stores all time tags to a file

### 8.17.2 Constructor & Destructor Documentation

#### 8.17.2.1 [FileWriter::FileWriter](#) ( [TimeTaggerBase](#) \* *tagger*, const std::string & *filename*, std::vector< [channel\\_t](#) > *channels* )

constructor of a [FileWriter](#)

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>filename</i>	name of the file to store to
<i>channels</i>	channels which are stored to the file

### 8.17.2.2 `FileWriter::~FileWriter ( )`

## 8.17.3 Member Function Documentation

### 8.17.3.1 `void FileWriter::clear_impl ( )` [override],[protected],[virtual]

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

### 8.17.3.2 `uint64_t FileWriter::getMaxFileSize ( )`

fetches the maximum file size. Please see `setMaxFileSize` for more details.

#### Returns

the maximum file size in bytes

### 8.17.3.3 `uint64_t FileWriter::getTotalEvents ( )`

queries the total amount of events stored in all files

#### Returns

the total amount of events stored

### 8.17.3.4 `uint64_t FileWriter::getTotalSize ( )`

queries the total amount of bytes stored in all files

#### Returns

the total amount of bytes stored

### 8.17.3.5 `bool FileWriter::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time )` [override],[protected],[virtual]

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

## Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

## Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

**8.17.3.6** void FileWriter::on\_start ( ) [override],[protected],[virtual]

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

**8.17.3.7** void FileWriter::on\_stop ( ) [override],[protected],[virtual]

callback when the measurement class is stopped

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

**8.17.3.8** void FileWriter::setMarker ( const std::string & marker )

writes a marker in the file. While parsing the file, the last marker can be extracted again.

## Parameters

<i>marker</i>	the marker to write into the file
---------------	-----------------------------------

**8.17.3.9** void FileWriter::setMaxFileSize ( uint64\_t max\_file\_size )

Set the maximum file size on disk and so when the automatical split happens. Note: This is a rough limit, the actual file might be larger by one block.

## Parameters

<i>max_file_size</i>	new maximum file size in bytes
----------------------	--------------------------------

8.17.3.10 `void FileWriter::split ( const std::string & new_filename = " " )`

Close the current file and create a new one

#### Parameters

<code>new_filename</code>	filename of the new file. If empty, the old one will be used.
---------------------------	---

## 8.17.4 Friends And Related Function Documentation

8.17.4.1 `friend class FileWriterImpl` [*friend*]

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.18 Flim Class Reference

Fluorescence lifetime imaging.

```
#include <Iterators.h>
```

Inheritance diagram for Flim:



## Public Member Functions

- [Flim](#) ([TimeTaggerBase](#) \*tagger, [channel\\_t](#) start\_channel, [channel\\_t](#) click\_channel, [channel\\_t](#) pixel\_begin\_↵  
channel, [uint32\\_t](#) n\_pixels, [uint32\\_t](#) n\_bins, [timestamp\\_t](#) binwidth, [channel\\_t](#) pixel\_end\_channel=CHANN↵  
EL\_UNUSED, [channel\\_t](#) frame\_begin\_channel=CHANNEL\_UNUSED, [uint32\\_t](#) finish\_after\_outputframe=0,  
[uint32\\_t](#) n\_frame\_average=1, [bool](#) pre\_initialize=true)  
  
*construct a [Flim](#) measurement with a variety of high-level functionality*
- [~Flim](#) ()
- [void initialize](#) ()  
  
*initializes and starts measuring this [Flim](#) measurement*
- [void getReadyFrame](#) ([std::function](#)< [uint32\\_t](#) \*([size\\_t](#), [size\\_t](#))> array\_out, [int32\\_t](#) index=-1)  
  
*obtain for each pixel the histogram for the given frame index*
- [void getReadyFrameIntensity](#) ([std::function](#)< [float](#) \*([size\\_t](#))> array\_out, [int32\\_t](#) index=-1)  
  
*obtain an array of the pixel intensity of the given frame index*
- [void getCurrentFrame](#) ([std::function](#)< [uint32\\_t](#) \*([size\\_t](#), [size\\_t](#))> array\_out)  
  
*obtain for each pixel the histogram for the frame currently active*
- [void getCurrentFrameIntensity](#) ([std::function](#)< [float](#) \*([size\\_t](#))> array\_out)  
  
*obtain the array of the pixel intensities of the frame currently active*
- [void getSummedFrames](#) ([std::function](#)< [uint32\\_t](#) \*([size\\_t](#), [size\\_t](#))> array\_out, [bool](#) only\_ready\_frames=true,  
[bool](#) clear\_summed=false)  
  
*obtain for each pixel the histogram from all frames acquired so far*
- [void getSummedFramesIntensity](#) ([std::function](#)< [float](#) \*([size\\_t](#))> array\_out, [bool](#) only\_ready\_frames=true,  
[bool](#) clear\_summed=false)  
  
*obtain the array of the pixel intensities from all frames acquired so far*
- [FlimFrameInfo getReadyFrameEx](#) ([int32\\_t](#) index=-1)  
  
*obtain a frame information object, for the given frame index*
- [FlimFrameInfo getCurrentFrameEx](#) ()  
  
*obtain a frame information object, for the currently active frame*
- [FlimFrameInfo getSummedFramesEx](#) ([bool](#) only\_ready\_frames=true, [bool](#) clear\_summed=false)  
  
*obtain a frame information object, that represents the sum of all frames acquired so far.*
- [uint32\\_t getFramesAcquired](#) () const  
  
*total number of frames completed so far*
- [void getIndex](#) ([std::function](#)< [long long](#) \*([size\\_t](#))> array\_out)  
  
*a vector of size n\_bins containing the time bins in ps*

## Protected Member Functions

- [void on\\_frame\\_end](#) () override
- [void clear\\_impl](#) () override  
  
*clear [Iterator](#) state.*
- [uint32\\_t get\\_ready\\_index](#) ([int32\\_t](#) index)
- [virtual void frameReady](#) ([uint32\\_t](#) frame\_number, [std::vector](#)< [uint32\\_t](#) > &data, [std::vector](#)< [timestamp\\_t](#)↵  
\_t > &pixel\_begin\_times, [std::vector](#)< [timestamp\\_t](#) > &pixel\_end\_times, [timestamp\\_t](#) frame\_begin\_time,  
[timestamp\\_t](#) frame\_end\_time)

## Protected Attributes

- `std::vector< std::vector< uint32_t > >` [back\\_frames](#)
- `std::vector< std::vector< timestamp_t > >` [frame\\_begins](#)
- `std::vector< std::vector< timestamp_t > >` [frame\\_ends](#)
- `std::vector< uint32_t >` [pixels\\_completed](#)
- `std::vector< uint32_t >` [summed\\_frames](#)
- `std::vector< timestamp_t >` [accum\\_diffs](#)
- `uint32_t` [captured\\_frames](#)
- `uint32_t` [total\\_frames](#)
- `int32_t` [last\\_frame](#)
- `std::mutex` [swap\\_chain\\_lock](#)

### 8.18.1 Detailed Description

Fluorescence lifetime imaging.

Successively acquires  $n$  histograms (one for each pixel in the image), where each histogram is determined by the number of bins and the binwidth. Clicks that fall outside the histogram range are ignored.

Fluorescence-lifetime imaging microscopy or [Flim](#) is an imaging technique for producing an image based on the differences in the exponential decay rate of the fluorescence from a fluorescent sample.

Fluorescence lifetimes can be determined in the time domain by using a pulsed source. When a population of fluorophores is excited by an ultrashort or delta pulse of light, the time-resolved fluorescence will decay exponentially.

### 8.18.2 Constructor & Destructor Documentation

**8.18.2.1** `Flim::Flim ( TimeTaggerBase * tagger, channel_t start_channel, channel_t click_channel, channel_t pixel_begin_channel, uint32_t n_pixels, uint32_t n_bins, timestamp_t binwidth, channel_t pixel_end_channel = CHANNEL_UNUSED, channel_t frame_begin_channel = CHANNEL_UNUSED, uint32_t finish_after_outputframe = 0, uint32_t n_frame_average = 1, bool pre_initialize = true )`

construct a [Flim](#) measurement with a variety of high-level functionality

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>start_channel</i>	channel on which start clicks are received for the time differences histogramming
<i>click_channel</i>	channel on which clicks are received for the time differences histogramming
<i>pixel_begin_channel</i>	start of a pixel (histogram)
<i>n_pixels</i>	number of pixels (histograms) of one frame
<i>n_bins</i>	number of histogram bins for each pixel
<i>binwidth</i>	bin size in picoseconds
<i>pixel_end_channel</i>	end marker of a pixel - incoming clicks on the click_channel will be ignored afterwards
<i>frame_begin_channel</i>	(optional) start the frame, or reset the pixel index
<i>finish_after_outputframe</i>	(optional) sets the number of frames stored within the measurement class. After reaching the number, the measurement will stop. If the number is 0 (default value), one frame is stored and the measurement runs continuously.
<i>n_frame_average</i>	(optional) average multiple input frames into one output frame, default: 1
<i>pre_initialize</i>	(optional) initializes the measurement on constructing.

## 8.18.2.2 Flim::~~Flim ( )

## 8.18.3 Member Function Documentation

## 8.18.3.1 void Flim::clear\_impl ( ) [override], [protected], [virtual]

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear\\_impl\(\)](#) method to reset its internal state. The [clear\\_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [FlimAbstract](#).

## 8.18.3.2 virtual void Flim::frameReady ( uint32\_t frame\_number, std::vector&lt; uint32\_t &gt; &amp; data, std::vector&lt; timestamp\_t &gt; &amp; pixel\_begin\_times, std::vector&lt; timestamp\_t &gt; &amp; pixel\_end\_times, timestamp\_t frame\_begin\_time, timestamp\_t frame\_end\_time ) [protected], [virtual]

## 8.18.3.3 uint32\_t Flim::get\_ready\_index ( int32\_t index ) [protected]

## 8.18.3.4 void Flim::getCurrentFrame ( std::function&lt; uint32\_t \*(size\_t, size\_t)&gt; array\_out )

obtain for each pixel the histogram for the frame currently active

This function returns the histograms for all pixels of the currently active frame

## 8.18.3.5 FlimFrameInfo Flim::getCurrentFrameEx ( )

obtain a frame information object, for the currently active frame

This function returns the frame information object for the currently active frame

## 8.18.3.6 void Flim::getCurrentFrameIntensity ( std::function&lt; float \*(size\_t)&gt; array\_out )

obtain the array of the pixel intensities of the frame currently active

This function returns the intensities of all pixels of the currently active frame

The pixel intensity is defined by the number of counts acquired within the pixel divided by the respective integration time.

## 8.18.3.7 uint32\_t Flim::getFramesAcquired ( ) const [inline]

total number of frames completed so far

This function returns the amount of frames that have been completed so far, since the creation / last clear of the object.

## 8.18.3.8 void Flim::getIndex ( std::function&lt; long long \*(size\_t)&gt; array\_out )

a vector of size n\_bins containing the time bins in ps

This function returns a vector of size n\_bins containing the time bins in ps.

## 8.18.3.9 void Flim::getReadyFrame ( std::function&lt; uint32\_t \*(size\_t, size\_t)&gt; array\_out, int32\_t index = -1 )

obtain for each pixel the histogram for the given frame index

This function returns the histograms for all pixels according to the frame index given. If the index is -1, it will return the last frame, which has been completed. When finish\_after\_outputframe is 0, the index value must be -1. If index >= finish\_after\_outputframe, it will throw an error.

## Parameters

<i>array_out</i>	callback for the array output allocation
<i>index</i>	index of the frame to be obtained. if -1, the last frame which has been completed is returned

#### 8.18.3.10 **FlimFrameInfo** **Flim::getReadyFrameEx** ( *int32\_t index* = -1 )

obtain a frame information object, for the given frame index

This function returns a frame information object according to the index given. If the index is -1, it will return the last completed frame. When `finish_after_outputframe` is 0, index must be -1. If `index >= finish_after_outputframe`, it will throw an error.

## Parameters

<i>index</i>	index of the frame to be obtained. if -1, last completed frame will be returned
--------------	---

#### 8.18.3.11 **void** **Flim::getReadyFrameIntensity** ( *std::function< float \*(size\_t)> array\_out*, *int32\_t index* = -1 )

obtain an array of the pixel intensity of the given frame index

This function returns the intensities according to the frame index given. If the index is -1, it will return the intensity of the last frame, which has been completed. When `finish_after_outputframe` is 0, the index value must be -1. If `index >= finish_after_outputframe`, it will throw an error.

The pixel intensity is defined by the number of counts acquired within the pixel divided by the respective integration time.

## Parameters

<i>array_out</i>	callback for the array output allocation
<i>index</i>	index of the frame to be obtained. if -1, the last frame which has been completed is returned

#### 8.18.3.12 **void** **Flim::getSummedFrames** ( *std::function< uint32\_t \*(size\_t, size\_t)> array\_out*, *bool only\_ready\_frames* = *true*, *bool clear\_summed* = *false* )

obtain for each pixel the histogram from all frames acquired so far

This function returns the histograms for all pixels. The counts within the histograms are integrated since the start or the last clear of the measurement.

## Parameters

<i>array_out</i>	callback for the array output allocation
<i>only_ready_frames</i>	if true, only the finished frames are added. On false, the currently active frame is aggregated.
<i>clear_summed</i>	if true, the summed frames memory will be cleared.



**8.18.3.13** `FlimFrameInfo Flim::getSummedFramesEx ( bool only_ready_frames = true, bool clear_summed = false )`

obtain a frame information object, that represents the sum of all frames acquired so far.

This function returns the frame information object that represents the sum of all acquired frames.

Parameters

<i>only_ready_frames</i>	if true only the finished frames are added. On false, the currently active is aggregated.
<i>clear_summed</i>	if true, the summed frames memory will be reset and all frames stored prior will be unaccounted in the future.

**8.18.3.14** `void Flim::getSummedFramesIntensity ( std::function< float *(size_t)> array_out, bool only_ready_frames = true, bool clear_summed = false )`

obtain the array of the pixel intensities from all frames acquired so far

The pixel intensity is the number of counts within the pixel divided by the integration time.

This function returns the intensities of all pixels summed over all acquired frames.

Parameters

<i>array_out</i>	callback for the array output allocation
<i>only_ready_frames</i>	if true only the finished frames are added. On false, the currently active frame is aggregated.
<i>clear_summed</i>	if true, the summed frames memory will be cleared.

**8.18.3.15** `void Flim::initialize ( )`

initializes and starts measuring this [Flim](#) measurement

This function initializes the [Flim](#) measurement and starts executing it. It does nothing if preinitialized in the constructor is set to true.

**8.18.3.16** `void Flim::on_frame_end ( )` `[override]`, `[protected]`, `[virtual]`

Implements [FlimAbstract](#).

## 8.18.4 Member Data Documentation

**8.18.4.1** `std::vector<timestamp_t> Flim::accum_diffs` `[protected]`

**8.18.4.2** `std::vector<std::vector<uint32_t>> Flim::back_frames` `[protected]`

**8.18.4.3** `uint32_t Flim::captured_frames` `[protected]`

8.18.4.4 `std::vector<std::vector<timestamp_t>> Flim::frame_begins` [protected]

8.18.4.5 `std::vector<std::vector<timestamp_t>> Flim::frame_ends` [protected]

8.18.4.6 `int32_t Flim::last_frame` [protected]

8.18.4.7 `std::vector<uint32_t> Flim::pixels_completed` [protected]

8.18.4.8 `std::vector<uint32_t> Flim::summed_frames` [protected]

8.18.4.9 `std::mutex Flim::swap_chain_lock` [protected]

8.18.4.10 `uint32_t Flim::total_frames` [protected]


The documentation for this class was generated from the following file:

- [Iterators.h](#)

## 8.19 FlimAbstract Class Reference

```
#include <Iterators.h>
```

Inheritance diagram for FlimAbstract:



### Public Member Functions

- `FlimAbstract (TimeTaggerBase *tagger, channel_t start_channel, channel_t click_channel, channel_t pixel_begin_channel, uint32_t n_pixels, uint32_t n_bins, timestamp_t binwidth, channel_t pixel_end_channel=CHANNEL_UNUSED, channel_t frame_begin_channel=CHANNEL_UNUSED, uint32_t finish_after_outputframe=0, uint32_t n_frame_average=1, bool pre_initialize=true)`  
*construct a **FlimAbstract** object, **Flim** and **FlimBase** classes inherit from it*
- `~FlimAbstract ()`
- `bool isAcquiring () const`  
*tells if the data aquisition has finished reaching finish\_after\_outputframe*

## Protected Member Functions

- `template<FastBinning::Mode bin_mode>`  
`void process_tags (const std::vector< Tag > &incoming_tags)`
- `bool next_impl (std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time)` override  
*update iterator state*
- `void clear_impl ()` override  
*clear iterator state.*
- `void on_start ()` override  
*callback when the measurement class is started*
- `virtual void on_frame_end ()=0`

## Protected Attributes

- `const channel_t start_channel`
- `const channel_t click_channel`
- `const channel_t pixel_begin_channel`
- `const uint32_t n_pixels`
- `const uint32_t n_bins`
- `const timestamp_t binwidth`
- `const channel_t pixel_end_channel`
- `const channel_t frame_begin_channel`
- `const uint32_t finish_after_outputframe`
- `const uint32_t n_frame_average`
- `const timestamp_t time_window`
- `timestamp_t current_frame_begin`
- `timestamp_t current_frame_end`
- `bool acquiring {}`
- `bool frame_acquisition {}`
- `bool pixel_acquisition {}`
- `uint32_t pixels_processed {}`
- `uint32_t frames_completed {}`
- `uint32_t ticks {}`
- `size_t data_base {}`
- `std::vector< uint32_t > frame`
- `std::vector< timestamp_t > pixel_begins`
- `std::vector< timestamp_t > pixel_ends`
- `std::deque< timestamp_t > previous_starts`
- `FastBinning binner`
- `std::recursive_mutex acquisition_lock`
- `bool initialized`

### 8.19.1 Constructor & Destructor Documentation

- 8.19.1.1 `FlimAbstract::FlimAbstract ( TimeTaggerBase * tagger, channel_t start_channel, channel_t click_channel, channel_t pixel_begin_channel, uint32_t n_pixels, uint32_t n_bins, timestamp_t binwidth, channel_t pixel_end_channel = CHANNEL_UNUSED, channel_t frame_begin_channel = CHANNEL_UNUSED, uint32_t finish_after_outputframe = 0, uint32_t n_frame_average = 1, bool pre_initialize = true )`

construct a `FlimAbstract` object, `Flim` and `FlimBase` classes inherit from it

## Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>start_channel</i>	channel on which start clicks are received for the time differences histogramming
<i>click_channel</i>	channel on which clicks are received for the time differences histogramming
<i>pixel_begin_channel</i>	start of a pixel (histogram)
<i>n_pixels</i>	number of pixels (histograms) of one frame
<i>n_bins</i>	number of histogram bins for each pixel
<i>binwidth</i>	bin size in picoseconds
<i>pixel_end_channel</i>	end marker of a pixel - incoming clicks on the click_channel will be ignored afterwards
<i>frame_begin_channel</i>	(optional) start the frame, or reset the pixel index
<i>finish_after_outputframe</i>	(optional) sets the number of frames stored within the measurement class. After reaching the number, the measurement will stop. If the number is 0 (default value), one frame is stored and the measurement runs continuously.
<i>n_frame_average</i>	(optional) average multiple input frames into one output frame, default: 1
<i>pre_initialize</i>	(optional) initializes the measurement on constructing.

## 8.19.1.2 FlimAbstract::~~FlimAbstract ( )

## 8.19.2 Member Function Documentation

## 8.19.2.1 void FlimAbstract::clear\_impl ( ) [override], [protected], [virtual]

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear\\_impl\(\)](#) method to reset its internal state. The [clear\\_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

Reimplemented in [Flim](#).

## 8.19.2.2 bool FlimAbstract::isAcquiring ( ) const [inline]

tells if the data aquisition has finished reaching finish\_after\_outputframe

This function returns a boolean which tells the user if the class is still acquiring data. It can only reach the false state for finish\_after\_outputframe > 0.

## Note

This can differ from isRunning. The return value of isRunning state depends only on start/startFor/stop.

## 8.19.2.3 bool FlimAbstract::next\_impl ( std::vector&lt; Tag &gt; &amp; incoming\_tags, timestamp\_t begin\_time, timestamp\_t end\_time ) [override], [protected], [virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

## Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

## Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

**8.19.2.4** `virtual void FlimAbstract::on_frame_end ( )` `[protected]`, `[pure virtual]`

Implemented in [Flim](#), and [FlimBase](#).

**8.19.2.5** `void FlimAbstract::on_start ( )` `[override]`, `[protected]`, `[virtual]`

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

**8.19.2.6** `template<FastBinning::Mode bin_mode> void FlimAbstract::process_tags ( const std::vector< Tag > & incoming_tags )` `[protected]`

### 8.19.3 Member Data Documentation

**8.19.3.1** `bool FlimAbstract::acquiring {}` `[protected]`

**8.19.3.2** `std::recursive_mutex FlimAbstract::acquisition_lock` `[protected]`

**8.19.3.3** `FastBinning FlimAbstract::binner` `[protected]`

**8.19.3.4** `const timestamp_t FlimAbstract::binwidth` `[protected]`

**8.19.3.5** `const channel_t FlimAbstract::click_channel` `[protected]`

**8.19.3.6** `timestamp_t FlimAbstract::current_frame_begin` `[protected]`

**8.19.3.7** `timestamp_t FlimAbstract::current_frame_end` `[protected]`

**8.19.3.8** `size_t FlimAbstract::data_base {}` `[protected]`

- 8.19.3.9 `const uint32_t FlimAbstract::finish_after_outputframe` [protected]
- 8.19.3.10 `std::vector<uint32_t> FlimAbstract::frame` [protected]
- 8.19.3.11 `bool FlimAbstract::frame_acquisition {}` [protected]
- 8.19.3.12 `const channel_t FlimAbstract::frame_begin_channel` [protected]
- 8.19.3.13 `uint32_t FlimAbstract::frames_completed {}` [protected]
- 8.19.3.14 `bool FlimAbstract::initialized` [protected]
- 8.19.3.15 `const uint32_t FlimAbstract::n_bins` [protected]
- 8.19.3.16 `const uint32_t FlimAbstract::n_frame_average` [protected]
- 8.19.3.17 `const uint32_t FlimAbstract::n_pixels` [protected]
- 8.19.3.18 `bool FlimAbstract::pixel_acquisition {}` [protected]
- 8.19.3.19 `const channel_t FlimAbstract::pixel_begin_channel` [protected]
- 8.19.3.20 `std::vector<timestamp_t> FlimAbstract::pixel_begins` [protected]
- 8.19.3.21 `const channel_t FlimAbstract::pixel_end_channel` [protected]
- 8.19.3.22 `std::vector<timestamp_t> FlimAbstract::pixel_ends` [protected]
- 8.19.3.23 `uint32_t FlimAbstract::pixels_processed {}` [protected]
- 8.19.3.24 `std::deque<timestamp_t> FlimAbstract::previous_starts` [protected]
- 8.19.3.25 `const channel_t FlimAbstract::start_channel` [protected]
- 8.19.3.26 `uint32_t FlimAbstract::ticks {}` [protected]
- 8.19.3.27 `const timestamp_t FlimAbstract::time_window` [protected]


The documentation for this class was generated from the following file:

- [Iterators.h](#)

## 8.20 FlimBase Class Reference

```
#include <Iterators.h>
```

Inheritance diagram for FlimBase:



### Public Member Functions

- `FlimBase (TimeTaggerBase *tagger, channel_t start_channel, channel_t click_channel, channel_t pixel_begin_channel, uint32_t n_pixels, uint32_t n_bins, timestamp_t binwidth, channel_t pixel_end_channel=CHANNEL_UNUSED, channel_t frame_begin_channel=CHANNEL_UNUSED, uint32_t finish_after_outputframe=0, uint32_t n_frame_average=1, bool pre_initialize=true)`  
*construct a basic [Flim](#) measurement, containing a minimum featureset for efficiently purposes*
- `~FlimBase ()`
- `void initialize ()`  
*initializes and starts measuring this [Flim](#) measurement*

### Protected Member Functions

- `void on_frame_end ()` override
- `virtual void frameReady (uint32_t frame_number, std::vector< uint32_t > &data, std::vector< timestamp_t > &pixel_begin_times, std::vector< timestamp_t > &pixel_end_times, timestamp_t frame_begin_time, timestamp_t frame_end_time)`

### Protected Attributes

- `uint32_t total_frames`

#### 8.20.1 Constructor & Destructor Documentation

- 8.20.1.1 `FlimBase::FlimBase ( TimeTaggerBase * tagger, channel_t start_channel, channel_t click_channel, channel_t pixel_begin_channel, uint32_t n_pixels, uint32_t n_bins, timestamp_t binwidth, channel_t pixel_end_channel = CHANNEL_UNUSED, channel_t frame_begin_channel = CHANNEL_UNUSED, uint32_t finish_after_outputframe = 0, uint32_t n_frame_average = 1, bool pre_initialize = true )`

construct a basic [Flim](#) measurement, containing a minimum featureset for efficiently purposes

## Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>start_channel</i>	channel on which start clicks are received for the time differences histogramming
<i>click_channel</i>	channel on which clicks are received for the time differences histogramming
<i>pixel_begin_channel</i>	start of a pixel (histogram)
<i>n_pixels</i>	number of pixels (histograms) of one frame
<i>n_bins</i>	number of histogram bins for each pixel
<i>binwidth</i>	bin size in picoseconds
<i>pixel_end_channel</i>	end marker of a pixel - incoming clicks on the click_channel will be ignored afterwards
<i>frame_begin_channel</i>	(optional) start the frame, or reset the pixel index
<i>finish_after_outputframe</i>	(optional) sets the number of frames stored within the measurement class. After reaching the number, the measurement will stop. If the number is 0 (default value), one frame is stored and the measurement runs continuously.
<i>n_frame_average</i>	(optional) average multiple input frames into one output frame, default: 1
<i>pre_initialize</i>	(optional) initializes the measurement on constructing.

8.20.1.2 `FlimBase::~~FlimBase ( )`

## 8.20.2 Member Function Documentation

8.20.2.1 `virtual void FlimBase::frameReady ( uint32_t frame_number, std::vector< uint32_t > & data, std::vector< timestamp_t > & pixel_begin_times, std::vector< timestamp_t > & pixel_end_times, timestamp_t frame_begin_time, timestamp_t frame_end_time )` `[protected]`, `[virtual]`

8.20.2.2 `void FlimBase::initialize ( )`

initializes and starts measuring this [Flim](#) measurement

This function initializes the [Flim](#) measurement and starts executing it. It does nothing if preinitialized in the constructor is set to true.

8.20.2.3 `void FlimBase::on_frame_end ( )` `[override]`, `[protected]`, `[virtual]`

Implements [FlimAbstract](#).

## 8.20.3 Member Data Documentation

8.20.3.1 `uint32_t FlimBase::total_frames` `[protected]`

The documentation for this class was generated from the following file:

- [Iterators.h](#)



## 8.21 FlimFrameInfo Class Reference

```
#include <Iterators.h>
```

### Public Member Functions

- `int32_t` [getFrameNumber](#) () const  
*index of this frame*
- `bool` [isValid](#) () const  
*tells if this frame is valid*
- `uint32_t` [getPixelPosition](#) () const  
*number of pixels aquired on this frame*
- `void` [getHistograms](#) (std::function< `uint32_t` \*(`size_t`, `size_t`)> `array_out`)
- `void` [getIntensities](#) (std::function< `float` \*(`size_t`)> `array_out`)
- `void` [getSummedCounts](#) (std::function< `uint64_t` \*(`size_t`)> `array_out`)
- `void` [getPixelBegins](#) (std::function< `long long` \*(`size_t`)> `array_out`)
- `void` [getPixelEnds](#) (std::function< `long long` \*(`size_t`)> `array_out`)

### Public Attributes

- `uint32_t` [pixels](#)
- `uint32_t` [bins](#)
- `int32_t` [frame\\_number](#)
- `uint32_t` [pixel\\_position](#)

### Friends

- class [Flim](#)

#### 8.21.1 Member Function Documentation

8.21.1.1 `int32_t` [FlimFrameInfo::getFrameNumber](#) ( ) const `[inline]`

index of this frame

This function returns the frame number, starting from 0 for the very first frame acquired. If the index is -1, it is an invalid frame which is returned on error.

8.21.1.2 `void` [FlimFrameInfo::getHistograms](#) ( std::function< `uint32_t` \*(`size_t`, `size_t`)> `array_out` )

8.21.1.3 `void` [FlimFrameInfo::getIntensities](#) ( std::function< `float` \*(`size_t`)> `array_out` )

8.21.1.4 `void` [FlimFrameInfo::getPixelBegins](#) ( std::function< `long long` \*(`size_t`)> `array_out` )

8.21.1.5 `void` [FlimFrameInfo::getPixelEnds](#) ( std::function< `long long` \*(`size_t`)> `array_out` )

8.21.1.6 `uint32_t` [FlimFrameInfo::getPixelPosition](#) ( ) const `[inline]`

number of pixels aquired on this frame

This function returns a value which tells how many pixels were processed for this frame.

8.21.1.7 `void FlimFrameInfo::getSummedCounts ( std::function< uint64_t *(size_t)> array_out )`

8.21.1.8 `bool FlimFrameInfo::isValid ( ) const [inline]`

tells if this frame is valid

This function returns a boolean which tells if this frame is valid or not. Invalid frames are possible on errors, such as asking for the last completed frame when no frame has been completed so far.

## 8.21.2 Friends And Related Function Documentation

8.21.2.1 `friend class Flim [friend]`

## 8.21.3 Member Data Documentation

8.21.3.1 `uint32_t FlimFrameInfo::bins`

8.21.3.2 `int32_t FlimFrameInfo::frame_number`

8.21.3.3 `uint32_t FlimFrameInfo::pixel_position`

8.21.3.4 `uint32_t FlimFrameInfo::pixels`

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.22 FrequencyMultiplier Class Reference

The signal of an input channel is scaled up to a higher frequency according to the multiplier passed as a parameter.

```
#include <Iterators.h>
```

Inheritance diagram for FrequencyMultiplier:



## Public Member Functions

- [FrequencyMultiplier](#) ([TimeTaggerBase](#) \**tagger*, [channel\\_t](#) *input\_channel*, [int32\\_t](#) *multiplier*)  
*constructor of a [FrequencyMultiplier](#)*
- [~FrequencyMultiplier](#) ()
- [channel\\_t](#) [getChannel](#) ()
- [int32\\_t](#) [getMultiplier](#) ()

## Protected Member Functions

- [bool](#) [next\\_impl](#) ([std::vector](#)< [Tag](#) > &*incoming\_tags*, [timestamp\\_t](#) *begin\_time*, [timestamp\\_t](#) *end\_time*) override  
*update iterator state*

## Friends

- class [FrequencyMultiplierImpl](#)

## Additional Inherited Members

### 8.22.1 Detailed Description

The signal of an input channel is scaled up to a higher frequency according to the multiplier passed as a parameter.

The [FrequencyMultiplier](#) inserts copies the original input events from the *input\_channel* and adds additional events to match the upscaling factor. The algorithm used assumes a constant frequency and calculates out of the last two incoming events linearly the intermediate timestamps to match the upscaled frequency given by the multiplier parameter.

The [FrequencyMultiplier](#) can be used to restore the actual frequency applied to an *input\_channel* which was reduces via the [EventDivider](#) to lower the effective data rate. For example a 80 MHz laser sync signal can be scaled down via [setEventDivider](#)(..., 80) to 1 MHz (hardware side) and an 80 MHz signal can be restored via [FrequencyMultiplier](#)(..., 80) on the software side with some loss in precision. The [FrequencyMultiplier](#) is an alternative way to reduce the data rate in comparison to the [EventFilter](#), which has a higher precision but can be more difficult to use.

### 8.22.2 Constructor & Destructor Documentation

#### 8.22.2.1 [FrequencyMultiplier::FrequencyMultiplier](#) ( [TimeTaggerBase](#) \* *tagger*, [channel\\_t](#) *input\_channel*, [int32\\_t](#) *multiplier* )

constructor of a [FrequencyMultiplier](#)

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>input_channel</i>	channel on which the upscaling of the frequency is based on
<i>multiplier</i>	frequency upscaling factor

#### 8.22.2.2 FrequencyMultiplier::~~FrequencyMultiplier ( )

### 8.22.3 Member Function Documentation

#### 8.22.3.1 channel\_t FrequencyMultiplier::getChannel ( )

#### 8.22.3.2 int32\_t FrequencyMultiplier::getMultiplier ( )

#### 8.22.3.3 bool FrequencyMultiplier::next\_impl ( std::vector< Tag > & incoming\_tags, timestamp\_t begin\_time, timestamp\_t end\_time ) [override], [protected], [virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

#### Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

#### Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

### 8.22.4 Friends And Related Function Documentation

#### 8.22.4.1 friend class FrequencyMultiplierImpl [friend]

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.23 GatedChannel Class Reference

An input channel is gated by a gate channel.

```
#include <Iterators.h>
```

Inheritance diagram for GatedChannel:



### Public Member Functions

- `GatedChannel (TimeTaggerBase *tagger, channel_t input_channel, channel_t gate_start_channel, channel_t gate_stop_channel)`  
*constructor of a [GatedChannel](#)*
- `~GatedChannel ()`
- `channel_t getChannel ()`  
*the new virtual channel*

### Protected Member Functions

- `bool next_impl (std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time)` override  
*update iterator state*

### Friends

- class `GatedChannelImpl`

### Additional Inherited Members

#### 8.23.1 Detailed Description

An input channel is gated by a gate channel.

Note: The gate is edge sensitive and not level sensitive. That means that the gate will transfer data only when an appropriate level change is detected on the `gate_start_channel`.

#### 8.23.2 Constructor & Destructor Documentation

8.23.2.1 `GatedChannel::GatedChannel ( TimeTaggerBase * tagger, channel_t input_channel, channel_t gate_start_channel, channel_t gate_stop_channel )`

constructor of a [GatedChannel](#)

## Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>input_channel</i>	channel which is gated
<i>gate_start_channel</i>	channel on which a signal detected will start the transmission of the input_channel through the gate
<i>gate_stop_channel</i>	channel on which a signal detected will stop the transmission of the input_channel through the gate

8.23.2.2 `GatedChannel::~GatedChannel ( )`

## 8.23.3 Member Function Documentation

8.23.3.1 `channel_t GatedChannel::getChannel ( )`

the new virtual channel

This function returns the new allocated virtual channel. It can be used now in any new iterator.

8.23.3.2 `bool GatedChannel::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time )` `[override]`, `[protected]`, `[virtual]`

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

## Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

## Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

## 8.23.4 Friends And Related Function Documentation

8.23.4.1 `friend class GatedChannelImpl` `[friend]`

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.24 Histogram Class Reference

Accumulate time differences into a histogram.

```
#include <Iterators.h>
```

Inheritance diagram for Histogram:



### Public Member Functions

- `Histogram` (`TimeTaggerBase` \*`tagger`, `channel_t` `click_channel`, `channel_t` `start_channel`=`CHANNEL_UNU`↵  
`SED`, `timestamp_t` `binwidth`=1000, `int32_t` `n_bins`=1000)  
*constructor of a `Histogram` measurement*
- `~Histogram` ()
- void `getData` (`std::function`< `int32_t` \*(`size_t`)> `array_out`)
- void `getIndex` (`std::function`< `long long` \*(`size_t`)> `array_out`)

### Protected Member Functions

- bool `next_impl` (`std::vector`< `Tag` > &`incoming_tags`, `timestamp_t` `begin_time`, `timestamp_t` `end_time`) override  
*update iterator state*
- void `clear_impl` () override  
*clear `Iterator` state.*
- void `on_start` () override  
*callback when the measurement class is started*

### Friends

- class `TimeDifferencesImpl`< `Histogram` >

## Additional Inherited Members

### 8.24.1 Detailed Description

Accumulate time differences into a histogram.

This is a simple multiple start, multiple stop measurement. This is a special case of the more general '[Time↔Differences](#)' measurement. Specifically, the thread waits for clicks on a first channel, the 'start channel', then measures the time difference between the last start click and all subsequent clicks on a second channel, the 'click channel', and stores them in a histogram. The histogram range and resolution is specified by the number of bins and the binwidth. Clicks that fall outside the histogram range are ignored. Data accumulation is performed independently for all start clicks. This type of measurement is frequently referred to as 'multiple start, multiple stop' measurement and corresponds to a full auto- or cross-correlation measurement.

### 8.24.2 Constructor & Destructor Documentation

**8.24.2.1** `Histogram::Histogram ( TimeTaggerBase * tagger, channel_t click_channel, channel_t start_channel = CHANNEL_UNUSED, timestamp_t binwidth = 1000, int32_t n_bins = 1000 )`

constructor of a [Histogram](#) measurement

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>click_channel</i>	channel that increments the count in a bin
<i>start_channel</i>	channel that sets start times relative to which clicks on the click channel are measured
<i>binwidth</i>	width of one histogram bin in ps
<i>n_bins</i>	number of bins in the histogram

**8.24.2.2** `Histogram::~~Histogram ( )`

### 8.24.3 Member Function Documentation

**8.24.3.1** `void Histogram::clear_impl ( )` `[override]`, `[protected]`, `[virtual]`

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented from [IteratorBase](#).



8.24.3.2 void Histogram::getData ( std::function< int32\_t \*(size\_t)> *array\_out* )

8.24.3.3 void Histogram::getIndex ( std::function< long long \*(size\_t)> *array\_out* )

8.24.3.4 bool Histogram::next\_impl ( std::vector< Tag > & *incoming\_tags*, timestamp\_t *begin\_time*, timestamp\_t *end\_time* ) [override],[protected],[virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

**Parameters**

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

**Returns**

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

**8.24.3.5** `void Histogram::on_start ( )` `[override]`, `[protected]`, `[virtual]`

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

**8.24.4 Friends And Related Function Documentation**

**8.24.4.1** `friend class TimeDifferencesImpl< Histogram >` `[friend]`

The documentation for this class was generated from the following file:


- [Iterators.h](#)

**8.25 Histogram2D Class Reference**

A 2-dimensional histogram of time differences. This can be used in measurements similar to 2D NRM spectroscopy.

```
#include <Iterators.h>
```

Inheritance diagram for Histogram2D:



## Public Member Functions

- [Histogram2D](#) ([TimeTaggerBase](#) \**tagger*, [channel\\_t](#) *start\_channel*, [channel\\_t](#) *stop\_channel\_1*, [channel\\_t](#) *stop\_channel\_2*, [timestamp\\_t](#) *binwidth\_1*, [timestamp\\_t](#) *binwidth\_2*, [int32\\_t](#) *n\_bins\_1*, [int32\\_t](#) *n\_bins\_2*)  
*constructor of a [Histogram2D](#) measurement*
- [~Histogram2D](#) ()
- void [getData](#) (std::function< [int32\\_t](#) \*([size\\_t](#), [size\\_t](#))> *array\_out*)
- void [getIndex](#) (std::function< long long \*([size\\_t](#), [size\\_t](#), [size\\_t](#))> *array\_out*)
- void [getIndex\\_1](#) (std::function< long long \*([size\\_t](#))> *array\_out*)
- void [getIndex\\_2](#) (std::function< long long \*([size\\_t](#))> *array\_out*)

## Protected Member Functions

- bool [next\\_impl](#) (std::vector< [Tag](#) > &*incoming\_tags*, [timestamp\\_t](#) *begin\_time*, [timestamp\\_t](#) *end\_time*) override  
*update iterator state*
- void [clear\\_impl](#) () override  
*clear [Iterator](#) state.*

## Friends

- class [Histogram2DImpl](#)

## Additional Inherited Members

### 8.25.1 Detailed Description

A 2-dimensional histogram of time differences. This can be used in measurements similar to 2D NRM spectroscopy.

This measurement is a 2-dimensional version of the [Histogram](#) measurement. The measurement accumulates two-dimensional histogram where stop signals from two separate channels define the bin coordinate. For instance, this kind of measurement is similar to that of typical 2D NMR spectroscopy.

### 8.25.2 Constructor & Destructor Documentation

- 8.25.2.1 [Histogram2D::Histogram2D](#) ( [TimeTaggerBase](#) \* *tagger*, [channel\\_t](#) *start\_channel*, [channel\\_t](#) *stop\_channel\_1*, [channel\\_t](#) *stop\_channel\_2*, [timestamp\\_t](#) *binwidth\_1*, [timestamp\\_t](#) *binwidth\_2*, [int32\\_t](#) *n\_bins\_1*, [int32\\_t](#) *n\_bins\_2* )

constructor of a [Histogram2D](#) measurement

#### Parameters

<i>tagger</i>	time tagger object
<i>start_channel</i>	channel on which start clicks are received
<i>stop_channel_1</i>	channel on which stop clicks for the time axis 1 are received
<i>stop_channel_2</i>	channel on which stop clicks for the time axis 2 are received
<i>binwidth_1</i>	bin width in ps for the time axis 1
<i>binwidth_2</i>	bin width in ps for the time axis 2
<i>n_bins_1</i>	the number of bins along the time axis 1
<i>n_bins_2</i>	the number of bins along the time axis 2

### 8.25.2.2 Histogram2D::~Histogram2D ( )

## 8.25.3 Member Function Documentation

### 8.25.3.1 void Histogram2D::clear\_impl ( ) [override],[protected],[virtual]

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear\\_impl\(\)](#) method to reset its internal state. The [clear\\_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

### 8.25.3.2 void Histogram2D::getData ( std::function< int32\_t \*(size\_t, size\_t)> array\_out )

Returns a two-dimensional array of size `n_bins_1` by `n_bins_2` containing the 2D histogram.

### 8.25.3.3 void Histogram2D::getIndex ( std::function< long long \*(size\_t, size\_t, size\_t)> array\_out )

Returns a 3D array containing two coordinate matrices (meshgrid) for time bins in ps for the time axes 1 and 2. For details on meshgrid please take a look at the respective documentation either for Matlab or Python NumPy

### 8.25.3.4 void Histogram2D::getIndex\_1 ( std::function< long long \*(size\_t)> array\_out )

Returns a vector of size `n_bins_1` containing the bin locations in ps for the time axis 1.

### 8.25.3.5 void Histogram2D::getIndex\_2 ( std::function< long long \*(size\_t)> array\_out )

Returns a vector of size `n_bins_2` containing the bin locations in ps for the time axis 2.

### 8.25.3.6 bool Histogram2D::next\_impl ( std::vector< Tag > &incoming\_tags, timestamp\_t begin\_time, timestamp\_t end\_time ) [override],[protected],[virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

#### Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

**Returns**

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

**8.25.4 Friends And Related Function Documentation****8.25.4.1 friend class Histogram2DImpl** [*friend*]

The documentation for this class was generated from the following file:


- [Iterators.h](#)

**8.26 HistogramLogBins Class Reference**

Accumulate time differences into a histogram with logarithmic increasing bin sizes.

```
#include <Iterators.h>
```

Inheritance diagram for HistogramLogBins:

**Public Member Functions**

- [HistogramLogBins](#) ([TimeTaggerBase](#) \*tagger, [channel\\_t](#) click\_channel, [channel\\_t](#) start\_channel, double exp\_start, double exp\_stop, int32\_t n\_bins)  
*constructor of a [HistogramLogBins](#) measurement*
- [~HistogramLogBins](#) ()
- void [getData](#) (std::function< uint64\_t \*(size\_t)> array\_out)
- void [getDataNormalizedCountsPerPs](#) (std::function< double \*(size\_t)> array\_out)
- void [getDataNormalizedG2](#) (std::function< double \*(size\_t)> array\_out)
- void [getBinEdges](#) (std::function< long long \*(size\_t)> array\_out)

## Protected Member Functions

- bool `next_impl` (std::vector< [Tag](#) > &incoming\_tags, [timestamp\\_t](#) begin\_time, [timestamp\\_t](#) end\_time) override  
*update iterator state*
- void `clear_impl` () override  
*clear [Iterator](#) state.*

## Friends

- class [HistogramLogBinsImpl](#)

## Additional Inherited Members

### 8.26.1 Detailed Description

Accumulate time differences into a histogram with logarithmic increasing bin sizes.

This is a multiple start, multiple stop measurement, and works the very same way as the histogram measurement but with logarithmic increasing bin widths. After initializing the measurement (or after an overflow) no data is accumulated in the histogram until the full histogram duration has passed to ensure a balanced count accumulation over the full histogram.

### 8.26.2 Constructor & Destructor Documentation

**8.26.2.1** `HistogramLogBins::HistogramLogBins ( TimeTaggerBase * tagger, channel\_t click_channel, channel\_t start_channel, double exp_start, double exp_stop, int32\_t n_bins )`

constructor of a [HistogramLogBins](#) measurement

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>click_channel</i>	channel that increments the count in a bin
<i>start_channel</i>	channel that sets start times relative to which clicks on the click channel are measured
<i>exp_start</i>	exponent for the lowest time differences in the histogram: $10^{\text{exp\_start}}$ s, lowest <i>exp_start</i> : -12 => 1ps
<i>exp_stop</i>	exponent for the highest time differences in the histogram: $10^{\text{exp\_stop}}$ s
<i>n_bins</i>	total number of bins in the histogram

**8.26.2.2** `HistogramLogBins::~~HistogramLogBins ( )`

### 8.26.3 Member Function Documentation

**8.26.3.1** `void HistogramLogBins::clear_impl ( )` [[override](#)], [[protected](#)], [[virtual](#)]

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear\\_impl\(\)](#) method to reset its internal state. The [clear\\_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

**8.26.3.2** `void HistogramLogBins::getBinEdges ( std::function< long long *(size_t)> array_out )`

returns the edges of the bins in ps

**8.26.3.3** `void HistogramLogBins::getData ( std::function< uint64_t *(size_t)> array_out )`

returns the absolute counts for the bins

**8.26.3.4** `void HistogramLogBins::getDataNormalizedCountsPerPs ( std::function< double *(size_t)> array_out )`

returns the counts normalized by the binwidth of each bin

**8.26.3.5** `void HistogramLogBins::getDataNormalizedG2 ( std::function< double *(size_t)> array_out )`

returns the counts normalized by the binwidth and the average count rate. This matches the implementation of [Correlation::getDataNormalized](#)

**8.26.3.6** `bool HistogramLogBins::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time ) [override], [protected], [virtual]`

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

#### Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

#### Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

## 8.26.4 Friends And Related Function Documentation

#### 8.26.4.1 friend class HistogramLogBinsImpl [friend]

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.27 Iterator Class Reference

a simple event queue

```
#include <Iterators.h>
```

Inheritance diagram for Iterator:



### Public Member Functions

- [Iterator](#) ([TimeTaggerBase](#) \*tagger, [channel\\_t](#) channel)  
*standard constructor*
- [~Iterator](#) ()
- [timestamp\\_t next](#) ()  
*get next timestamp*
- [uint64\\_t size](#) ()  
*get queue size*

### Protected Member Functions

- [bool next\\_impl](#) (std::vector< [Tag](#) > &incoming\_tags, [timestamp\\_t](#) begin\_time, [timestamp\\_t](#) end\_time) override  
*update iterator state*
- [void clear\\_impl](#) () override  
*clear [Iterator](#) state.*

### Friends

- class [IteratorImpl](#)



## Additional Inherited Members

### 8.27.1 Detailed Description

a simple event queue

A simple [Iterator](#), just keeping a first-in first-out queue of event timestamps.

**Deprecated** use [TimeTagStream](#)

### 8.27.2 Constructor & Destructor Documentation

#### 8.27.2.1 `Iterator::Iterator ( TimeTaggerBase * tagger, channel_t channel )`

standard constructor

Parameters

<i>tagger</i>	the backend
<i>channel</i>	the channel to get events from

#### 8.27.2.2 `Iterator::~~Iterator ( )`

### 8.27.3 Member Function Documentation

#### 8.27.3.1 `void Iterator::clear_impl ( ) [override], [protected], [virtual]`

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear\\_impl\(\)](#) method to reset its internal state. The [clear\\_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

#### 8.27.3.2 `timestamp_t Iterator::next ( )`

get next timestamp

get the next timestamp from the queue.

#### 8.27.3.3 `bool Iterator::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time ) [override], [protected], [virtual]`

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

**Parameters**

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

**Returns**

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

8.27.3.4 `uint64_t Iterator::size ( )`

get queue size

## 8.27.4 Friends And Related Function Documentation

8.27.4.1 `friend class IteratorImpl [friend]`

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.28 IteratorBase Class Reference

Base class for all iterators.

```
#include <TimeTagger.h>
```

Inheritance diagram for IteratorBase:



## Public Member Functions

- virtual `~IteratorBase()`  
*destructor*
- void `start()`  
*start the iterator*
- void `startFor(timestamp_t capture_duration, bool clear=true)`

- *start the iterator, and stops it after the capture\_duration*
- bool `waitUntilFinished` (int64\_t timeout=-1)  
*wait until the iterator has finished running.*
- void `stop` ()  
*stop the iterator*
- void `clear` ()  
*clear *Iterator* state.*
- bool `isRunning` ()  
*query the *Iterator* state.*
- `timestamp_t` `getCaptureDuration` ()  
*query the evaluation time*

## Protected Member Functions

- `IteratorBase` (`TimeTaggerBase` \*tagger, std::string base\_type\_="IteratorBase", std::string extra\_info\_="")  
*standard constructor*
- void `registerChannel` (`channel_t` channel)  
*register a channel*
- void `unregisterChannel` (`channel_t` channel)  
*unregister a channel*
- `channel_t` `getNewVirtualChannel` ()  
*allocate a new virtual output channel for this iterator*
- void `finishInitialization` ()  
*method to call after finishing the initialization of the measurement*
- virtual void `clear_impl` ()  
*clear *Iterator* state.*
- virtual void `on_start` ()  
*callback when the measurement class is started*
- virtual void `on_stop` ()  
*callback when the measurement class is stopped*
- void `lock` ()  
*acquire update lock*
- void `unlock` ()  
*release update lock*
- `OrderedBarrier::OrderInstance` `parallelize` (`OrderedPipeline` &pipeline)  
*release lock and continue work in parallel*
- std::unique\_lock< std::mutex > `getLock` ()  
*acquire update lock*
- virtual bool `next_impl` (std::vector< `Tag` > &incoming\_tags, `timestamp_t` begin\_time, `timestamp_t` end\_time)=0  
*update iterator state*
- void `finish_running` ()

## Protected Attributes

- std::set< `channel_t` > `channels_registered`  
*list of channels used by the iterator*
- bool `running`  
*running state of the iterator*
- bool `autostart`
- `TimeTaggerBase` \* `tagger`
- `timestamp_t` `capture_duration`

## Friends

- class [TimeTaggerRunner](#)
- class [TimeTaggerProxy](#)
- class [SynchronizedMeasurements](#)

### 8.28.1 Detailed Description

Base class for all iterators.

### 8.28.2 Constructor & Destructor Documentation

**8.28.2.1** `IteratorBase::IteratorBase ( TimeTaggerBase * tagger, std::string base_type_ = "IteratorBase", std::string extra_info_ = " " )` `[protected]`

standard constructor

will register with the [TimeTagger](#) backend.

**8.28.2.2** `virtual IteratorBase::~~IteratorBase ( )` `[virtual]`

destructor

will stop and unregister prior finalization.

### 8.28.3 Member Function Documentation

**8.28.3.1** `void IteratorBase::clear ( )`

clear [Iterator](#) state.

**8.28.3.2** `virtual void IteratorBase::clear_impl ( )` `[inline]`, `[protected]`, `[virtual]`

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented in [Flim](#), [FlimAbstract](#), [CustomMeasurementBase](#), [EventGenerator](#), [FileWriter](#), [Scope](#), [Correlation](#), [HistogramLogBins](#), [Histogram](#), [TimeDifferencesND](#), [Histogram2D](#), [TimeDifferences](#), [StartStop](#), [Dump](#), [TimeTagStream](#), [Iterator](#), [Countrate](#), [Counter](#), [CountBetweenMarkers](#), and [Combiner](#).

**8.28.3.3** `void IteratorBase::finish_running ( )` `[protected]`

**8.28.3.4** `void IteratorBase::finishInitialization ( )` `[protected]`

method to call after finishing the initialization of the measurement

#### 8.28.3.5 `timestamp_t` `IteratorBase::getCaptureDuration ( )`

query the evaluation time

Query the total capture duration since the last call to `clear`. This might have a wrong amount of time if there were some overflows within this range.

##### Returns

capture duration of the data

#### 8.28.3.6 `std::unique_lock<std::mutex>` `IteratorBase::getLock ( )` `[protected]`

acquire update lock

All mutable operations on a iterator are guarded with an update mutex. Implementers are advised to lock an iterator, whenever internal state is queried or changed.

##### Returns

a lock object, which releases the lock when this instance is freed

#### 8.28.3.7 `channel_t` `IteratorBase::getNewVirtualChannel ( )` `[protected]`

allocate a new virtual output channel for this iterator

#### 8.28.3.8 `bool` `IteratorBase::isRunning ( )`

query the [Iterator](#) state.

Fetches if this iterator is running.

#### 8.28.3.9 `void` `IteratorBase::lock ( )` `[protected]`

acquire update lock

All mutable operations on a iterator are guarded with an update mutex. Implementers are advised to [lock\(\)](#) an iterator, whenever internal state is queried or changed.

**Deprecated** use `getLock`

#### 8.28.3.10 `virtual bool` `IteratorBase::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time )` `[protected]`, `[pure virtual]`

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

## Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

## Returns

true if the content of this block was modified, false otherwise

Implemented in [SyntheticSingleTag](#), [FlimAbstract](#), [CustomMeasurementBase](#), [EventGenerator](#), [FileWriter](#), [ConstantFractionDiscriminator](#), [Scope](#), [Correlation](#), [HistogramLogBins](#), [Histogram](#), [TimeDifferencesND](#), [Histogram2D](#), [TimeDifferences](#), [StartStop](#), [Dump](#), [TimeTagStream](#), [Iterator](#), [FrequencyMultiplier](#), [GatedChannel](#), [TriggerOnCountrate](#), [DelayedChannel](#), [Countrate](#), [Coincidences](#), [Counter](#), [CountBetweenMarkers](#), and [Combiner](#).

**8.28.3.11** `virtual void IteratorBase::on_start ( ) [inline],[protected],[virtual]`

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented in [FlimAbstract](#), [CustomMeasurementBase](#), [EventGenerator](#), [FileWriter](#), [ConstantFractionDiscriminator](#), [Histogram](#), [TimeDifferencesND](#), [TimeDifferences](#), [StartStop](#), [Dump](#), [TriggerOnCountrate](#), [DelayedChannel](#), [Countrate](#), and [Counter](#).

**8.28.3.12** `virtual void IteratorBase::on_stop ( ) [inline],[protected],[virtual]`

callback when the measurement class is stopped

This function is guarded by the update lock.

Reimplemented in [CustomMeasurementBase](#), [FileWriter](#), and [Dump](#).

**8.28.3.13** `OrderedBarrier::OrderInstance IteratorBase::parallelize ( OrderedPipeline & pipeline ) [protected]`

release lock and continue work in parallel

The measurement's lock is released, allowing this measurement to continue, while still executing work in parallel.

## Returns

a ordered barrier instance that can be synced afterwards.

**8.28.3.14** `void IteratorBase::registerChannel ( channel_t channel ) [protected]`

register a channel

Only channels registered by any iterator attached to a backend are delivered over the usb.

## Parameters

<i>channel</i>	the channel
----------------	-------------

## 8.28.3.15 void IteratorBase::start ( )

start the iterator

The default behavior for iterators is to start automatically on creation.

## 8.28.3.16 void IteratorBase::startFor ( timestamp\_t capture\_duration, bool clear = true )

start the iterator, and stops it after the capture\_duration

## Parameters

<i>capture_duration</i>	capture duration until the measurement is stopped
<i>clear</i>	resets the data aquired

When the startFor is called before the previous measurement has ended and the clear parameter is set to false, then the passed capture\_duration will be added on top to the current max\_capture\_duration

## 8.28.3.17 void IteratorBase::stop ( )

stop the iterator

The iterator is put into the STOPPED state, but will still be registered with the backend.

## 8.28.3.18 void IteratorBase::unlock ( ) [protected]

release update lock

see [lock\(\)](#)

**Deprecated** use getLock

## 8.28.3.19 void IteratorBase::unregisterChannel ( channel\_t channel ) [protected]

unregister a channel

## Parameters

<i>channel</i>	the channel
----------------	-------------



8.28.3.20 `bool IteratorBase::waitUntilFinished ( int64_t timeout = -1 )`

wait until the iterator has finished running.

#### Parameters

<i>timeout</i>	time in milliseconds to wait for the measurements. If negative, wait until finished.
----------------	--

`waitUntilFinished` will wait according to the timeout and return true if the iterator finished or false if not. Furthermore, when `waitUntilFinished` is called on a iterator running indefinitely, it will log an error and return immediately.

## 8.28.4 Friends And Related Function Documentation

8.28.4.1 `friend class SynchronizedMeasurements` [friend]

8.28.4.2 `friend class TimeTaggerProxy` [friend]

8.28.4.3 `friend class TimeTaggerRunner` [friend]

## 8.28.5 Member Data Documentation

8.28.5.1 `bool IteratorBase::autostart` [protected]

8.28.5.2 `timestamp_t IteratorBase::capture_duration` [protected]

8.28.5.3 `std::set<channel_t> IteratorBase::channels_registered` [protected]

list of channels used by the iterator

8.28.5.4 `bool IteratorBase::running` [protected]

running state of the iterator

8.28.5.5 `TimeTaggerBase* IteratorBase::tagger` [protected]

The documentation for this class was generated from the following file:

- [TimeTagger.h](#)

## 8.29 OrderedBarrier Class Reference

```
#include <TimeTagger.h>
```

## Classes

- class [OrderInstance](#)

## Public Member Functions

- [OrderedBarrier](#) ()
- [~OrderedBarrier](#) ()
- [OrderInstance](#) queue ()
- void [waitUntilFinished](#) ()

## Friends

- class [OrderInstance](#)

## 8.29.1 Constructor & Destructor Documentation

8.29.1.1 `OrderedBarrier::OrderedBarrier ( )`

8.29.1.2 `OrderedBarrier::~~OrderedBarrier ( )`

## 8.29.2 Member Function Documentation

8.29.2.1 `OrderInstance OrderedBarrier::queue ( )`

8.29.2.2 `void OrderedBarrier::waitUntilFinished ( )`

## 8.29.3 Friends And Related Function Documentation

8.29.3.1 `friend class OrderInstance` [`friend`]

The documentation for this class was generated from the following file:

- [TimeTagger.h](#)

## 8.30 OrderedPipeline Class Reference

```
#include <TimeTagger.h>
```

## Public Member Functions

- [OrderedPipeline](#) ()
- [~OrderedPipeline](#) ()

## Friends

- class [IteratorBase](#)

### 8.30.1 Constructor & Destructor Documentation

8.30.1.1 `OrderedPipeline::OrderedPipeline ( )`

8.30.1.2 `OrderedPipeline::~~OrderedPipeline ( )`

### 8.30.2 Friends And Related Function Documentation

8.30.2.1 `friend class IteratorBase [friend]`

The documentation for this class was generated from the following file:

- [TimeTagger.h](#)

## 8.31 OrderedBarrier::OrderInstance Class Reference

```
#include <TimeTagger.h>
```

## Public Member Functions

- [OrderInstance](#) ( )
- [OrderInstance](#) ([OrderedBarrier](#) \*parent, uint64\_t instance\_id)
- [~OrderInstance](#) ( )
- void [sync](#) ( )
- void [release](#) ( )

## Friends

- class [OrderedBarrier](#)

### 8.31.1 Constructor & Destructor Documentation

8.31.1.1 `OrderedBarrier::OrderInstance::OrderInstance ( )`

8.31.1.2 `OrderedBarrier::OrderInstance::OrderInstance ( OrderedBarrier * parent, uint64_t instance_id )`

8.31.1.3 `OrderedBarrier::OrderInstance::~~OrderInstance ( )`

### 8.31.2 Member Function Documentation

8.31.2.1 `void OrderedBarrier::OrderInstance::release ( )`

8.31.2.2 `void OrderedBarrier::OrderInstance::sync ( )`

### 8.31.3 Friends And Related Function Documentation

8.31.3.1 `friend class OrderedBarrier [friend]`


The documentation for this class was generated from the following file:

- [TimeTagger.h](#)

## 8.32 Scope Class Reference

```
#include <Iterators.h>
```

Inheritance diagram for Scope:



### Public Member Functions

- `Scope (TimeTaggerBase *tagger, std::vector< channel_t > event_channels, channel_t trigger_channel, timestamp_t window_size=1000000000, int32_t n_traces=1, int32_t n_max_events=1000)`  
*constructor of a [Scope](#) measurement*
- `~Scope ()`
- `bool ready ()`
- `int32_t triggered ()`
- `std::vector< std::vector< Event > > getData ()`
- `timestamp_t getWindowSize ()`

### Protected Member Functions

- `bool next_impl (std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time)` override  
*update iterator state*
- `void clear_impl ()` override  
*clear [Iterator](#) state.*

### Friends

- class [ScopeImpl](#)

### Additional Inherited Members

#### 8.32.1 Constructor & Destructor Documentation

8.32.1.1 `Scope::Scope ( TimeTaggerBase * tagger, std::vector< channel_t > event_channels, channel_t trigger_channel, timestamp_t window_size = 1000000000, int32_t n_traces = 1, int32_t n_max_events = 1000 )`

constructor of a [Scope](#) measurement

## Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>event_channels</i>	channels which are captured
<i>trigger_channel</i>	channel that starts a new trace
<i>window_size</i>	window time of each trace
<i>n_traces</i>	amount of traces ( $n\_traces < 1$ , automatic retrigger)
<i>n_max_events</i>	maximum number of tags in each trace

## 8.32.1.2 Scope::~~Scope ( )

## 8.32.2 Member Function Documentation

## 8.32.2.1 void Scope::clear\_impl ( ) [override],[protected],[virtual]

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear\\_impl\(\)](#) method to reset its internal state. The [clear\\_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

## 8.32.2.2 std::vector&lt;std::vector&lt;Event&gt; &gt; Scope::getData ( )

## 8.32.2.3 timestamp\_t Scope::getWindowSize ( )

## 8.32.2.4 bool Scope::next\_impl ( std::vector&lt; Tag &gt; &amp; incoming\_tags, timestamp\_t begin\_time, timestamp\_t end\_time ) [override],[protected],[virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

## Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

## Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

8.32.2.5 `bool Scope::ready ( )`

8.32.2.6 `int32_t Scope::triggered ( )`

### 8.32.3 Friends And Related Function Documentation

8.32.3.1 `friend class ScopeImpl [friend]`

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.33 StartStop Class Reference

simple start-stop measurement

```
#include <Iterators.h>
```

Inheritance diagram for StartStop:



### Public Member Functions

- [StartStop](#) ([TimeTaggerBase](#) \*tagger, [channel\\_t](#) click\_channel, [channel\\_t](#) start\_channel=CHANNEL\_UNU↵  
SED, [timestamp\\_t](#) binwidth=1000)  
*constructor of StartStop*
- [~StartStop](#) ()
- void [getData](#) (std::function< long long \*(size\_t, size\_t)> array\_out)

### Protected Member Functions

- bool [next\\_impl](#) (std::vector< [Tag](#) > &incoming\_tags, [timestamp\\_t](#) begin\_time, [timestamp\\_t](#) end\_time) over-  
ride  
*update iterator state*
- void [clear\\_impl](#) () override  
*clear Iterator state.*
- void [on\\_start](#) () override  
*callback when the measurement class is started*

## Friends

- class [StartStopImpl](#)

## Additional Inherited Members

### 8.33.1 Detailed Description

simple start-stop measurement

This class performs a start-stop measurement between two channels and stores the time differences in a histogram. The histogram resolution is specified beforehand (binwidth) but the histogram range is unlimited. It is adapted to the largest time difference that was detected. Thus all pairs of subsequent clicks are registered.

Be aware, on long-running measurements this may considerably slow down system performance and even crash the system entirely when attached to an unsuitable signal source.

### 8.33.2 Constructor & Destructor Documentation

**8.33.2.1** `StartStop::StartStop ( TimeTaggerBase * tagger, channel_t click_channel, channel_t start_channel = CHANNEL_UNUSED, timestamp_t binwidth = 1000 )`

constructor of [StartStop](#)

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>click_channel</i>	channel for stop clicks
<i>start_channel</i>	channel for start clicks
<i>binwidth</i>	width of one histogram bin in ps

**8.33.2.2** `StartStop::~StartStop ( )`

### 8.33.3 Member Function Documentation

**8.33.3.1** `void StartStop::clear_impl ( )` `[override]`, `[protected]`, `[virtual]`

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

8.33.3.2 `void StartStop::getData ( std::function< long long *(size_t, size_t)> array_out )`

8.33.3.3 `bool StartStop::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time )` `[override]`, `[protected]`, `[virtual]`

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

#### Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

#### Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

8.33.3.4 `void StartStop::on_start ( )` `[override]`, `[protected]`, `[virtual]`

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

## 8.33.4 Friends And Related Function Documentation

8.33.4.1 `friend class StartStopImpl` `[friend]`

The documentation for this class was generated from the following file:

- [Iterators.h](#)

## 8.34 SynchronizedMeasurements Class Reference

start, stop and clear several measurements synchronized

```
#include <Iterators.h>
```



## Public Member Functions

- [SynchronizedMeasurements](#) ([TimeTaggerBase](#) \*tagger)  
*construct a [SynchronizedMeasurements](#) object*
- [~SynchronizedMeasurements](#) ()
- void [registerMeasurement](#) ([IteratorBase](#) \*measurement)  
*register a measurement (iterator) to the SynchronizedMeasurements-group.*
- void [unregisterMeasurement](#) ([IteratorBase](#) \*measurement)  
*unregister a measurement (iterator) from the SynchronizedMeasurements-group.*
- void [clear](#) ()  
*clear all registered measurements synchronously*
- void [start](#) ()  
*start all registered measurements synchronously*
- void [stop](#) ()  
*stop all registered measurements synchronously*
- void [startFor](#) ([timestamp\\_t](#) capture\_duration, bool [clear](#)=true)  
*start all registered measurements synchronously, and stops them after the capture\_duration*
- bool [waitUntilFinished](#) ([int64\\_t](#) timeout=-1)  
*wait until all registered measurements have finished running.*
- bool [isRunning](#) ()  
*check if any iterator is running*
- [TimeTaggerBase](#) \* [getTagger](#) ()

## Protected Member Functions

- void [runCallback](#) ([TimeTaggerBase::IteratorCallback](#) callback, [std::unique\\_lock](#)< [std::mutex](#) > &lk, bool block=true)  
*run a callback on all registered measurements synchronously*

## Friends

- class [TimeTaggerProxy](#)

### 8.34.1 Detailed Description

start, stop and clear several measurements synchronized

For the case that several measurements should be started, stopped or cleared at the very same time, a [SynchronizedMeasurements](#) object can be create to which all the measurements (also called iterators) can be registered with [.registerMeasurement\(measurement\)](#). Calling [.stop\(\)](#), [.start\(\)](#) or [.clear\(\)](#) on the [SynchronizedMeasurements](#) object will call the respective method on each of the registered measurements at the very same time. That means that all measurements taking part will have processed the very same time tags.

### 8.34.2 Constructor & Destructor Documentation

#### 8.34.2.1 SynchronizedMeasurements::SynchronizedMeasurements ( [TimeTaggerBase](#) \* tagger )

construct a [SynchronizedMeasurements](#) object

## Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
---------------	---

8.34.2.2 `SynchronizedMeasurements::~~SynchronizedMeasurements ( )`

## 8.34.3 Member Function Documentation

8.34.3.1 `void SynchronizedMeasurements::clear ( )`

clear all registered measurements synchronously

8.34.3.2 `TimeTaggerBase* SynchronizedMeasurements::getTagger ( )`

Returns a proxy tagger object, which shall be used to create immediately registered measurements. Those measurements will not start automatically.

8.34.3.3 `bool SynchronizedMeasurements::isRunning ( )`

check if any iterator is running

8.34.3.4 `void SynchronizedMeasurements::registerMeasurement ( IteratorBase * measurement )`

register a measurement (iterator) to the SynchronizedMeasurements-group.

All available methods called on the [SynchronizedMeasurements](#) will happen at the very same time for all the registered measurements.

8.34.3.5 `void SynchronizedMeasurements::runCallback ( TimeTaggerBase::IteratorCallback callback, std::unique_lock< std::mutex > & lk, bool block = true ) [protected]`

run a callback on all registered measurements synchronously

Please keep in mind that the callback is copied for each measurement. So please avoid big captures.

8.34.3.6 `void SynchronizedMeasurements::start ( )`

start all registered measurements synchronously

8.34.3.7 `void SynchronizedMeasurements::startFor ( timestamp_t capture_duration, bool clear = true )`

start all registered measurements synchronously, and stops them after the *capture\_duration*

**8.34.3.8 void SynchronizedMeasurements::stop ( )**

stop all registered measurements synchronously

**8.34.3.9 void SynchronizedMeasurements::unregisterMeasurement ( IteratorBase \* *measurement* )**

unregister a measurement (iterator) from the SynchronizedMeasurements-group.

Stops synchronizing calls on the selected measurement, if the measurement is not within this synchronized group, the method does nothing.

**8.34.3.10 bool SynchronizedMeasurements::waitUntilFinished ( int64\_t *timeout* = -1 )**

wait until all registered measurements have finished running.

**Parameters**

<i>timeout</i>	time in milliseconds to wait for the measurements. If negative, wait until finished.
----------------	--

waitUntilFinished will wait according to the timeout and return true if all measurements finished or false if not. Furthermore, when waitUntilFinished is called on a set running indefinitely, it will log an error and return immediately.

**8.34.4 Friends And Related Function Documentation****8.34.4.1 friend class TimeTaggerProxy [friend]**

The documentation for this class was generated from the following file:


- [Iterators.h](#)

**8.35 SyntheticSingleTag Class Reference**

synthetic trigger timetag generator.

```
#include <Iterators.h>
```

Inheritance diagram for SyntheticSingleTag:



## Public Member Functions

- [SyntheticSingleTag](#) ([TimeTaggerBase](#) \**tagger*, [channel\\_t](#) *base\_channel*=[CHANNEL\\_UNUSED](#))  
*Construct a pulse event generator.*
- [~SyntheticSingleTag](#) ()
- void [trigger](#) ()  
*Generate a timetag for each call of this method.*
- [channel\\_t](#) [getChannel](#) () const

## Protected Member Functions

- bool [next\\_impl](#) (std::vector< [Tag](#) > &*incoming\_tags*, [timestamp\\_t](#) *begin\_time*, [timestamp\\_t](#) *end\_time*) override  
*update iterator state*

## Friends

- class [SyntheticSingleTagImpl](#)

## Additional Inherited Members

### 8.35.1 Detailed Description

synthetic trigger timetag generator.

Creates timetags based on a trigger method. Whenever the user calls the 'trigger' method, a timetag will be added to the *base\_channel*.

This synthetic channel can inject timetags into an existing channel or create a new virtual channel.

### 8.35.2 Constructor & Destructor Documentation

#### 8.35.2.1 [SyntheticSingleTag::SyntheticSingleTag](#) ( [TimeTaggerBase](#) \* *tagger*, [channel\\_t](#) *base\_channel* = [CHANNEL\\_UNUSED](#) )

Construct a pulse event generator.

##### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>base_channel</i>	base channel to which this signal will be added. If unused, a new channel will be created.

#### 8.35.2.2 [SyntheticSingleTag::~~SyntheticSingleTag](#) ( )

### 8.35.3 Member Function Documentation

8.35.3.1 `channel_t SyntheticSingleTag::getChannel ( ) const`

8.35.3.2 `bool SyntheticSingleTag::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time ) [override], [protected], [virtual]`

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

#### Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

#### Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

8.35.3.3 `void SyntheticSingleTag::trigger ( )`

Generate a timetag for each call of this method.

## 8.35.4 Friends And Related Function Documentation

8.35.4.1 `friend class SyntheticSingleTagImpl [friend]`

The documentation for this class was generated from the following file:

- [Iterators.h](#)

## 8.36 Tag Struct Reference

a single event on a channel

```
#include <TimeTagger.h>
```

### Public Types

- enum [Type](#) : unsigned char {  
[Type::TimeTag](#) = 0, [Type::Error](#) = 1, [Type::OverflowBegin](#) = 2, [Type::OverflowEnd](#) = 3,  
[Type::MissedEvents](#) = 4 }

## Public Attributes

- enum [Tag::Type](#) type
- char [reserved](#)
- unsigned short [missed\\_events](#)
- [channel\\_t](#) channel
- [timestamp\\_t](#) time

### 8.36.1 Detailed Description

a single event on a channel

Channel events are passed from the backend to registered iterators by the `IteratorBase::next()` callback function.

A [Tag](#) describes a single event on a channel.

### 8.36.2 Member Enumeration Documentation

#### 8.36.2.1 enum `Tag::Type` : unsigned char `[strong]`

This enum marks what kind of event this object represents: `TimeTag`: a normal event from any input channel  
`Error`: an error in the internal data processing, e.g. on plugging the external clock. This invalidates the global time  
`OverflowBegin`: this marks the begin of an interval with incomplete data because of too high data rates `Overflow`↔  
`End`: this marks the end of the interval. All events, which were lost in this interval, have been handled `Missed`↔  
`Events`: this virtual event signals the amount of lost events per channel within an overflow interval. Repeated usage for higher amounts of events

Enumerator

***TimeTag***  
***Error***  
***OverflowBegin***  
***OverflowEnd***  
***MissedEvents***

### 8.36.3 Member Data Documentation

#### 8.36.3.1 `channel_t` `Tag::channel`

#### 8.36.3.2 unsigned short `Tag::missed_events`

#### 8.36.3.3 char `Tag::reserved`

#### 8.36.3.4 `timestamp_t` `Tag::time`

#### 8.36.3.5 enum `Tag::Type` `Tag::type`

The documentation for this struct was generated from the following file:


- [TimeTagger.h](#)

## 8.37 TimeDifferences Class Reference

Accumulates the time differences between clicks on two channels in one or more histograms.

```
#include <Iterators.h>
```

Inheritance diagram for TimeDifferences:



### Public Member Functions

- [TimeDifferences](#) ([TimeTaggerBase](#) \*tagger, [channel\\_t](#) click\_channel, [channel\\_t](#) start\_channel=[CHANNEL\\_UNUSED](#), [channel\\_t](#) next\_channel=[CHANNEL\\_UNUSED](#), [channel\\_t](#) sync\_channel=[CHANNEL\\_UNUSED](#), [timestamp\\_t](#) binwidth=1000, [int32\\_t](#) n\_bins=1000, [int32\\_t](#) n\_histograms=1)  
*constructor of a [TimeDifferences](#) measurement*
- [~TimeDifferences](#) ()
- void [getData](#) (std::function< [int32\\_t](#) \*([size\\_t](#), [size\\_t](#))> array\_out)  
*returns a two-dimensional array of size 'n\_bins' by 'n\_histograms' containing the histograms*
- void [getIndex](#) (std::function< long long \*([size\\_t](#))> array\_out)  
*returns a vector of size 'n\_bins' containing the time bins in ps*
- void [setMaxCounts](#) ([uint64\\_t](#) max\_counts)  
*set the number of rollovers at which the measurement stops integrating*
- [uint64\\_t](#) [getCounts](#) ()  
*returns the number of rollovers (histogram index resets)*
- bool [ready](#) ()  
*returns 'true' when the required number of rollovers set by 'setMaxCounts' has been reached*

### Protected Member Functions

- bool [next\\_impl](#) (std::vector< [Tag](#) > &incoming\_tags, [timestamp\\_t](#) begin\_time, [timestamp\\_t](#) end\_time) override  
*update iterator state*
- void [clear\\_impl](#) () override  
*clear [Iterator](#) state.*
- void [on\\_start](#) () override  
*callback when the measurement class is started*

## Friends

- class [TimeDifferencesImpl](#) < [TimeDifferences](#) >

## Additional Inherited Members

### 8.37.1 Detailed Description

Accumulates the time differences between clicks on two channels in one or more histograms.

A multidimensional histogram measurement with the option up to include three additional channels that control how to step through the indices of the histogram array. This is a very powerful and generic measurement. You can use it to record cross-correlation, lifetime measurements, fluorescence lifetime imaging and many more measurements based on pulsed excitation. Specifically, the measurement waits for a tag on the 'start\_channel', then measures the time difference between the start tag and all subsequent tags on the 'click\_channel' and stores them in a histogram. If no 'start\_channel' is specified, the 'click\_channel' is used as 'start\_channel' corresponding to an auto-correlation measurement. The histogram has a number 'n\_bins' of bins of bin width 'binwidth'. Clicks that fall outside the histogram range are discarded. Data accumulation is performed independently for all start tags. This type of measurement is frequently referred to as 'multiple start, multiple stop' measurement and corresponds to a full auto- or cross-correlation measurement.

The data obtained from subsequent start tags can be accumulated into the same histogram (one- dimensional measurement) or into different histograms (two-dimensional measurement). In this way, you can perform more general two-dimensional time-difference measurements. The parameter 'n\_histograms' specifies the number of histograms. After each tag on the 'next\_channel', the histogram index is incremented by one and reset to zero after reaching the last valid index. The measurement starts with the first tag on the 'next\_channel'.

You can also provide a synchronization trigger that resets the histogram index by specifying a 'sync\_channel'. The measurement starts when a tag on the 'sync\_channel' arrives with a subsequent tag on 'next\_channel'. When a rollover occurs, the accumulation is stopped until the next sync and subsequent next signal. A sync signal before a rollover will stop the accumulation, reset the histogram index and a subsequent signal on the 'next\_channel' starts the accumulation again.

Typically, you will run the measurement indefinitely until stopped by the user. However, it is also possible to specify the maximum number of rollovers of the histogram index. In this case the measurement stops when the number of rollovers has reached the specified value. This means that for both a one-dimensional and for a two-dimensional measurement, it will measure until the measurement went through the specified number of rollovers / sync tags.

### 8.37.2 Constructor & Destructor Documentation

**8.37.2.1** `TimeDifferences::TimeDifferences ( TimeTaggerBase * tagger, channel_t click_channel, channel_t start_channel = CHANNEL_UNUSED, channel_t next_channel = CHANNEL_UNUSED, channel_t sync_channel = CHANNEL_UNUSED, timestamp_t binwidth = 1000, int32_t n_bins = 1000, int32_t n_histograms = 1 )`

constructor of a [TimeDifferences](#) measurement

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>click_channel</i>	channel that increments the count in a bin
<i>start_channel</i>	channel that sets start times relative to which clicks on the click channel are measured
<i>next_channel</i>	channel that increments the histogram index
<i>sync_channel</i>	channel that resets the histogram index to zero
<i>binwidth</i>	width of one histogram bin in ps
<i>n_bins</i>	number of bins in each histogram
<i>n_histograms</i>	number of histograms



## 8.37.2.2 TimeDifferences::~~TimeDifferences ( )

## 8.37.3 Member Function Documentation

## 8.37.3.1 void TimeDifferences::clear\_impl ( ) [override],[protected],[virtual]

clear [Iterator](#) state.

Each [Iterator](#) should implement the [clear\\_impl\(\)](#) method to reset its internal state. The [clear\\_impl\(\)](#) function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

## 8.37.3.2 uint64\_t TimeDifferences::getCounts ( )

returns the number of rollovers (histogram index resets)

## 8.37.3.3 void TimeDifferences::getData ( std::function&lt; int32\_t \*(size\_t, size\_t)&gt; array\_out )

returns a two-dimensional array of size 'n\_bins' by 'n\_histograms' containing the histograms

## 8.37.3.4 void TimeDifferences::getIndex ( std::function&lt; long long \*(size\_t)&gt; array\_out )

returns a vector of size 'n\_bins' containing the time bins in ps

## 8.37.3.5 bool TimeDifferences::next\_impl ( std::vector&lt; Tag &gt; &amp; incoming\_tags, timestamp\_t begin\_time, timestamp\_t end\_time ) [override],[protected],[virtual]

update iterator state

Each [Iterator](#) must implement the [next\\_impl\(\)](#) method. The [next\\_impl\(\)](#) function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

## Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

## Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

**8.37.3.6** `void TimeDifferences::on_start ( )` `[override]`, `[protected]`, `[virtual]`

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

**8.37.3.7** `bool TimeDifferences::ready ( )`

returns 'true' when the required number of rollovers set by 'setMaxCounts' has been reached

**8.37.3.8** `void TimeDifferences::setMaxCounts ( uint64_t max_counts )`

set the number of rollovers at which the measurement stops integrating

Parameters

<i>max_counts</i>	maximum number of sync/next clicks
-------------------	------------------------------------

## 8.37.4 Friends And Related Function Documentation

**8.37.4.1** `friend class TimeDifferencesImpl< TimeDifferences >` `[friend]`

The documentation for this class was generated from the following file:

- [Iterators.h](#)

## 8.38 TimeDifferencesImpl< T > Class Template Reference

```
#include <Iterators.h>
```

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.39 TimeDifferencesND Class Reference

Accumulates the time differences between clicks on two channels in a multi-dimensional histogram.

```
#include <Iterators.h>
```

Inheritance diagram for TimeDifferencesND:



### Public Member Functions

- `TimeDifferencesND` (`TimeTaggerBase` \*tagger, `channel_t` click\_channel, `channel_t` start\_channel, `std::vector`< `channel_t` > next\_channels, `std::vector`< `channel_t` > sync\_channels, `std::vector`< `int32_t` > n\_hists, `timestamp_t` binwidth, `int32_t` n\_bins)  
*constructor of a TimeDifferencesND measurement*
- `~TimeDifferencesND` ()
- void `getData` (`std::function`< `int32_t` \*(`size_t`, `size_t`)> array\_out)  
*returns a two-dimensional array of size n\_bins by all n\_hists containing the histograms*
- void `getIndex` (`std::function`< `long long` \*(`size_t`)> array\_out)  
*returns a vector of size n\_bins containing the time bins in ps*

### Protected Member Functions

- bool `next_impl` (`std::vector`< `Tag` > &incoming\_tags, `timestamp_t` begin\_time, `timestamp_t` end\_time) override  
*update iterator state*
- void `clear_impl` () override  
*clear Iterator state.*
- void `on_start` () override  
*callback when the measurement class is started*

### Friends

- class `TimeDifferencesNDImpl`

## Additional Inherited Members

### 8.39.1 Detailed Description

Accumulates the time differences between clicks on two channels in a multi-dimensional histogram.

This is a multidimensional implementation of the [TimeDifferences](#) measurement class. Please read their documentation first.

This measurement class extends the [TimeDifferences](#) interface for a multidimensional amount of histograms. It captures many multiple start - multiple stop histograms, but with many asynchronous next\_channel triggers. After each tag on each next\_channel, the histogram index of the associated dimension is incremented by one and reset to zero after reaching the last valid index. The elements of the parameter n\_histograms specifies the number of histograms per dimension. The accumulation starts when next\_channel has been triggered on all dimensions.

You should provide a synchronization trigger by specifying a sync\_channel per dimension. It will stop the accumulation when an associated histogram index rollover occurs. A sync event will also stop the accumulation, reset the histogram index of the associated dimension, and a subsequent event on the corresponding next\_channel starts the accumulation again. The synchronization is done asynchronous, so an event on the next\_channel increases the histogram index even if the accumulation is stopped. The accumulation starts when a tag on the sync\_channel arrives with a subsequent tag on next\_channel for all dimensions.

Please use setInputDelay to adjust the latency of all channels. In general, the order of the provided triggers including maximum jitter should be: old start trigger – all sync triggers – all next triggers – new start trigger

### 8.39.2 Constructor & Destructor Documentation

**8.39.2.1** `TimeDifferencesND::TimeDifferencesND ( TimeTaggerBase * tagger, channel_t click_channel, channel_t start_channel, std::vector< channel_t > next_channels, std::vector< channel_t > sync_channels, std::vector< int32_t > n_histograms, timestamp_t binwidth, int32_t n_bins )`

constructor of a [TimeDifferencesND](#) measurement

#### Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>click_channel</i>	channel that increments the count in a bin
<i>start_channel</i>	channel that sets start times relative to which clicks on the click channel are measured
<i>next_channels</i>	vector of channels that increments the histogram index
<i>sync_channels</i>	vector of channels that resets the histogram index to zero
<i>n_histograms</i>	vector of numbers of histograms per dimension.
<i>binwidth</i>	width of one histogram bin in ps
<i>n_bins</i>	number of bins in each histogram

**8.39.2.2** `TimeDifferencesND::~~TimeDifferencesND ( )`

### 8.39.3 Member Function Documentation

8.39.3.1 `void TimeDifferencesND::clear_impl ( ) [override],[protected],[virtual]`

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

8.39.3.2 `void TimeDifferencesND::getData ( std::function< int32_t *(size_t, size_t)> array_out )`

returns a two-dimensional array of size `n_bins` by all `n_histograms` containing the histograms

8.39.3.3 `void TimeDifferencesND::getIndex ( std::function< long long *(size_t)> array_out )`

returns a vector of size `n_bins` containing the time bins in ps

8.39.3.4 `bool TimeDifferencesND::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time ) [override],[protected],[virtual]`

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

#### Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

#### Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

8.39.3.5 `void TimeDifferencesND::on_start ( ) [override],[protected],[virtual]`

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

### 8.39.4 Friends And Related Function Documentation

#### 8.39.4.1 friend class TimeDifferencesNDImpl [friend]

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.40 TimeTagger Class Reference

backend for the [TimeTagger](#).

```
#include <TimeTagger.h>
```

Inheritance diagram for TimeTagger:



### Public Member Functions

- virtual void [reset](#) ()=0  
*reset the [TimeTagger](#) object to default settings and detach all iterators*
- virtual void [setTestSignalDivider](#) (int divider)=0  
*set the divider for the frequency of the test signal*
- virtual int [getTestSignalDivider](#) ()=0  
*get the divider for the frequency of the test signal*
- virtual void [setTriggerLevel](#) ([channel\\_t](#) channel, double voltage)=0  
*set the trigger voltage threshold of a channel*
- virtual double [getTriggerLevel](#) ([channel\\_t](#) channel)=0  
*get the trigger voltage threshold of a channel*
- virtual [timestamp\\_t](#) [getHardwareDelayCompensation](#) ([channel\\_t](#) channel)=0  
*get hardware delay compensation of a channel*
- virtual void [setInputMux](#) ([channel\\_t](#) channel, int mux\_mode)=0  
*configures the input multiplexer*
- virtual int [getInputMux](#) ([channel\\_t](#) channel)=0  
*fetches the configuration of the input multiplexer*

- virtual void [setConditionalFilter](#) (std::vector< [channel\\_t](#) > trigger, std::vector< [channel\\_t](#) > filtered, bool hardwareDelayCompensation=true)=0  
*configures the conditional filter*
- virtual void [clearConditionalFilter](#) ()=0  
*deactivates the conditional filter*
- virtual std::vector< [channel\\_t](#) > [getConditionalFilterTrigger](#) ()=0  
*fetches the configuration of the conditional filter*
- virtual std::vector< [channel\\_t](#) > [getConditionalFilterFiltered](#) ()=0  
*fetches the configuration of the conditional filter*
- virtual void [setNormalization](#) (std::vector< [channel\\_t](#) > channel, bool state)=0  
*enables or disables the normalization of the distribution.*
- virtual bool [getNormalization](#) ([channel\\_t](#) channel)=0  
*returns the the normalization of the distribution.*
- virtual void [setHardwareBufferSize](#) (int size)=0  
*sets the maximum USB buffer size*
- virtual int [getHardwareBufferSize](#) ()=0  
*queries the size of the USB queue*
- virtual void [setStreamBlockSize](#) (int max\_events, int max\_latency)=0  
*sets the maximum events and latency for the stream block size*
- virtual int [getStreamBlockSizeEvents](#) ()=0
- virtual int [getStreamBlockSizeLatency](#) ()=0
- virtual void [setEventDivider](#) ([channel\\_t](#) channel, unsigned int divider)=0  
*Divides the amount of transmitted edge per channel.*
- virtual unsigned int [getEventDivider](#) ([channel\\_t](#) channel)=0  
*Returns the factor of the dividing filter.*
- virtual void [autoCalibration](#) (std::function< double \*(size\_t)> array\_out)=0  
*runs a calibrations based on the on-chip uncorrelated signal generator.*
- virtual std::string [getSerial](#) ()=0  
*identifies the hardware by serial number*
- virtual std::string [getModel](#) ()=0  
*identifies the hardware by Time Tagger Model*
- virtual int [getChannelNumberScheme](#) ()=0  
*Fetch the configured numbering scheme for this [TimeTagger](#) object.*
- virtual std::vector< double > [getDACRange](#) ()=0  
*returns the minumum and the maximum voltage of the DACs as a trigger reference*
- virtual void [getDistributionCount](#) (std::function< uint64\_t \*(size\_t, size\_t)> array\_out)=0  
*get internal calibration data*
- virtual void [getDistributionPSEcs](#) (std::function< long long \*(size\_t, size\_t)> array\_out)=0  
*get internal calibration data This method is not supported any more on the Time Tagger Ultra series*
- virtual std::vector< [channel\\_t](#) > [getChannelList](#) ([ChannelEdge](#) type=[ChannelEdge::All](#))=0  
*fetch a vector of all physical input channel ids*
- virtual [timestamp\\_t](#) [getPsPerClock](#) ()=0  
*fetch the duration of each clock cycle in picoseconds*
- virtual std::string [getPcbVersion](#) ()=0  
*Return the hardware version of the PCB board. Version 0 is everything before mid 2018 and with the channel configuration ZERO. version >= 1 is channel configuration ONE.*
- virtual std::string [getFirmwareVersion](#) ()=0  
*Return an unique identifier for the applied firmware.*
- virtual std::string [getSensorData](#) ()=0  
*Show the status of the sensor data from the FPGA and peripherals on the console.*
- virtual void [setLED](#) (uint32\_t bitmask)=0

Enforce a state to the LEDs 0: led\_status[R] 16: led\_status[R] - mux 1: led\_status[G] 17: led\_status[G] - mux 2: led\_status[B] 18: led\_status[B] - mux 3: led\_power[R] 19: led\_power[R] - mux 4: led\_power[G] 20: led\_power[G] - mux 5: led\_power[B] 21: led\_power[B] - mux 6: led\_clock[R] 22: led\_clock[R] - mux 7: led\_clock[G] 23: led\_clock[G] - mux 8: led\_clock[B] 24: led\_clock[B] - mux.

- virtual std::string [getLicenseInfo](#) ()=0
- virtual uint32\_t [factoryAccess](#) (uint32\_t pw, uint32\_t addr, uint32\_t data, uint32\_t mask)=0  
Direct read/write access to WireIn/WireOuts in FPGA (mask==0 for readonly)
- virtual void [setSoundFrequency](#) (uint32\_t freq\_hz)=0  
Set the Time Taggers internal buzzer to a frequency in Hz (freq\_hz==0 to disable)

## Additional Inherited Members

### 8.40.1 Detailed Description

backend for the [TimeTagger](#).

The [TimeTagger](#) class connects to the hardware, and handles the communication over the usb. There may be only one instance of the backend per physical device.

### 8.40.2 Member Function Documentation

**8.40.2.1** virtual void [TimeTagger::autoCalibration](#) ( std::function< double \*(size\_t)> *array\_out* ) [pure virtual]

runs a calibrations based on the on-chip uncorrelated signal generator.

**8.40.2.2** virtual void [TimeTagger::clearConditionalFilter](#) ( ) [pure virtual]

deactivates the conditional filter

equivilent to setConditionalFilter({},{})

**8.40.2.3** virtual uint32\_t [TimeTagger::factoryAccess](#) ( uint32\_t *pw*, uint32\_t *addr*, uint32\_t *data*, uint32\_t *mask* ) [pure virtual]

Direct read/write access to WireIn/WireOuts in FPGA (mask==0 for readonly)

DO NOT USE. Only for internal debug purposes.

**8.40.2.4** virtual std::vector<channel\_t> [TimeTagger::getChannelList](#) ( [ChannelEdge](#) *type* = [ChannelEdge::All](#) ) [pure virtual]

fetch a vector of all physical input channel ids

The function returns the channel of all rising and falling edges. For example for the Time Tagger 20 (8 input channels) TT\_CHANNEL\_NUMBER\_SCHEME\_ZERO: {0,1,2,3,4,5,6,7,8,9,10,11,12,13,14,15} and for TT\_CHANNEL\_NUMBER\_SCHEME\_ONE: {-8,-7,-6,-5,-4,-3,-2,-1,1,2,3,4,5,6,7,8}

TT\_CHANNEL\_RISING\_EDGES returns only the rising edges SCHEME\_ONE: {1,2,3,4,5,6,7,8} and TT\_CHANNEL\_FALLING\_EDGES return only the falling edges SCHEME\_ONE: {-1,-2,-3,-4,-5,-6,-7,-8} which are the invertedChannels of the rising edges.



8.40.2.5 `virtual int TimeTagger::getChannelNumberScheme ( ) [pure virtual]`

Fetch the configured numbering scheme for this [TimeTagger](#) object.

Please see [setTimeTaggerChannelNumberScheme\(\)](#) for details.

8.40.2.6 `virtual std::vector<channel_t> TimeTagger::getConditionalFilterFiltered ( ) [pure virtual]`

fetches the configuration of the conditional filter

see [setConditionalFilter](#)

8.40.2.7 `virtual std::vector<channel_t> TimeTagger::getConditionalFilterTrigger ( ) [pure virtual]`

fetches the configuration of the conditional filter

see [setConditionalFilter](#)

8.40.2.8 `virtual std::vector<double> TimeTagger::getDACRange ( ) [pure virtual]`

returns the minumum and the maximum voltage of the DACs as a trigger reference

8.40.2.9 `virtual void TimeTagger::getDistributionCount ( std::function< uint64_t *(size_t, size_t)> array_out ) [pure virtual]`

get internal calibration data

8.40.2.10 `virtual void TimeTagger::getDistributionPSecs ( std::function< long long *(size_t, size_t)> array_out ) [pure virtual]`

get internal calibration data This method is not supported any more on the Time Tagger Ultra series

### Deprecated

8.40.2.11 `virtual unsigned int TimeTagger::getEventDivider ( channel_t channel ) [pure virtual]`

Returns the factor of the dividing filter.

See [setEventDivider](#) for further details.

#### Parameters

<i>channel</i>	channel to be queried
----------------	-----------------------

**Returns**

the configured divider

**8.40.2.12** `virtual std::string TimeTagger::getFirmwareVersion ( ) [pure virtual]`

Return an unique identifier for the applied firmware.

This function returns a comma separated list of the firmware version with

- the device identifier: TT-20 or TT-Ultra
- the firmware identifier: FW 3
- optional the timestamp of the assembling of the firmware
- the firmware identifier of the USB chip: OK 1.30 eg "TT-Ultra, FW 3, TS 2018-11-13 22:57:32, OK 1.30"

**8.40.2.13** `virtual int TimeTagger::getHardwareBufferSize ( ) [pure virtual]`

queries the size of the USB queue

See setHardwareBufferSize for more information.

**Returns**

the actual size of the USB queue in events

**8.40.2.14** `virtual timestamp_t TimeTagger::getHardwareDelayCompensation ( channel_t channel ) [pure virtual]`

get hardware delay compensation of a channel

The physical input delays are calibrated and compensated. However this compensation is implemented after the conditional filter and so affects its result. This function queries the effective input delay, which compensates the hardware delay.

**Parameters**

<i>channel</i>	the channel
----------------	-------------

**Returns**

the hardware delay compensation in picoseconds

**8.40.2.15** `virtual int TimeTagger::getInputMux ( channel_t channel ) [pure virtual]`

fetches the configuration of the input multiplexer

## Parameters

<i>channel</i>	the physcal channel of the input multiplexer
----------------	--

## Returns

the configuration mode of the input multiplexer

**8.40.2.16** `virtual std::string TimeTagger::getLicenseInfo ( ) [pure virtual]`

Fetches and parses the current installed license on this device

## Returns

a human readable string containing all information about the license on this device

**8.40.2.17** `virtual std::string TimeTagger::getModel ( ) [pure virtual]`

identifies the hardware by Time Tagger Model

**8.40.2.18** `virtual bool TimeTagger::getNormalization ( channel_t channel ) [pure virtual]`

returns the the normalization of the distribution.

Refer the Manual for a description of this function.

## Parameters

<i>channel</i>	the channel to query
----------------	----------------------

## Returns

if the normalization is enabled

**8.40.2.19** `virtual std::string TimeTagger::getPcbVersion ( ) [pure virtual]`

Return the hardware version of the PCB board. Version 0 is everything before mid 2018 and with the channel configuration ZERO. version  $\geq 1$  is channel configuration ONE.

**8.40.2.20** `virtual timestamp_t TimeTagger::getPsPerClock ( ) [pure virtual]`

fetch the duration of each clock cycle in picoseconds

8.40.2.21 `virtual std::string TimeTagger::getSensorData ( ) [pure virtual]`

Show the status of the sensor data from the FPGA and peripherals on the console.

8.40.2.22 `virtual std::string TimeTagger::getSerial ( ) [pure virtual]`

identifies the hardware by serial number

8.40.2.23 `virtual int TimeTagger::getStreamBlockSizeEvents ( ) [pure virtual]`

8.40.2.24 `virtual int TimeTagger::getStreamBlockSizeLatency ( ) [pure virtual]`

8.40.2.25 `virtual int TimeTagger::getTestSignalDivider ( ) [pure virtual]`

get the divider for the frequency of the test signal

8.40.2.26 `virtual double TimeTagger::getTriggerLevel ( channel_t channel ) [pure virtual]`

get the trigger voltage threshold of a channel

#### Parameters

<i>channel</i>	the channel
----------------	-------------

8.40.2.27 `virtual void TimeTagger::reset ( ) [pure virtual]`

reset the [TimeTagger](#) object to default settings and detach all iterators

8.40.2.28 `virtual void TimeTagger::setConditionalFilter ( std::vector< channel_t > trigger, std::vector< channel_t > filtered, bool hardwareDelayCompensation = true ) [pure virtual]`

configures the conditional filter

After each event on the trigger channels, one event per filtered channel will pass afterwards. This is handled in a very early stage in the pipeline, so all event limitations but the deadtime are suppressed. But the accuracy of the order of those events is low.

Refer the Manual for a description of this function.

#### Parameters

<i>trigger</i>	the channels that sets the condition
<i>filtered</i>	the channels that are filtered by the condition
<i>hardwareDelayCompensation</i>	if false, the physical hardware delay will not be compensated

**8.40.2.29** `virtual void TimeTagger::setEventDivider ( channel_t channel, unsigned int divider )` `[pure virtual]`

Divides the amount of transmitted edge per channel.

This filter decimates the events on a given channel by a specified factor. So for a divider  $n$ , every  $n$ th event is transmitted through the filter and  $n-1$  events are skipped between consecutive transmitted events. If a conditional filter is also active, the event divider is applied after the conditional filter, so the conditional is applied to the complete event stream and only events which pass the conditional filter are forwarded to the divider.

As it is a hardware filter, it reduces the required USB bandwidth and CPU processing power, but it cannot be configured for virtual channels.

#### Parameters

<i>channel</i>	channel to be configured
<i>divider</i>	new divider, must be smaller than 65536

**8.40.2.30** `virtual void TimeTagger::setHardwareBufferSize ( int size )` `[pure virtual]`

sets the maximum USB buffer size

This option controls the maximum buffer size of the USB connection. This can be used to balance low input latency vs high (peak) throughput.

#### Parameters

<i>size</i>	the maximum buffer size in events
-------------	-----------------------------------

**8.40.2.31** `virtual void TimeTagger::setInputMux ( channel_t channel, int mux_mode )` `[pure virtual]`

configures the input multiplexer

Every physical input channel has an input multiplexer with 4 modes: 0: normal input mode 1: use the input from channel -1 (left) 2: use the input from channel +1 (right) 3: use the reference oscillator

Mode 1 and 2 cascades, so many inputs can be configured to get the same input events.

#### Parameters

<i>channel</i>	the physical channel of the input multiplexer
<i>mux_mode</i>	the configuration mode of the input multiplexer

**8.40.2.32** `virtual void TimeTagger::setLED ( uint32_t bitmask )` `[pure virtual]`

Enforce a state to the LEDs 0: led\_status[R] 16: led\_status[R] - mux 1: led\_status[G] 17: led\_status[G] - mux 2: led\_status[B] 18: led\_status[B] - mux 3: led\_power[R] 19: led\_power[R] - mux 4: led\_power[G] 20: led\_power[G] - mux 5: led\_power[B] 21: led\_power[B] - mux 6: led\_clock[R] 22: led\_clock[R] - mux 7: led\_clock[G] 23: led\_clock[G] - mux 8: led\_clock[B] 24: led\_clock[B] - mux.

**8.40.2.33** `virtual void TimeTagger::setNormalization ( std::vector< channel_t > channel, bool state )` [pure virtual]

enables or disables the normalization of the distribution.

Refer the Manual for a description of this function.

#### Parameters

<i>channel</i>	list of channels to modify
<i>state</i>	the new state

**8.40.2.34** `virtual void TimeTagger::setSoundFrequency ( uint32_t freq_hz )` [pure virtual]

Set the Time Taggers internal buzzer to a frequency in Hz (freq\_hz==0 to disable)

#### Parameters

<i>freq_hz</i>	the generated audio frequency
----------------	-------------------------------

**8.40.2.35** `virtual void TimeTagger::setStreamBlockSize ( int max_events, int max_latency )` [pure virtual]

sets the maximum events and latency for the stream block size

This option controls the latency and the block size of the data stream. The default values are max\_events = 131072 events and max\_latency = 20 ms. Depending on which of the two parameters is exceeded first, the block stream size is adjusted accordingly. The block size will be reduced automatically for blocks when no signal is arriving for 512 ns on the Time Tagger Ultra and 1536 ns for the Time Tagger 20. \*

#### Parameters

<i>max_events</i>	maximum number of events
<i>max_latency</i>	maximum latency in ms

**8.40.2.36** `virtual void TimeTagger::setTestSignalDivider ( int divider )` [pure virtual]

set the divider for the frequency of the test signal

The base clock of the test signal oscillator for the Time Tagger Ultra is running at 100.8 MHz sampled down by an factor of 2 to have a similar base clock as the Time Tagger 20 (~50 MHz). The default divider is 63 -> ~800 kEvents/s

#### Parameters

<i>divider</i>	frequency divisor of the oscillator
----------------	-------------------------------------

8.40.2.37 `virtual void TimeTagger::setTriggerLevel ( channel_t channel, double voltage ) [pure virtual]`

set the trigger voltage threshold of a channel

#### Parameters

<i>channel</i>	the channel to set
<i>voltage</i>	voltage level.. [0..1]


The documentation for this class was generated from the following file:

- [TimeTagger.h](#)

## 8.41 TimeTaggerBase Class Reference

```
#include <TimeTagger.h>
```

Inheritance diagram for TimeTaggerBase:



### Public Types

- `typedef std::function< void(IteratorBase *)> IteratorCallback`
- `typedef std::map< IteratorBase *, IteratorCallback > IteratorCallbackMap`

### Public Member Functions

- `virtual unsigned int getFence (bool alloc_fence=true)=0`
- `virtual bool waitForFence (unsigned int fence, int64_t timeout=-1)=0`
- `virtual bool sync (int64_t timeout=-1)=0`
- `virtual channel_t getInvertedChannel (channel_t channel)=0`  
*get the falling channel id for a raising channel and vice versa*
- `virtual bool isUnusedChannel (channel_t channel)=0`  
*compares the provided channel with CHANNEL\_UNUSED*
- `virtual void runSynchronized (const IteratorCallbackMap &callbacks, bool block=true)=0`

*Run synchronized callbacks for a list of iterators.*

- virtual std::string `getConfiguration` ()=0
- virtual void `setInputDelay` (channel\_t channel, timestamp\_t delay)=0  
*set time delay on a channel*
- virtual void `setDelayHardware` (channel\_t channel, timestamp\_t delay)=0  
*set time delay on a channel*
- virtual void `setDelaySoftware` (channel\_t channel, timestamp\_t delay)=0  
*set time delay on a channel*
- virtual timestamp\_t `getInputDelay` (channel\_t channel)=0  
*get time delay of a channel*
- virtual timestamp\_t `getDelaySoftware` (channel\_t channel)=0  
*get time delay of a channel*
- virtual timestamp\_t `getDelayHardware` (channel\_t channel)=0  
*get time delay of a channel*
- virtual timestamp\_t `setDeadtime` (channel\_t channel, timestamp\_t deadtime)=0  
*set the deadtime between two edges on the same channel.*
- virtual timestamp\_t `getDeadtime` (channel\_t channel)=0  
*get the deadtime between two edges on the same channel.*
- virtual void `setTestSignal` (channel\_t channel, bool enabled)=0  
*enable the calibration on a channel.*
- virtual void `setTestSignal` (std::vector< channel\_t > channel, bool enabled)=0
- virtual bool `getTestSignal` (channel\_t channel)=0  
*fetch the status of the test signal generator*
- virtual long long `getOverflows` ()=0  
*get overflow count*
- virtual void `clearOverflows` ()=0  
*clear overflow counter*
- virtual long long `getOverflowsAndClear` ()=0  
*get and clear overflow counter*

## Protected Member Functions

- `TimeTaggerBase` ()  
*abstract interface class*
- virtual `~TimeTaggerBase` ()
- `TimeTaggerBase` (const `TimeTaggerBase` &)=delete
- `TimeTaggerBase` & `operator=` (const `TimeTaggerBase` &)=delete
- virtual std::shared\_ptr< IteratorBaseListNode > `addIterator` (IteratorBase \*it)=0
- virtual void `freeIterator` (IteratorBase \*it)=0
- virtual channel\_t `getNewVirtualChannel` ()=0
- virtual void `freeVirtualChannel` (channel\_t channel)=0
- virtual void `registerChannel` (channel\_t channel)=0  
*register a FPGA channel.*
- virtual void `unregisterChannel` (channel\_t channel)=0  
*release a previously registered channel.*
- virtual void `addChild` (`TimeTaggerBase` \*child)=0
- virtual void `removeChild` (`TimeTaggerBase` \*child)=0
- virtual void `release` ()=0



## Friends

- class [IteratorBase](#)
- class [TimeTaggerProxy](#)
- class [TimeTaggerRunner](#)

## 8.41.1 Member Typedef Documentation

8.41.1.1 `typedef std::function<void(IteratorBase *)> TimeTaggerBase::IteratorCallback`

8.41.1.2 `typedef std::map<IteratorBase *, IteratorCallback> TimeTaggerBase::IteratorCallbackMap`

## 8.41.2 Constructor & Destructor Documentation

8.41.2.1 `TimeTaggerBase::TimeTaggerBase ( ) [inline], [protected]`

abstract interface class

8.41.2.2 `virtual TimeTaggerBase::~~TimeTaggerBase ( ) [inline], [protected], [virtual]`

destructor

8.41.2.3 `TimeTaggerBase::TimeTaggerBase ( const TimeTaggerBase & ) [protected], [delete]`

## 8.41.3 Member Function Documentation

8.41.3.1 `virtual void TimeTaggerBase::addChild ( TimeTaggerBase * child ) [protected], [pure virtual]`

8.41.3.2 `virtual std::shared_ptr<IteratorBaseListNode> TimeTaggerBase::addIterator ( IteratorBase * it ) [protected], [pure virtual]`

8.41.3.3 `virtual void TimeTaggerBase::clearOverflows ( ) [pure virtual]`

clear overflow counter

Sets the overflow counter to zero

8.41.3.4 `virtual void TimeTaggerBase::freeIterator ( IteratorBase * it ) [protected], [pure virtual]`

8.41.3.5 `virtual void TimeTaggerBase::freeVirtualChannel ( channel_t channel ) [protected], [pure virtual]`

8.41.3.6 `virtual std::string TimeTaggerBase::getConfiguration ( ) [pure virtual]`

Fetches the overall configuration status of the Time Tagger object.

Returns

a JSON serialized string with all configuration and status flags.

8.41.3.7 `virtual timestamp_t TimeTaggerBase::getDeadtime ( channel_t channel ) [pure virtual]`

get the deadtime between two edges on the same channel.

This function gets the user configureable deadtime.

## Parameters

<i>channel</i>	channel to be queried
----------------	-----------------------

## Returns

the real configured deadtime

**8.41.3.8** `virtual timestamp_t TimeTaggerBase::getDelayHardware ( channel_t channel )` [pure virtual]

get time delay of a channel

see setDelayHardware

## Parameters

<i>channel</i>	the channel
----------------	-------------

**8.41.3.9** `virtual timestamp_t TimeTaggerBase::getDelaySoftware ( channel_t channel )` [pure virtual]

get time delay of a channel

see setDelaySoftware

## Parameters

<i>channel</i>	the channel
----------------	-------------

**8.41.3.10** `virtual unsigned int TimeTaggerBase::getFence ( bool alloc_fence = true )` [pure virtual]

Generate a new fence object, which validates the current configuration and the current time. This fence is uploaded to the earliest pipeline stage of the Time Tagger. Waiting on this fence ensures that all hardware settings such as trigger levels, channel registrations, etc., have propagated to the FPGA and are physically active. Synchronizes the Time Tagger internal memory, so that all tags arriving after the waitFence call were actually produced after the getFence call. The waitFence function waits until all tags, which are present at the time of the function call within the internal memory of the Time Tagger, are processed. This call might block to limit the amount of active fences.

## Parameters

<i>alloc_fence</i>	if false, a reference to the most recently created fence will be returned instead
--------------------	---

## Returns

the allocated fence

8.41.3.11 `virtual timestamp_t TimeTaggerBase::getInputDelay ( channel_t channel ) [pure virtual]`

get time delay of a channel

see setInputDelay

Parameters

<i>channel</i>	the channel
----------------	-------------

8.41.3.12 `virtual channel_t TimeTaggerBase::getInvertedChannel ( channel_t channel ) [pure virtual]`

get the falling channel id for a raising channel and vice versa

8.41.3.13 `virtual channel_t TimeTaggerBase::getNewVirtualChannel ( ) [protected],[pure virtual]`

8.41.3.14 `virtual long long TimeTaggerBase::getOverflows ( ) [pure virtual]`

get overflow count

Get the number of communication overflows occurred

8.41.3.15 `virtual long long TimeTaggerBase::getOverflowsAndClear ( ) [pure virtual]`

get and clear overflow counter

Get the number of communication overflows occurred and sets them to zero

8.41.3.16 `virtual bool TimeTaggerBase::getTestSignal ( channel_t channel ) [pure virtual]`

fetch the status of the test signal generator

Parameters

<i>channel</i>	the channel
----------------	-------------

8.41.3.17 `virtual bool TimeTaggerBase::isUnusedChannel ( channel_t channel ) [pure virtual]`

compares the provided channel with CHANNEL\_UNUSED

But also keeps care about the channel number scheme and selects either CHANNEL\_UNUSED or CHANNEL\_UNUSED\_OLD

**8.41.3.18** `TimeTaggerBase& TimeTaggerBase::operator= ( const TimeTaggerBase & )` [protected], [delete]

**8.41.3.19** `virtual void TimeTaggerBase::registerChannel ( channel_t channel )` [protected], [pure virtual]

register a FPGA channel.

Only events on previously registered channels will be transfered over the communication channel.

Parameters

<i>channel</i>	the channel
----------------	-------------

**8.41.3.20** `virtual void TimeTaggerBase::release ( )` [protected], [pure virtual]

**8.41.3.21** `virtual void TimeTaggerBase::removeChild ( TimeTaggerBase * child )` [protected], [pure virtual]

**8.41.3.22** `virtual void TimeTaggerBase::runSynchronized ( const IteratorCallbackMap & callbacks, bool block=true )` [pure virtual]

Run synchronized callbacks for a list of iterators.

This method has a list of callbacks for a list of iterators. Those callbacks are called for a synchronized data set, but in parallel. They are called from an internal worker thread. As the data set is synchronized, this creates a bottleneck for one worker thread, so only fast and non-blocking callbacks are allowed.

Parameters

<i>callbacks</i>	Map of callbacks per iterator
<i>block</i>	Shall this method block until all callbacks are finished

**8.41.3.23** `virtual timestamp_t TimeTaggerBase::setDeadtime ( channel_t channel, timestamp_t deadtime )` [pure virtual]

set the deadtime between two edges on the same channel.

This function sets the user configureable deadtime. The requested time will be rounded to the nearest multiple of the clock time. The deadtime will also be clamped to device specific limitations.

As the actual deadtime will be altered, the real value will be returned.

Parameters

<i>channel</i>	channel to be configured
<i>deadtime</i>	new deadtime

**Returns**

the real configured deadtime

**8.41.3.24** `virtual void TimeTaggerBase::setDelayHardware ( channel_t channel, timestamp_t delay )` [pure virtual]

set time delay on a channel

When set, every event on this physical input channel is delayed by the given delay in picoseconds. This delay is implemented on the hardware before any filter with no performance overhead. The maximum delay on the Time Tagger Ultra series is 2 us. This affects both the rising and the falling event at the same time.

**Parameters**

<i>channel</i>	the channel to set
<i>delay</i>	the delay in picoseconds

**8.41.3.25** `virtual void TimeTaggerBase::setDelaySoftware ( channel_t channel, timestamp_t delay )` [pure virtual]

set time delay on a channel

When set, every event on this channel is delayed by the given delay in picoseconds. This happens on the computer and so after the on-device filters. Please use setDelayHardware instead for better performance. This affects either the the rising or the falling event only.

This method has the best performance with "small delays". The delay is considered "small" when less than 100 events arrive within the time of the largest delay set. For example, if the total event-rate over all channels used is 10 Mevent/s, the signal can be delayed efficiently up to 10 microseconds. For large delays, please use [DelayedChannel](#) instead.

**Parameters**

<i>channel</i>	the channel to set
<i>delay</i>	the delay in picoseconds

**8.41.3.26** `virtual void TimeTaggerBase::setInputDelay ( channel_t channel, timestamp_t delay )` [pure virtual]

set time delay on a channel

When set, every event on this channel is delayed by the given delay in picoseconds.

This method has the best performance with "small delays". The delay is considered "small" when less than 100 events arrive within the time of the largest delay set. For example, if the total event-rate over all channels used is 10 Mevent/s, the signal can be delayed efficiently up to 10 microseconds. For large delays, please use [DelayedChannel](#) instead.

## Parameters

<i>channel</i>	the channel to set
<i>delay</i>	the delay in picoseconds

8.41.3.27 `virtual void TimeTaggerBase::setTestSignal ( channel_t channel, bool enabled )` [pure virtual]

enable the calibration on a channel.

This will connect or disconnect the channel with the on-chip uncorrelated signal generator.

## Parameters

<i>channel</i>	the channel
<i>enabled</i>	enabled / disabled flag

8.41.3.28 `virtual void TimeTaggerBase::setTestSignal ( std::vector< channel_t > channel, bool enabled )` [pure virtual]

8.41.3.29 `virtual bool TimeTaggerBase::sync ( int64_t timeout = -1 )` [pure virtual]

Sync the timetagger pipeline, so that all started iterators and their enabled channels are ready This is a shortcut for calling `getFence` and `waitForFence` at once. See `getFence` for more details.

## Parameters

<i>timeout</i>	timeout in milliseconds. Negative means no timeout, zero returns immediately.
----------------	---

## Returns

true on success, false on timeout

8.41.3.30 `virtual void TimeTaggerBase::unregisterChannel ( channel_t channel )` [protected], [pure virtual]

release a previously registered channel.

## Parameters

<i>channel</i>	the channel
----------------	-------------

8.41.3.31 `virtual bool TimeTaggerBase::waitForFence ( unsigned int fence, int64_t timeout = -1 )` [pure virtual]

Wait for a fence in the data stream. See `getFence` for more details.

## Parameters

<i>fence</i>	fence object, which shall be waited on
<i>timeout</i>	timeout in milliseconds. Negative means no timeout, zero returns immediately.

## Returns

true if the fence has passed, false on timeout

### 8.41.4 Friends And Related Function Documentation

8.41.4.1 friend class `IteratorBase` `[friend]`

8.41.4.2 friend class `TimeTaggerProxy` `[friend]`

8.41.4.3 friend class `TimeTaggerRunner` `[friend]`

The documentation for this class was generated from the following file:


- [TimeTagger.h](#)

## 8.42 TimeTaggerVirtual Class Reference

virtual [TimeTagger](#) based on dump files

```
#include <TimeTagger.h>
```

Inheritance diagram for TimeTaggerVirtual:



## Public Member Functions

- virtual uint64\_t **replay** (const std::string &file, timestamp\_t begin=0, timestamp\_t duration=-1, bool queue=true)=0  
*replay a given dump file on the disc*
- virtual void **stop** ()=0  
*stops the current and all queued files.*
- virtual void **reset** ()=0  
*stops the all queued files and resets the *TimeTaggerVirtual* to its default settings*
- virtual bool **waitForCompletion** (uint64\_t ID=0, int64\_t timeout=-1)=0  
*block the current thread until the replay finish*
- virtual void **setReplaySpeed** (double speed)=0  
*configures the speed factor for the virtual tagger.*
- virtual double **getReplaySpeed** ()=0  
*fetches the speed factor*
- virtual void **setConditionalFilter** (std::vector< channel\_t > trigger, std::vector< channel\_t > filtered)=0  
*configures the conditional filter*
- virtual void **clearConditionalFilter** ()=0  
*deactivates the conditional filter*
- virtual std::vector< channel\_t > **getConditionalFilterTrigger** ()=0  
*fetches the configuration of the conditional filter*
- virtual std::vector< channel\_t > **getConditionalFilterFiltered** ()=0  
*fetches the configuration of the conditional filter*

## Additional Inherited Members

### 8.42.1 Detailed Description

virtual [TimeTagger](#) based on dump files

The [TimeTaggerVirtual](#) class represents a virtual Time Tagger. But instead of connecting to Swabians hardware, it replays all tags from a recorded file.

### 8.42.2 Member Function Documentation

#### 8.42.2.1 virtual void TimeTaggerVirtual::clearConditionalFilter ( ) [pure virtual]

deactivates the conditional filter

equivilent to setConditionalFilter({},{})

#### 8.42.2.2 virtual std::vector<channel\_t> TimeTaggerVirtual::getConditionalFilterFiltered ( ) [pure virtual]

fetches the configuration of the conditional filter

see setConditionalFilter



**8.42.2.3** `virtual std::vector<channel_t> TimeTaggerVirtual::getConditionalFilterTrigger ( ) [pure virtual]`

fetches the configuration of the conditional filter

see `setConditionalFilter`

**8.42.2.4** `virtual double TimeTaggerVirtual::getReplaySpeed ( ) [pure virtual]`

fetches the speed factor

Please see `setReplaySpeed` for more details.

#### Returns

the speed factor

**8.42.2.5** `virtual uint64_t TimeTaggerVirtual::replay ( const std::string & file, timestamp_t begin = 0, timestamp_t duration = -1, bool queue = true ) [pure virtual]`

replay a given dump file on the disc

This method adds the file to the replay queue. If the flag 'queue' is false, the current queue will be flushed and this file will be replayed immediatelly.

#### Parameters

<i>file</i>	the file to be replayed
<i>begin</i>	amount of ps to skip at the begin of the file. A negativ time will generate a pause in the replay
<i>duration</i>	time period in ps of the file. -1 replays till the last tag
<i>queue</i>	flag if this file shall be queued

#### Returns

ID of the queued file

**8.42.2.6** `virtual void TimeTaggerVirtual::reset ( ) [pure virtual]`

stops the all queued files and resets the [TimeTaggerVirtual](#) to its default settings

This method stops the current file, clears the replay queue and resets the [TimeTaggerVirtual](#) to its default settings.

**8.42.2.7** `virtual void TimeTaggerVirtual::setConditionalFilter ( std::vector< channel_t > trigger, std::vector< channel_t > filtered ) [pure virtual]`

configures the conditional filter

After each event on the trigger channels, one event per filtered channel will pass afterwards. This is handled in a very early stage in the pipeline, so all event limitations but the deadtime are suppressed. But the accuracy of the order of those events is low.

Refer the Manual for a description of this function.

## Parameters

<i>trigger</i>	the channels that sets the condition
<i>filtered</i>	the channels that are filtered by the condition

8.42.2.8 `virtual void TimeTaggerVirtual::setReplaySpeed ( double speed ) [pure virtual]`

configures the speed factor for the virtual tagger.

This method configures the speed factor of this virtual Time Tagger. A value of 1.0 will replay in real time. All values < 0.0 will replay the data as fast as possible, but stops at the end of all data. This is the default value.

## Parameters

<i>speed</i>	ratio of the replay speed and the real time
--------------	---

8.42.2.9 `virtual void TimeTaggerVirtual::stop ( ) [pure virtual]`

stops the current and all queued files.

This method stops the current file and clears the replay queue.

8.42.2.10 `virtual bool TimeTaggerVirtual::waitForCompletion ( uint64_t ID = 0, int64_t timeout = -1 ) [pure virtual]`

block the current thread until the replay finish

This method blocks the current execution and waits till the given file has finished its replay. If no ID is provided, it waits until all queued files are replayed.

This function does not block on a zero timeout. Negative timeouts are interpreted as infinite timeouts.

## Parameters

<i>ID</i>	selects which file to wait for
<i>timeout</i>	timeout in milliseconds

## Returns

true if the file is complete, false on timeout

The documentation for this class was generated from the following file:


- [TimeTagger.h](#)

## 8.43 TimeTagStream Class Reference

access the time tag stream

```
#include <Iterators.h>
```

Inheritance diagram for TimeTagStream:



### Public Member Functions

- `TimeTagStream (TimeTaggerBase *tagger, uint64_t n_max_events, std::vector< channel_t > channels=std::vector< channel_t >())`  
*constructor of a `TimeTagStream` thread*
- `~TimeTagStream ()`  
*tbd*
- `uint64_t getCounts ()`  
*get incoming time tags*
- `TimeTagStreamBuffer getData ()`  
*fetches all stored tags and clears the internal state*

### Protected Member Functions

- `bool next_impl (std::vector< Tag > &incoming_tags, timestamp_t begin_time, timestamp_t end_time)` override  
*update iterator state*
- `void clear_impl ()` override  
*clear `Iterator` state.*

### Friends

- class `TimeTagStreamImpl`

### Additional Inherited Members

#### 8.43.1 Detailed Description

access the time tag stream

### 8.43.2 Constructor & Destructor Documentation

8.43.2.1 `TimeTagStream::TimeTagStream ( TimeTaggerBase * tagger, uint64_t n_max_events, std::vector< channel_t > channels = std::vector< channel_t > () )`

constructor of a [TimeTagStream](#) thread

Gives access to the time tag stream

Parameters

<i>tagger</i>	reference to a <a href="#">TimeTagger</a>
<i>n_max_events</i>	maximum number of tags stored
<i>channels</i>	channels which are dumped to the file (when empty or not passed all active channels are dumped)

8.43.2.2 `TimeTagStream::~~TimeTagStream ( )`

tbd

### 8.43.3 Member Function Documentation

8.43.3.1 `void TimeTagStream::clear_impl ( )` `[override]`, `[protected]`, `[virtual]`

clear [Iterator](#) state.

Each [Iterator](#) should implement the `clear_impl()` method to reset its internal state. The `clear_impl()` function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

8.43.3.2 `uint64_t TimeTagStream::getCounts ( )`

get incoming time tags

All incoming time tags are stored in a buffer (max size: `max_tags`). The buffer is cleared after retrieving the data with `getData()` return the number of stored tags

8.43.3.3 `TimeTagStreamBuffer TimeTagStream::getData ( )`

fetches all stored tags and clears the internal state

8.43.3.4 `bool TimeTagStream::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time )` `[override]`, `[protected]`, `[virtual]`

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

## Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

## Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

## 8.43.4 Friends And Related Function Documentation

## 8.43.4.1 friend class TimeTagStreamImpl [friend]

The documentation for this class was generated from the following file:

- [Iterators.h](#)

## 8.44 TimeTagStreamBuffer Class Reference

```
#include <Iterators.h>
```

## Public Member Functions

- void [getOverflows](#) (std::function< unsigned char \*(size\_t)> array\_out)
- void [getChannels](#) (std::function< int \*(size\_t)> array\_out)
- void [getTimestamps](#) (std::function< long long \*(size\_t)> array\_out)
- void [getMissedEvents](#) (std::function< unsigned short \*(size\_t)> array\_out)
- void [getEventTypes](#) (std::function< unsigned char \*(size\_t)> array\_out)

## Public Attributes

- uint64\_t [size](#)
- bool [hasOverflows](#)
- timestamp\_t [tStart](#)
- timestamp\_t [tGetData](#)

## Friends

- class [TimeTagStreamImpl](#)
- class [FileReaderImpl](#)

### 8.44.1 Member Function Documentation

8.44.1.1 void TimeTagStreamBuffer::getChannels ( std::function< int \*(size\_t)> array\_out )

8.44.1.2 void TimeTagStreamBuffer::getEventTypes ( std::function< unsigned char \*(size\_t)> array\_out )

8.44.1.3 void TimeTagStreamBuffer::getMissedEvents ( std::function< unsigned short \*(size\_t)> array\_out )

8.44.1.4 void TimeTagStreamBuffer::getOverflows ( std::function< unsigned char \*(size\_t)> array\_out )

8.44.1.5 void TimeTagStreamBuffer::getTimestamps ( std::function< long long \*(size\_t)> array\_out )

### 8.44.2 Friends And Related Function Documentation

8.44.2.1 friend class FileReaderImpl [friend]

8.44.2.2 friend class TimeTagStreamImpl [friend]

### 8.44.3 Member Data Documentation

8.44.3.1 bool TimeTagStreamBuffer::hasOverflows

8.44.3.2 uint64\_t TimeTagStreamBuffer::size

8.44.3.3 timestamp\_t TimeTagStreamBuffer::tGetData

8.44.3.4 timestamp\_t TimeTagStreamBuffer::tStart

The documentation for this class was generated from the following file:


- [Iterators.h](#)

## 8.45 TriggerOnCountrate Class Reference

Inject trigger events when exceeding or falling below a given count rate within a rolling time window.

```
#include <Iterators.h>
```

Inheritance diagram for TriggerOnCountrate:



## Public Member Functions

- [TriggerOnCountrate](#) ([TimeTaggerBase](#) \*tagger, [channel\\_t](#) input\_channel, double reference\_countrate, double hysteresis, [timestamp\\_t](#) time\_window)  
*constructor of a [TriggerOnCountrate](#)*
- [~TriggerOnCountrate](#) ()
- [channel\\_t](#) [getChannelAbove](#) ()  
*Get the channel number of the above channel.*
- [channel\\_t](#) [getChannelBelow](#) ()  
*Get the channel number of the below channel.*
- [std::vector](#)< [channel\\_t](#) > [getChannels](#) ()  
*Get both virtual channel numbers: [[getChannelAbove](#) (), [getChannelBelow](#) ()].*
- bool [isAbove](#) ()  
*Returns whether the Virtual Channel is currently in the above state.*
- bool [isBelow](#) ()  
*Returns whether the Virtual Channel is currently in the below state.*
- double [getCurrentCountrate](#) ()  
*Get the current count rate averaged within the time\_window.*
- bool [injectCurrentState](#) ()  
*Emit a time-tag into the respective channel according to the current state.*

## Protected Member Functions

- bool [next\\_impl](#) ([std::vector](#)< [Tag](#) > &incoming\_tags, [timestamp\\_t](#) begin\_time, [timestamp\\_t](#) end\_time) override  
*update iterator state*
- void [on\\_start](#) () override  
*callback when the measurement class is started*

## Friends

- class [TriggerOnCountrateImpl](#)

## Additional Inherited Members

### 8.45.1 Detailed Description

Inject trigger events when exceeding or falling below a given count rate within a rolling time window.

Measures the count rate inside a rolling time window and emits tags when a given reference\_countrate is crossed. A [TriggerOnCountrate](#) object provides two virtual channels: The above channel is triggered when the count rate exceeds the threshold (transition from below to above). The below channel is triggered when the count rate falls below the threshold (transition from above to below). To avoid the emission of multiple trigger tags in the transition area, the hysteresis count rate modifies the threshold with respect to the transition direction: An event in the above channel will be triggered when the channel is in the below state and rises to reference\_countrate + hysteresis or above. Vice versa, the below channel fires when the channel is in the above state and falls to the limit of reference\_countrate - hysteresis or below.

The time-tags are always injected at the end of the integration window. You can use the [DelayedChannel](#) to adjust the temporal position of the trigger tags with respect to the integration time window.

The very first tag of the virtual channel will be emitted time\_window after the instantiation of the object and will reflect the current state, so either above or below.

## 8.45.2 Constructor & Destructor Documentation

8.45.2.1 `TriggerOnCountrate::TriggerOnCountrate ( TimeTaggerBase * tagger, channel_t input_channel, double reference_countrate, double hysteresis, timestamp_t time_window )`

constructor of a [TriggerOnCountrate](#)

### Parameters

<i>tagger</i>	Reference to a <a href="#">TimeTagger</a> object.
<i>input_channel</i>	Channel number of the channel whose count rate will control the trigger channels.
<i>reference_countrate</i>	The reference count rate in Hz that separates the <code>above</code> range from the <code>below</code> range.
<i>hysteresis</i>	The threshold count rate in Hz for transitioning to the <code>above</code> threshold state is <code>countrate &gt;= reference_countrate + hysteresis</code> , whereas it is <code>countrate &lt;= reference_countrate - hysteresis</code> for transitioning to the <code>below</code> threshold state. The hysteresis avoids the emission of multiple trigger tags upon a single transition.
<i>time_window</i>	Rolling time window size in ps. The count rate is analyzed within this time window and compared to the threshold count rate.

8.45.2.2 `TriggerOnCountrate::~~TriggerOnCountrate ( )`

## 8.45.3 Member Function Documentation

8.45.3.1 `channel_t TriggerOnCountrate::getChannelAbove ( )`

Get the channel number of the `above` channel.

8.45.3.2 `channel_t TriggerOnCountrate::getChannelBelow ( )`

Get the channel number of the `below` channel.

8.45.3.3 `std::vector<channel_t> TriggerOnCountrate::getChannels ( )`

Get both virtual channel numbers: [[getChannelAbove\(\)](#), [getChannelBelow\(\)](#)].

8.45.3.4 `double TriggerOnCountrate::getCurrentCountrate ( )`

Get the current count rate averaged within the `time_window`.

8.45.3.5 `bool TriggerOnCountrate::injectCurrentState ( )`

Emit a time-tag into the respective channel according to the current state.

Emit a time-tag into the respective channel according to the current state. This is useful if you start a new measurement that requires the information. The function returns whether it was possible to inject the event. The injection is not possible if the Time Tagger is in overflow mode or the time window has not passed yet. The function call is non-blocking.



8.45.3.6 `bool TriggerOnCountrate::isAbove ( )`

Returns whether the Virtual Channel is currently in the `above` state.

8.45.3.7 `bool TriggerOnCountrate::isBelow ( )`

Returns whether the Virtual Channel is currently in the `below` state.

8.45.3.8 `bool TriggerOnCountrate::next_impl ( std::vector< Tag > & incoming_tags, timestamp_t begin_time, timestamp_t end_time ) [override], [protected], [virtual]`

update iterator state

Each [Iterator](#) must implement the `next_impl()` method. The `next_impl()` function is guarded by the update lock.

The backend delivers each [Tag](#) on each registered channel to this callback function.

## Parameters

<i>incoming_tags</i>	block of events
<i>begin_time</i>	earliest event in the block
<i>end_time</i>	begin_time of the next block, not including in this block

## Returns

true if the content of this block was modified, false otherwise

Implements [IteratorBase](#).

8.45.3.9 `void TriggerOnCountrate::on_start ( ) [override], [protected], [virtual]`

callback when the measurement class is started

This function is guarded by the update lock.

Reimplemented from [IteratorBase](#).

## 8.45.4 Friends And Related Function Documentation

8.45.4.1 `friend class TriggerOnCountrateImpl [friend]`

The documentation for this class was generated from the following file:

- [Iterators.h](#)




## Chapter 9

# File Documentation

### 9.1 Iterators.h File Reference

```
#include <algorithm>
#include <array>
#include <assert.h>
#include <deque>
#include <fstream>
#include <iostream>
#include <limits>
#include <list>
#include <map>
#include <memory>
#include <mutex>
#include <queue>
#include <set>
#include <stdint.h>
#include <stdio.h>
#include <unordered_map>
#include <vector>
#include "TimeTagger.h"
```

Include dependency graph for Iterators.h:



### Classes

- class [FastBinning](#)
- class [Combiner](#)

*Combine some channels in a virtual channel which has a tick for each tick in the input channels.*

- class [CountBetweenMarkers](#)

*a simple counter where external marker signals determine the bins*

- class [Counter](#)

- a simple counter on one or more channels*
- class [Coincidences](#)
  - a coincidence monitor for one or more channel groups*
- class [Coincidence](#)
  - a coincidence monitor for one or more channel groups*
- class [Countrate](#)
  - count rate on one or more channels*
- class [DelayedChannel](#)
  - a simple delayed queue*
- class [TriggerOnCountrate](#)
  - Inject trigger events when exceeding or falling below a given count rate within a rolling time window.*
- class [GatedChannel](#)
  - An input channel is gated by a gate channel.*
- class [FrequencyMultiplier](#)
  - The signal of an input channel is scaled up to a higher frequency according to the multiplier passed as a parameter.*
- class [Iterator](#)
  - a simple event queue*
- class [TimeTagStreamBuffer](#)
- class [TimeTagStream](#)
  - access the time tag stream*
- class [Dump](#)
  - dump all time tags to a file*
- class [StartStop](#)
  - simple start-stop measurement*
- class [TimeDifferencesImpl< T >](#)
- class [TimeDifferences](#)
  - Accumulates the time differences between clicks on two channels in one or more histograms.*
- class [Histogram2D](#)
  - A 2-dimensional histogram of time differences. This can be used in measurements similar to 2D NRM spectroscopy.*
- class [TimeDifferencesND](#)
  - Accumulates the time differences between clicks on two channels in a multi-dimensional histogram.*
- class [Histogram](#)
  - Accumulate time differences into a histogram.*
- class [HistogramLogBins](#)
  - Accumulate time differences into a histogram with logarithmic increasing bin sizes.*
- class [Correlation](#)
  - cross-correlation between two channels*
- struct [Event](#)
- class [Scope](#)
- class [SynchronizedMeasurements](#)
  - start, stop and clear several measurements synchronized*
- class [ConstantFractionDiscriminator](#)
  - a virtual CFD implementation which returns the mean time between a raising and a falling pair of edges*
- class [FileWriter](#)
  - compresses and stores all time tags to a file*
- class [FileReader](#)
- class [EventGenerator](#)
  - Generate predefined events in a virtual channel relative to a trigger event.*
- class [CustomMeasurementBase](#)
- class [FlimAbstract](#)
- class [FlimBase](#)

- class [FlimFrameInfo](#)
- class [Flim](#)  
*Fluorescence lifetime imaging.*
- class [SyntheticSingleTag](#)  
*synthetic trigger timetag generator.*

## Macros

- `#define BINNING\_TEMPLATE\_HELPER(fun_name, binner, ...)`  
*FastBinning caller helper.*

## Enumerations

- enum [CoincidenceTimestamp](#) : `uint32_t` { [CoincidenceTimestamp::Last](#) = 0, [CoincidenceTimestamp::Average](#) = 1, [CoincidenceTimestamp::First](#) = 2, [CoincidenceTimestamp::ListedFirst](#) = 3 }
- enum [State](#) { [UNKNOWN](#), [HIGH](#), [LOW](#) }

### 9.1.1 Macro Definition Documentation

#### 9.1.1.1 `#define BINNING\_TEMPLATE\_HELPER( fun_name, binner, ... )`

##### Value:

```
switch (binner.getMode()) {
case FastBinning::Mode::ConstZero:
    fun_name<FastBinning::Mode::ConstZero>(__VA_ARGS__);
    break;
case FastBinning::Mode::Dividend:
    fun_name<FastBinning::Mode::Dividend>(__VA_ARGS__);
    break;
case FastBinning::Mode::PowerOfTwo:
    fun_name<FastBinning::Mode::PowerOfTwo>(__VA_ARGS__);
    break;
case FastBinning::Mode::FixedPoint\_32:
    fun_name<FastBinning::Mode::FixedPoint\_32>(__VA_ARGS__);
    break;
case FastBinning::Mode::FixedPoint\_64:
    fun_name<FastBinning::Mode::FixedPoint\_64>(__VA_ARGS__);
    break;
case FastBinning::Mode::Divide\_32:
    fun_name<FastBinning::Mode::Divide\_32>(__VA_ARGS__);
    break;
case FastBinning::Mode::Divide\_64:
    fun_name<FastBinning::Mode::Divide\_64>(__VA_ARGS__);
    break;
}
```

[FastBinning](#) caller helper.

## 9.1.2 Enumeration Type Documentation

### 9.1.2.1 enum CoincidenceTimestamp : uint32\_t [strong]

type of timestamp for the [Coincidence](#) virtual channel (Last, Average, First, ListedFirst)

Enumerator

***Last***

***Average***

***First***

***ListedFirst***

### 9.1.2.2 enum State

Enumerator

***UNKNOWN***

***HIGH***

***LOW***


## 9.2 TimeTagger.h File Reference

```
#include <condition_variable>
#include <cstdint>
#include <functional>
#include <limits>
#include <list>
#include <map>
#include <memory>
#include <mutex>
#include <set>
#include <stdint.h>
#include <string>
#include <unordered_set>
#include <vector>
```

Include dependency graph for TimeTagger.h:



This graph shows which files directly or indirectly include this file:



## Classes

- class [CustomLogger](#)
- class [TimeTaggerBase](#)
- class [TimeTaggerVirtual](#)  
*virtual [TimeTagger](#) based on dump files*
- class [TimeTagger](#)  
*backend for the [TimeTagger](#).*
- struct [Tag](#)  
*a single event on a channel*
- class [OrderedBarrier](#)
- class [OrderedBarrier::OrderInstance](#)
- class [OrderedPipeline](#)
- class [IteratorBase](#)  
*Base class for all iterators.*

## Macros

- `#define TT\_API __declspec(dllimport)`
- `#define timestamp\_t long long`
- `#define channel\_t int`
- `#define TIMETAGGER\_VERSION "2.9.0"`
- `#define GET\_DATA\_1D(function_name, type, argout, attribute) attribute void function_name(std::function<type *(size_t)> argout)`
- `#define GET\_DATA\_1D\_OP1(function_name, type, argout, optional_type, optional_name, optional_default, attribute) attribute void function_name(std::function<type *(size_t)> argout, optional_type optional_name = optional_default)`
- `#define GET\_DATA\_1D\_OP2(function_name, type, argout, optional_type, optional_name, optional_default, optional_type2, optional_name2, optional_default2, attribute)`
- `#define GET\_DATA\_2D(function_name, type, argout, attribute) attribute void function_name(std::function<type *(size_t, size_t)> argout)`
- `#define GET\_DATA\_2D\_OP1(function_name, type, argout, optional_type, optional_name, optional_default, attribute)`
- `#define GET\_DATA\_2D\_OP2(function_name, type, argout, optional_type, optional_name, optional_default, optional_type2, optional_name2, optional_default2, attribute)`

- `#define GET_DATA_3D(function_name, type, argout, attribute) attribute void function_name(std::function<type *(size_t, size_t, size_t)> argout)`
- `#define LogMessage(level, ...) LogBase(level, __FILE__, __LINE__, false, __VA_ARGS__);`
- `#define ErrorLog(...) LogMessage(LOGGER_ERROR, __VA_ARGS__);`
- `#define WarningLog(...) LogMessage(LOGGER_WARNING, __VA_ARGS__);`
- `#define InfoLog(...) LogMessage(LOGGER_INFO, __VA_ARGS__);`
- `#define LogMessageCensored(level, ...) LogBase(level, __FILE__, __LINE__, true, __VA_ARGS__);`
- `#define ErrorLogCensored(...) LogMessage(LOGGER_ERROR, __VA_ARGS__);`
- `#define WarningLogCensored(...) LogMessage(LOGGER_WARNING, __VA_ARGS__);`
- `#define InfoLogCensored(...) LogMessage(LOGGER_INFO, __VA_ARGS__);`

## Typedefs

- `typedef void(* logger_callback) (LogLevel level, std::string msg)`
- `using _Iterator = IteratorBase`

## Enumerations

- `enum Resolution { Resolution::Standard = 0, Resolution::HighResA = 1, Resolution::HighResB = 2, Resolution::HighResC = 3 }`
- `enum ChannelEdge : int32_t { ChannelEdge::NoFalling = 1 << 0, ChannelEdge::NoRising = 1 << 1, ChannelEdge::NoStandard = 1 << 2, ChannelEdge::NoHighRes = 1 << 3, ChannelEdge::All = 0, ChannelEdge::Rising = 1, ChannelEdge::Falling = 2, ChannelEdge::HighResAll = 4, ChannelEdge::HighResRising = 4 | 1, ChannelEdge::HighResFalling = 4 | 2, ChannelEdge::StandardAll = 8, ChannelEdge::StandardRising = 8 | 1, ChannelEdge::StandardFalling = 8 | 2 }`
- `enum LogLevel { LOGGER_ERROR = 40, LOGGER_WARNING = 30, LOGGER_INFO = 10 }`
- `enum LanguageUsed : std::uint32_t { LanguageUsed::Cpp = 0, LanguageUsed::Python, LanguageUsed::Csharp, LanguageUsed::Matlab, LanguageUsed::Labview, LanguageUsed::Mathematica, LanguageUsed::Unknown = 255 }`
- `enum FrontendType : std::uint32_t { FrontendType::Undefined = 0, FrontendType::WebApp, FrontendType::Firefly, FrontendType::Pyro5RPC, FrontendType::UserFrontend }`
- `enum UsageStatisticsStatus { UsageStatisticsStatus::Disabled, UsageStatisticsStatus::Collecting, UsageStatisticsStatus::CollectingAndUploading }`

## Functions

- `TT_API std::string getVersion ()`
- `TT_API TimeTagger * createTimeTagger (std::string serial="", Resolution resolution=Resolution::Standard)`  
*default constructor factory.*
- `TT_API TimeTaggerVirtual * createTimeTaggerVirtual ()`  
*default constructor factory for the createTimeTaggerVirtual class.*
- `TT_API void setCustomBitFileName (const std::string &bitFileName)`  
*set path and filename of the bitfile to be loaded into the FPGA*
- `TT_API bool freeTimeTagger (TimeTaggerBase *tagger)`  
*free a copy of a TimeTagger reference.*
- `TT_API std::vector< std::string > scanTimeTagger ()`  
*fetches a list of all available TimeTagger serials.*
- `TT_API std::string getTimeTaggerModel (const std::string &serial)`
- `TT_API void setTimeTaggerChannelNumberScheme (int scheme)`



- Configure the numbering scheme for new *TimeTagger* objects.

  - **TT\_API** int `getTimeTaggerChannelNumberScheme` ()  
Fetch the currently configured global numbering scheme.
  - **TT\_API** bool `hasTimeTaggerVirtualLicense` ()  
Check if a license for the *TimeTaggerVirtual* is available.
  - **TT\_API** void `flashLicense` (const std::string &serial, const std::string &license)
  - **TT\_API** std::string `extractLicenseInfo` (const std::string &license)
  - **TT\_API** logger\_callback `setLogger` (logger\_callback callback)  
Sets the notifier callback which is called for each log message.
  - **TT\_API** void `LogBase` (LogLevel level, const char \*file, int line, bool censored, const char \*fmt,...)  
Raise a new log message. Please use the `XXXLog` macro instead.
  - **TT\_API** void `setLanguageInfo` (std::uint32\_t pw, LanguageUsed language, std::string version)  
sets the language being used currently for usage statistics system.
  - **TT\_API** void `setFrontend` (FrontendType frontend)  
sets the frontend being used currently for usage statistics system.
  - **TT\_API** void `setUsageStatisticsStatus` (UsageStatisticsStatus new\_status)  
sets the status of the usage statistics system.
  - **TT\_API** UsageStatisticsStatus `getUsageStatisticsStatus` ()  
gets the status of the usage statistics system.
  - **TT\_API** std::string `getUsageStatisticsReport` ()  
gets the current recorded data by the usage statistics system.

## Variables

- constexpr `channel_t CHANNEL_UNUSED` = -134217728  
Constant for unused channel. Magic `channel_t` value to indicate an unused channel. So the iterators either have to disable this channel, or to choose a default one.
- constexpr `channel_t CHANNEL_UNUSED_OLD` = -1
- constexpr int `TT_CHANNEL_NUMBER_SCHEME_AUTO` = 0
- constexpr int `TT_CHANNEL_NUMBER_SCHEME_ZERO` = 1
- constexpr int `TT_CHANNEL_NUMBER_SCHEME_ONE` = 2
- constexpr `ChannelEdge TT_CHANNEL_RISING_AND_FALLING_EDGES` = `ChannelEdge::All`
- constexpr `ChannelEdge TT_CHANNEL_RISING_EDGES` = `ChannelEdge::Rising`
- constexpr `ChannelEdge TT_CHANNEL_FALLING_EDGES` = `ChannelEdge::Falling`

### 9.2.1 Macro Definition Documentation

9.2.1.1 `#define channel_t int`

9.2.1.2 `#define ErrorLog( ... ) LogMessage(LOGGER_ERROR, __VA_ARGS__);`

9.2.1.3 `#define ErrorLogCensored( ... ) LogMessage(LOGGER_ERROR, __VA_ARGS__);`

9.2.1.4 `#define GET_DATA_1D( function_name, type, argout, attribute ) attribute void function_name(std::function<type *(size_t)> argout)`

This are the default wrapper functions without any overloadings.

9.2.1.5 **#define** GET\_DATA\_1D\_OP1( *function\_name*, *type*, *argout*, *optional\_type*, *optional\_name*, *optional\_default*, *attribute* ) **attribute** void *function\_name*(std::function<type \*(size\_t)> *argout*, *optional\_type* *optional\_name* = *optional\_default*)

9.2.1.6 **#define** GET\_DATA\_1D\_OP2( *function\_name*, *type*, *argout*, *optional\_type*, *optional\_name*, *optional\_default*, *optional\_type2*, *optional\_name2*, *optional\_default2*, *attribute* )

**Value:**

```
attribute void function_name(std::function<type *(size_t)> argout, optional_type optional_name =
    optional_default, \
        optional_type2 optional_name2 = optional_default2)
```

9.2.1.7 **#define** GET\_DATA\_2D( *function\_name*, *type*, *argout*, *attribute* ) **attribute** void *function\_name*(std::function<type \*(size\_t, size\_t)> *argout*)

9.2.1.8 **#define** GET\_DATA\_2D\_OP1( *function\_name*, *type*, *argout*, *optional\_type*, *optional\_name*, *optional\_default*, *attribute* )

**Value:**

```
attribute void function_name(std::function<type *(size_t, size_t)> argout,
    \
        optional_type optional_name = optional_default)
```

9.2.1.9 **#define** GET\_DATA\_2D\_OP2( *function\_name*, *type*, *argout*, *optional\_type*, *optional\_name*, *optional\_default*, *optional\_type2*, *optional\_name2*, *optional\_default2*, *attribute* )

**Value:**

```
attribute void function_name(std::function<type *(size_t, size_t)> argout,
    \
        optional_type optional_name = optional_default,
    \
        optional_type2 optional_name2 = optional_default2)
```

9.2.1.10 **#define** GET\_DATA\_3D( *function\_name*, *type*, *argout*, *attribute* ) **attribute** void *function\_name*(std::function<type \*(size\_t, size\_t, size\_t)> *argout*)

9.2.1.11 **#define** InfoLog( ... ) LogMessage(LOGGER\_INFO, \_\_VA\_ARGS\_\_);

9.2.1.12 **#define** InfoLogCensored( ... ) LogMessage(LOGGER\_INFO, \_\_VA\_ARGS\_\_);

9.2.1.13 **#define** LogMessage( *level*, ... ) LogBase(level, \_\_FILE\_\_, \_\_LINE\_\_, false, \_\_VA\_ARGS\_\_);

9.2.1.14 **#define** LogMessageCensored( *level*, ... ) LogBase(level, \_\_FILE\_\_, \_\_LINE\_\_, true, \_\_VA\_ARGS\_\_);

9.2.1.15 **#define** timestamp\_t long long

9.2.1.16 `#define TIMETAGGER_VERSION "2.9.0"`

9.2.1.17 `#define TT_API __declspec(dllimport)`

9.2.1.18 `#define WarningLog( ... ) LogMessage(LOGGER_WARNING, __VA_ARGS__);`

9.2.1.19 `#define WarningLogCensored( ... ) LogMessage(LOGGER_WARNING, __VA_ARGS__);`

## 9.2.2 Typedef Documentation

9.2.2.1 `using _Iterator = IteratorBase`

9.2.2.2 `typedef void(* logger_callback) (LogLevel level, std::string msg)`

## 9.2.3 Enumeration Type Documentation

9.2.3.1 `enum ChannelEdge : int32_t [strong]`

Enum for filtering the channel list returned by `getChannelList`.

Enumerator

***NoFalling***  
***NoRising***  
***NoStandard***  
***NoHighRes***  
***All***  
***Rising***  
***Falling***  
***HighResAll***  
***HighResRising***  
***HighResFalling***  
***StandardAll***  
***StandardRising***  
***StandardFalling***

9.2.3.2 `enum FrontendType : std::uint32_t [strong]`

Enumerator

***Undefined***  
***WebApp***  
***Firefly***  
***Pyro5RPC***  
***UserFrontend***

### 9.2.3.3 enum LanguageUsed : std::uint32\_t [strong]

Enumerator

***Cpp***  
***Python***  
***Csharp***  
***Matlab***  
***Labview***  
***Mathematica***  
***Unknown***

### 9.2.3.4 enum LogLevel

Enumerator

***LOGGER\_ERROR***  
***LOGGER\_WARNING***  
***LOGGER\_INFO***

### 9.2.3.5 enum Resolution [strong]

This enum selects the high resolution mode of the Time Tagger series. If any high resolution mode is selected, the hardware will combine 2, 4 or even 8 input channels and average their timestamps. This results in a discretization jitter improvement of factor  $\sqrt{N}$  for  $N$  combined channels. The averaging is implemented before any filter, buffer or USB transmission. So all of those features are available with the averaged timestamps. Because of hardware limitations, only fixed combinations of channels are supported:

- HighResA: 1 : [1,2], 3 : [3,4], 5 : [5,6], 7 : [7,8], 10 : [10,11], 12 : [12,13], 14 : [14,15], 16 : [16,17], 9, 18
- HighResB: 1 : [1,2,3,4], 5 : [5,6,7,8], 10 : [10,11,12,13], 14 : [14,15,16,17], 9, 18
- HighResC: 5 : [1,2,3,4,5,6,7,8], 14 : [10,11,12,13,14,15,16,17], 9, 18 The inputs 9 and 18 are always available without averaging. The number of channels available will be limited to the number of channels licensed.

Enumerator

***Standard***  
***HighResA***  
***HighResB***  
***HighResC***

### 9.2.3.6 enum UsageStatisticsStatus [strong]

Enumerator

***Disabled***  
***Collecting***  
***CollectingAndUploading***

## 9.2.4 Function Documentation

### 9.2.4.1 TT\_API TimeTagger\* createTimeTagger ( std::string serial = " ", Resolution resolution = Resolution::Standard )

default constructor factory.

## Parameters

<i>serial</i>	serial number of FPGA board to use. if empty, the first board found is used.
<i>resolution</i>	enum for how many channels shall be grouped.

## See also

[Resolution](#) for details

## 9.2.4.2 TT\_API TimeTaggerVirtual\* createTimeTaggerVirtual ( )

default constructor factory for the createTimeTaggerVirtual class.

## 9.2.4.3 TT\_API std::string extractLicenseInfo ( const std::string &amp; license )

Parses the binary license and return a human readable information about this license.

## Parameters

<i>license</i>	the binary license, encoded as a hexadecimal string
----------------	---

## Returns

a human readable string containing all information about this license

## 9.2.4.4 TT\_API void flashLicense ( const std::string &amp; serial, const std::string &amp; license )

Update the license on the device. Updated license may be fetched by getRemoteLicense. The Time Tagger must not be instancated while updating the license.

## Parameters

<i>serial</i>	the serial of the device to update the license. Must not be empty
<i>license</i>	the binary license, encoded as a hexadecimal string

## 9.2.4.5 TT\_API bool freeTimeTagger ( TimeTaggerBase \* tagger )

free a copy of a [TimeTagger](#) reference.

## Parameters

<i>tagger</i>	the <a href="#">TimeTagger</a> reference to free
---------------	--

#### 9.2.4.6 TT\_API int getTimeTaggerChannelNumberScheme ( )

Fetch the currently configured global numbering scheme.

Please see [setTimeTaggerChannelNumberScheme\(\)](#) for details. Please use [TimeTagger::getChannelNumberScheme\(\)](#) to query the actual used numbering scheme, this function here will just return the scheme a newly created [TimeTagger](#) object will use.

#### 9.2.4.7 TT\_API std::string getTimeTaggerModel ( const std::string & *serial* )

#### 9.2.4.8 TT\_API std::string getUsageStatisticsReport ( )

gets the current recorded data by the usage statistics system.

Use this function to see what data has been collected so far and what will be sent to Swabian Instruments if 'CollectingAndUploading' is enabled. All data is pseudonymize.

#### Note

if no data has been collected or due to a system error, the database was corrupted, it will return an error. else it will be a database in json format.

#### Returns

the current recorded data by the usage statistics system.

#### 9.2.4.9 TT\_API UsageStatisticsStatus getUsageStatisticsStatus ( )

gets the status of the usage statistics system.

#### Returns

the current status of the usage statistics system.

#### 9.2.4.10 TT\_API std::string getVersion ( )

#### 9.2.4.11 TT\_API bool hasTimeTaggerVirtualLicense ( )

Check if a license for the [TimeTaggerVirtual](#) is available.

#### 9.2.4.12 TT\_API void LogBase ( LogLevel *level*, const char \* *file*, int *line*, bool *censored*, const char \* *fmt*, ... )

Raise a new log message. Please use the XXXLog macro instead.

#### 9.2.4.13 TT\_API std::vector<std::string> scanTimeTagger ( )

fetches a list of all available [TimeTagger](#) serials.

This function may return serials blocked by other processes or already disconnected some milliseconds later.

#### 9.2.4.14 TT\_API void setCustomBitFileName ( const std::string & *bitFileName* )

set path and filename of the bitfile to be loaded into the FPGA

For debugging/development purposes the firmware loaded into the FPGA can be set manually with this function. To load the default bitfile set bitFileName = ""

## Parameters

<i>bitFileName</i>	custom bitfile to use for the FPGA.
--------------------	-------------------------------------

9.2.4.15 TT\_API void setFrontend ( FrontendType *frontend* )

sets the frontend being used currently for usage statistics system.

## Parameters

<i>frontend</i>	the frontend currently being used.
-----------------	------------------------------------

9.2.4.16 TT\_API void setLanguageInfo ( std::uint32\_t *pw*, LanguageUsed *language*, std::string *version* )

sets the language being used currently for usage statistics system.

## Parameters

<i>pw</i>	password for authporization to change the language.
<i>language</i>	programming language being used.
<i>version</i>	version of the programming langugae being used.

9.2.4.17 TT\_API logger\_callback setLogger ( logger\_callback *callback* )

Sets the notifier callback which is called for each log message.

## Returns

The old callback

If this function is called with nullptr, the default callback will be used.

9.2.4.18 TT\_API void setTimeTaggerChannelNumberScheme ( int *scheme* )

Configure the numbering scheme for new [TimeTagger](#) objects.

## Parameters

<i>scheme</i>	new numbering scheme, must be TT_CHANNEL_NUMBER_SCHEME_AUTO, TT_CHANNEL_NUMBER_SCHEME_ZERO or TT_CHANNEL_NUMBER_SCHEME_ONE
---------------	--

This function sets the numbering scheme for newly created [TimeTagger](#) objects. The default value is `_AUTO`.

Note: [TimeTagger](#) objects are cached internally, so the scheme should be set before the first call of [createTimeTagger\(\)](#).

`_ZERO` will typically allocate the channel numbers 0 to 7 for the 8 input channels. 8 to 15 will be allocated for the corresponding falling events.

`_ONE` will typically allocate the channel numbers 1 to 8 for the 8 input channels. -1 to -8 will be allocated for the corresponding falling events.

`_AUTO` will choose the scheme based on the hardware revision and so based on the printed label.

#### 9.2.4.19 `TT_API void setUsageStatisticsStatus ( UsageStatisticsStatus new_status )`

sets the status of the usage statistics system.

This functionality allows configuring the usage statistics system.

##### Parameters

<code>new_status</code>	new status of the usage statistics system.
-------------------------	--

### 9.2.5 Variable Documentation

#### 9.2.5.1 `constexpr channel_t CHANNEL_UNUSED = -134217728`

Constant for unused channel. Magic `channel_t` value to indicate an unused channel. So the iterators either have to disable this channel, or to choose a default one.

This value changed in version 2.1. The old value -1 aliases with falling events. The old value will still be accepted for now if the old numbering scheme is active.

#### 9.2.5.2 `constexpr channel_t CHANNEL_UNUSED_OLD = -1`

#### 9.2.5.3 `constexpr ChannelEdge TT_CHANNEL_FALLING_EDGES = ChannelEdge::Falling`

#### 9.2.5.4 `constexpr int TT_CHANNEL_NUMBER_SCHEME_AUTO = 0`

Allowed values for [setTimeTaggerChannelNumberScheme\(\)](#).

`_ZERO` will typically allocate the channel numbers 0 to 7 for the 8 input channels. 8 to 15 will be allocated for the corresponding falling events.

`_ONE` will typically allocate the channel numbers 1 to 8 for the 8 input channels. -1 to -8 will be allocated for the corresponding falling events.

`_AUTO` will choose the scheme based on the hardware revision and so based on the printed label.

#### 9.2.5.5 `constexpr int TT_CHANNEL_NUMBER_SCHEME_ONE = 2`

#### 9.2.5.6 `constexpr int TT_CHANNEL_NUMBER_SCHEME_ZERO = 1`

#### 9.2.5.7 `constexpr ChannelEdge TT_CHANNEL_RISING_AND_FALLING_EDGES = ChannelEdge::All`

#### 9.2.5.8 `constexpr ChannelEdge TT_CHANNEL_RISING_EDGES = ChannelEdge::Rising`